

Oct 23, 2019

## iPSC editing with TALENs

DOI

[dx.doi.org/10.17504/protocols.io.8dahs2e](https://dx.doi.org/10.17504/protocols.io.8dahs2e)

iPSC editing  
with TALENs

Ruilin Tian<sup>1</sup>, Jason Hong<sup>1</sup>, Martin Kampmann<sup>1</sup>

<sup>1</sup>University of California, San Francisco

Neurodegeneration Met...

KampmannLab



**Martin Kampmann**

University of California, San Francisco

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.8dahs2e](https://dx.doi.org/10.17504/protocols.io.8dahs2e)

**Protocol Citation:** Ruilin Tian, Jason Hong, Martin Kampmann 2019. iPSC editing with TALENs. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.8dahs2e>

**Manuscript citation:**

Tian et al (2019). CRISPR Interference-Based Platform for Multimodal Genetic Screens in Human iPSC-Derived Neurons. Neuron pii: S0896-6273(19)30640-3. [Epub ahead of print] PubMed PMID: 31422865.

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 17, 2019

**Last Modified:** October 23, 2019

**Protocol Integer ID:** 28802

**Keywords:** iPSC, TALENs, CRISPR interference, CRISPRi, CROP-seq, Perturb-Seq, essential genes, functional genomics, high-content microscopy, neuron, single-cell RNA sequencing, stem cell

## Attachments



Editing\_iPSCs\_with\_T...

58KB

## Materials

### MATERIALS

⊗ pZT-C13-L1 **addgene Catalog #62196**

⊗ pZT-C13-R1 **addgene Catalog #62197**

⊗ DPBS (no Ca, no Mg) **Thermofisher Catalog #14190144**

⊗ Essential 8™ Medium **Gibco - Thermo Fisher Scientific Catalog #A1517001**

⊗ StemPro™ Accutase™ Cell Dissociation Reagent **Thermo Fisher Scientific Catalog #A1110501**

⊗ Lipofectamine™ Stem Transfection Reagent **Thermo Fisher Scientific Catalog #STEM00008**

⊗ StemFlex™ Medium **Thermo Fisher Scientific Catalog #A3349401**

⊗ Matrigel **Corning Catalog #356231**

⊗ KnockOut™ DMEM **Thermo Fisher Scientific Catalog #10829018**

⊗ Opti-MEM™ I Reduced Serum Medium **Thermo Fisher Scientific Catalog #31985070**

⊗ Y-27632 dihydrochloride (Rock Inhibitor) **R&D Systems Catalog #1254/10**



- BCL-XL plasmid (pEF1-BCL-XL-wpre-polyA P1102): Gift from Xiaobing Zhang, described in PMID: 30239926

## Safety warnings

! Please see SDS (Safety Data Sheet) for hazards and safety warnings.







## Pre-coating 6 well plate

- 1 Pre-coat 6 well plate with Matrigel (**diluted 1:50** with Knockout DMEM) by adding  1 mL to each well.
- 2 Let it sit in incubator for at least  00:20:00 .



## Passaging and Plating iPSC Cells

- 3 Grow iPSCs until they are ~ **85 %** confluent in one well of a 6 well plate or larger format.
- 4 Remove old medium and wash iPSCs with DPBS 1x.
- 5 Add Accutase for  00:03:00 –  00:05:00 to lift the cells.
- 6 Singularize the cells by gently pipetting them up and down several times.
- 7 In a 15 ml conical tube, add DPBS and then add the lifted cells in Accutase.
- 8 Spin down at  200 x g for  00:05:00 .
- 9 Aspirate supernatant and resuspend cells in Stemflex.
- 10 Perform cell count.
- 11 Remove the Matrigel in the pre-coated plate and add appropriate amount of Stemflex medium with Rock inhibitor (1000x).
- 12 Re-seed 0.5 M cells into a 6 well plate so that it can reach ~ **60 %** confluency the next day.





13 Put the plate in incubator and culture overnight.



## Transfection

14 The following day, change the media with new E8 media.

15 Prepare the following mixes:

*Tube 1:*

- 100  $\mu$ L Opti-Mem
- 10  $\mu$ L Lipofectamine Stem Reagent:

*Tube 2:*

- 100  $\mu$ L Opti-Mem
- DNA mix ( 0.5  $\mu$ g – 5.0  $\mu$ g total)

Optimized TALENS ratio

- 1.5  $\mu$ g Your Plasmid of choice
- 0.75  $\mu$ g TALENS L
- 0.75  $\mu$ g TALENS R
- 0.3  $\mu$ g BCL-XL

16 Add tube 2 to tube 1 and mix well.



17 Incubate the mixture for 00:10:00 .



18 Add the whole mixture to your cells in the 6 well plate.

19 Incubate overnight and change the media to Stemflex the next day.



20 If confluent, passage and expand them into 10 cm plate.



21 When confluent in 10 cm plate, continue to selection protocol.