

Feb 27, 2019

iPSC Freezing

 In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.x8zfrx6

Celeste M M. Karch¹, Rita Martinez¹, Jacob Marsh¹


¹Washington University in St Louis

Neurodegeneration Method Development Community
Tech. support email: ndcn-help@chanzuckerberg.com



Celeste M M. Karch

Washington University in St Louis



OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.x8zfrx6

Protocol Citation: Celeste M M. Karch, Rita Martinez, Jacob Marsh 2019. iPSC Freezing. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.x8zfrx6>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 18, 2019

Last Modified: February 27, 2019

Protocol Integer ID: 20473

Guidelines

This protocols is part of the **Screening Edited iPSC Clones collection**.

Materials

STEP MATERIALS


 Accutase™ Cell Dissociation Reagent **Gibco - Thermo Fisher Scientific Catalog #A1110501**

Protocol materials

 Accutase™ Cell Dissociation Reagent **Gibco - Thermo Fisher Scientific Catalog #A1110501**

 Accutase™ Cell Dissociation Reagent **Gibco - Thermo Fisher Scientific Catalog #A1110501**

Safety warnings

 Please refer to the SDS (Safety Data Sheet) for information about hazards, and to obtain advice on safety precautions.





1 Aspirate media

2 Gently wash cells with 1x PBS.

Note

Use 2-3 mL per well in 6 well plate.

3 Add Accutase (Gibco A11105-01) directly to the cells and incubate at  37 °C for 3-4 minutes.  00:03:00

Note

Individual donor cell lines exhibit variable sensitivity to accutase-mediated dissociation. Thus, monitor cells closely to determine when single cell dissociation is achieved.

Note

For a 6 well plate, add 0.75-1 mL per well.
For a 10cm² dish, add 3 mL.









Accutase™ Cell Dissociation Reagent **Gibco - Thermo Fisher**
Scientific Catalog #A1110501

4 Tap dish to aid in dislocation of cells.

5 Add DMEM/F12 directly to cells.

**Note**

- For a 6 well plate, add  2 mL -  4 mL per well.
- For a 10cm² dish, add  9 mL
- If cells remain attached, use a cell scraper to gently dislodge cells (apply gentle pressure and use 1-2 passes to remove cells)

- 6 Collect cells in conical tube (15mL/50mL depending on volume).
- 7 Add  2 mL -  5 mL DMEM/F12 to dish to remove all cells from the dish and add to conical tube.
- 8 Centrifuge cells at 750 rpm for  00:03:00 at room temperature.
- 9 Carefully aspirate supernatant.

Note



To avoid aspirating cell pellet, it is OK to leave a small amount of media (0.5-1mL).

- 10 Resuspend cell pellet with mTesR1 (No Rock Inhibitor).
 - Use volume appropriate for freezing.
 - Assume 1 mL per cryovial total and add ½ total volume of mTesR1.
 - Pipet cells 1-2 times only to preserve cell clumps.

Note

Example: to freeze 10 tubes, you will need 10 mL total and will add 5 mL mTesR1 to cell pellet (and 5 mL of 2x Freezing Media below).



- 11 Add an equal volume of cold 2x Freezing Media (20% DMSO, FBS). Pipet cells 1 time only to preserve cell clumps.
- 12 Transfer cell suspension to pre-labeled cryovials (1 mL per cryovial).
 - Ensure that cryovials are labeled with the following
 - Cell Type
 - Line Name
 - Passage #
 - Date
 - Your Name
- 13 Freeze vials at  -80 °C in foam racks for 48-72 hours.  48:00:00
- 14 Transfer vials to liquid nitrogen for long-term storage.