

Feb 26, 2019

iPSC gDNA Extraction: For Screening Edited Clones

 In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.x8afrse

Celeste M M. Karch¹, Rita Martinez¹, Jacob Marsh¹


¹Washington University in St Louis

Neurodegeneration Method Development Community
Tech. support email: ndcn-help@chanzuckerberg.com



Celeste M M. Karch

Washington University in St Louis



OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.x8afrse

Protocol Citation: Celeste M M. Karch, Rita Martinez, Jacob Marsh 2019. iPSC gDNA Extraction: For Screening Edited Clones. [protocols.io](https://dx.doi.org/10.17504/protocols.io.x8afrse) <https://dx.doi.org/10.17504/protocols.io.x8afrse>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 16, 2019

Last Modified: February 26, 2019

Protocol Integer ID: 20450

Attachments



Comprehensive


Genomi...

31KB


Guidelines

This protocols is part of the Screening Edited iPSC Clones collection.

Safety warnings


 Please refer to the SDS (Safety Data Sheet) for information about hazards, and to obtain advice on safety precautions.





1 Spin 96 well PCR plate at 3800 rpm for  00:30:00 at  4 °C .

2 Pipet off supernatant.


Note


If worried about removing cells, transfer supernatant with multichannel pipet into new PCR plates and store at  -20 °C until DNA extraction is complete.


3 To cell pellet (often not visible), add  50 µL QuickExtract DNA Solution (**Epicentre Technologies QE09050**).

4 Vortex plate for  00:00:15 .

5 Incubate plate at  65 °C for  00:06:00 .

6 Vortex plate for  00:00:15 .

7 Incubate plate in thermocycler at  98 °C for  00:02:00 .

8 Place plate in  -20 °C for storage, until ready to use for further screening experiments .