

Feb 17, 2025

IPSC Neuron Development

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl4wrzovo5/v1>

Chelsie Steele¹, Yujie Fan¹

¹Caltech



Chelsie Steele

Caltech

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4wrzovo5/v1>

Protocol Citation: Chelsie Steele, Yujie Fan 2025. IPSC Neuron Development . protocols.io

<https://dx.doi.org/10.17504/protocols.io.3byl4wrzovo5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 14, 2025



Last Modified: February 17, 2025

Protocol Integer ID: 118317

Keywords: ipsc neuron development, ipsc, development

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-020495

Abstract

IPSC Neuron Development

Materials

E8 medium (Thermo Fisher, A1517001)

Vitronectin (Thermo Fisher, A14700)

Matrigel (Corning, 354262)

Accutase (Fisher Scientific, A1110501)

E6 media - (ThermoFisher, A1516401)

Y-27632 - (Bio-Techne, 1254/10)

SB (R&D #1614)

LDN (Stemgent, 04-0074)

XAV (R&D , 3748)

N2 supplement (Stem Cell Technologies, 07156)

B27 (Fisher Scientific,17-504-044)

NAEE (Fisher Scientific,11-140-050)

Glutamax (ThermoFisher, 35050061)

PenStrep (ThermoFisher, 15140122)

Poly-L-ornithine (Sigma Aldrich, P4957-50ML)

GDNF (Recombinant Human Glial Derived Neurotrophic Factor) (Peprotech, 450-10)

BDNF (Recombinant Human/Mouse BDNF) (R&D, 248-BDB)

AA (Ascorbic Acid) (Sigma, 4034-100g)

dbCAMP(Cyclic Monophosphate Sodium Salt)(Sigma,D0627)



before day -1

- 1 hESCs were cultured on 6-well coated with Vitronectin (1:100 diluted in DPBS and coated in cold room overnight) and maintained in E8 medium (Thermo Fisher, A1517001).
- 2 The E8 medium was changed daily, and the cells were split every 3-5 days at 70-85% confluence. then seeded into the matrigel plates for day -1

Day-1 (Cortical differentiation (H9 /E8-E6))

1d

- 3 Thaw a 10mL Matrigel container overnight on ice, then aliquot & freeze. 12h
- 4 Dilute Matrigel 1:100 in cold DMEM/F12 and coat a 24 well with 1/2ml per well overnight at 4°C 12h
- 5 Aspirate E8 medium from hESCs and wash 1X with PBS. Add sufficient Accutase to cover the surface of the cells. In a 6-well plate, this is about 1mL of Accutase.
- 6 Incubate the plate at 37°C for 10 minutes. Then, triturate the cells in Accutase 10X using a P1000 to break them into single cells. Collect the cells into a 15mL tube, leaving one well for maintaining.
- 7 Spin down the 15mL canonical tube at 300g for 5 minutes. Aspirate the supernatant, and re-suspend the pellet in 1mL of E6 medium.
- 8 Count the cells by taking out 10uL of the re-suspended cells and mixing it thoroughly with 10uL of trypan blue in a separate 1.5mL Eppendorf tube. Load 10uL of this mixture onto a counting plate to count the cells. When you count, it will give you a concentration. Based on this concentration, plate approximately 100,000 cells/cm². For instance, in a 6-well plate, this is approximately 1,000,000 cells per well, and so you would need $N = 1,000,000 \text{ cells/well} \times 6 \text{ wells} = 6 \text{ million cells}$.
- 9 Using $V = N/C$, calculate the volume, V , needed from the 15mL tube of re-suspended cells for seeding. In a 50mL falcon tube, **E6 media+ Y-27632** (10 μ M, 1:1000) + the appropriate volume (V) from the 15mL tube.
- 10 Add 2mL of the mixture from the 50mL plate onto the geltrex or matrigel new well plates.



Day 0 to day 3

- 11 Check the cell confluence. The cells need to be 100% confluent before initiating differentiation.
- 12 Feed cells daily with E6 media supplemented as follows with: SB (10 μ M - 1:1000 of the stock, stock at 10mM); LDN (100nM - 1:5000 of the stock, stock at 500 μ M); XAV (500nM 1:2000 of the stock, stock at 1mM)

Day 4 to day 10

- 13 Feed cells daily with E6 media supplemented with: SB (10 μ M - 1:1000 of the stock, stock at 10mM);
LDN (100nM - 1:5000 of the stock, stock at 500 μ M)

Day 10 to day 20

- 14 Feed cells daily with Neurobasal media supplemented with the following: N2(1:100) + B27(1:50)+NAEE+GluMax (1:1000)+P/S (1:1000)

D18 to D19 Prepare PO/LM/FN plates

1d 16h

- 15 Prep the plates with a poly-L-ornithine solution diluted in DPBS for 24 hours. 1d
- 16 Aspirate the poly-L-ornithine solution, and wash each well twice with 0.5 mL DPBS.
- 17 Add 12 μ L Laminin (LM, 1 μ g/mL) and 12 μ L Fibronectin (FN, 1 μ g /mL) to 12 mL DPBS. Mix and add 0.5 mL of the mixture to each well. Incubate at 37°C for $\geq$ 16 h. 16h

Day 20 - Passage Cells

20m

- 18 Wash each well with DPBS then add Accutase for 20 min / 37° C in the incubator. Follow the steps from 6-9 to re-plate the cells. However, plate the cells at a density of 150, 000 cells/cm² in the pre-prepared PO/LM/FN plates, and use the following media mixture to feed the cells instead of E6: Neurobasal media supplemented with N2(1:100) + B27(1:50)+NAEE+GluMax (1:1000)+P/S (1:1000)+ ROCK inhibitor (1:1000, Y-27632). 20m



Day 21

- 19 Feed the cells on with Neurobasal basic media supplemented with: N2(1:100) + B27(1:50)+NAEE+GluMax (1:1000)+P/S (1:1000)+ ROCK inhibitor (1:1000, Y-27632) + GDNF (10ng/ml, 1X1000 dilution from 10ug/ml stock)+BDNF (10ng/ml, 1X1000 dilution from 10ug/ml stock) + AA (100uM, 1X1000 dilution from 100mM stock) + dbCAMP (100uM, 1×1000 dilution from 100mM stock).

Day 30 - end point

- 20 Maintain in the same media mentioned in day 21 changing ½ of the media every 5 days approx.