

Feb 27, 2019

# iPSC Restriction Digest: For Screening Edited Clones

 In 1 collection

DOI

[dx.doi.org/10.17504/protocols.io.x8sfrwe](https://dx.doi.org/10.17504/protocols.io.x8sfrwe)

Celeste M M. Karch<sup>1</sup>, Rita Martinez<sup>1</sup>, Jacob Marsh<sup>1</sup>


<sup>1</sup>Washington University in St Louis

Neurodegeneration Method Development Community  
Tech. support email: [ndcn-help@chanzuckerberg.com](mailto:ndcn-help@chanzuckerberg.com)



**Celeste M M. Karch**

Washington University in St Louis



## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.x8sfrwe>

**Protocol Citation:** Celeste M M. Karch, Rita Martinez, Jacob Marsh 2019. iPSC Restriction Digest: For Screening Edited Clones. [protocols.io https://dx.doi.org/10.17504/protocols.io.x8sfrwe](https://dx.doi.org/10.17504/protocols.io.x8sfrwe)

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** February 17, 2019

**Last Modified:** February 27, 2019

**Protocol Integer ID:** 20466

**Keywords:** ipsc restriction digest, screening edited clone, ipsc, edited clone

## Attachments



Comprehensive

Genomi...

31KB

## Guidelines

This protocols is part of the [Screening Edited iPSC Clones collection](#).

## Materials

### STEP MATERIALS

 CutSmart Buffer - 5.0 ml **New England Biolabs Catalog #B7204S**


## Protocol materials

 CutSmart Buffer - 5.0 ml **New England Biolabs Catalog #B7204S**

 CutSmart Buffer - 5.0 ml **New England Biolabs Catalog #B7204S**

## Troubleshooting

## Safety warnings

 Please refer to the SDS (Safety Data Sheet) for information about hazards, and to obtain advice on safety precautions.



## Before start

Identify the appropriate enzyme for the cell line/mutation in question.

a. A mutation or correction may either introduce or destroy a cut site.

i. We recommend a screening strategy where the editing event introduces a new enzyme cut site as this tends to be more specific

ii. An editing event that destroys an enzyme cut site can be used but tends to result in a high rate of false positives.

If no cut sites exist for mutation/correction in question, T7 assay should be performed.


After confirming the presence of the appropriate sized PCR product, you can move on to performing a restriction digest.



- 1 Prepare the following reagents in a strip cap tube with the previously identified enzyme.

Restriction Digest Protocol		
Reagent	Volum e	# of rxns
Buffer (Cut Smart/ NEB)	2 $\mu$ l	
DNA from PCR	17.75 $\mu$ l	
Enzyme (specific)	0.25 $\mu$ l	
Total	20 $\mu$ l	




 CutSmart Buffer - 5.0 ml **New England Biolabs Catalog #B7204S**

- 2 Incubate reaction at temperature ideal for the enzyme being used (e.g. 37°C , 42°C, 25°C, etc.)
- 3 Incubate at appropriate temperature for 2-3 hours.  02:00:00
- 4 Run the enzyme reaction on a gel to visualize product.
- 5 To make the gel, combine an appropriate amount of Agarose, TBE and Ethidium Bromide using the following guidelines.
  - a. The 2% gel will be cast in one of the following ways:

	15×15 cast	15×25 cast
Agarose	1.5 g	3.0 g
TBE	75 mL	150 mL
Ethidiu m Bromide	3.75 uL	7.5 uL

- b. Combine Agarose and TBE in an appropriately sized flask and microwave until Agarose is completely dissolved. Swirl intermittently during heating.



- c. Once completely dissolved add appropriate amount of Ethidium Bromide to flask and swirl until dispersed evenly.
- 6 Pour gel from flask into casting tray (be sure to add appropriate amount of combs to casting tray).
- 7 Let sit for 30-40 minutes, or until firm.  00:30:00
- 8 Place gel cast into the gel rig apparatus.
- 9 Load samples.
- 10 Load 50 bp ladder.
- 11 Place lid on gel rig apparatus.
- 12 Run gel at 150 volts for  01:30:00 (checking at  01:00:00 to ensure samples have not run too far or off the gel).
- 13 Turn off gel rig apparatus and remove cast.
- 14 Blot off excess TBE from cast.
- 15 Analyze gel images and select potentially edited clones based on banding patterns.