

Nov 13, 2023

Version 1

iPSCs Maintenance and Banking V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov1qd5ogr2/v1>

Mahmoud ElAchwah¹, Sol Diaz de Leon¹, Zhiping Pang¹

¹Rutgers University



Mahmoud ElAchwah

Rutgers University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov1qd5ogr2/v1>

Protocol Citation: Mahmoud ElAchwah, Sol Diaz de Leon, Zhiping Pang 2023. iPSCs Maintenance and Banking. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.ewov1qd5ogr2/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 15, 2023



Last Modified: November 13, 2023

Protocol Integer ID: 87845

Keywords: pluripotent stem cell, induced pluripotent stem cell, ipscs maintenance, ipsc, cell

Abstract

This protocols offers a thorough description of the maintenance and banking of induced pluripotent stem cells.

Materials

Matrigel:

Matrigel catalog number: **354234**

Matrigel stock concentration: **9-11 µg/µL**

Matrigel Stock aliquot volume: **400 µL**

Matrigel LoT number: **13823002**

Matrigel working concentration: **90-110 µg/mL**

Culture Antibiotic:

Primocin catalog number: **ant-pm-1**

Primocin stock concentration: **50 mg/mL**

Primocin stock aliquot volume: **1 mL**

Primocin LoT number: **NA**

Primocin working concentration: **100 µg/mL**

CEPT:

Chroman catalog number: **HY-15392**

Chroman stock concentration: **0.5 mM**

Chroman working concentration: **50 nM**

Emricasan catalog number: **S7775**

Emricasan stock concentration: **50 mM**

Emricasan working concentration: **5000 nM**

Polyamine supplement catalog number: **P8483**

Polyamine supplement stock concentration: **NA**

Polyamine supplement working concentration: **NA**

Trans-ISRIB catalog number: **5284**

Trans-ISRIB stock concentration: **7 mM**

Trans-IRIB working concentration: **700 nM**

Culture Media:

MEM catalog number: **51200-038**

MEM LoT number: **2451155**

Accutase catalog number: **AT-104**

Acutase LoT number: **2B2527A**

mTesr+ basal medium catalog number: **100-0274**

mTesr+ 10x supplement: **100-0275**

mTesr Plus kit: **100-0276**

mTesr+ LoT number: **AH29845535**

mTESR+ supplemented media

- Add 100ml of mTESR 10x supplement into 400ml of mTESR+ media



- Filter and aliquot 45ml into conical tubes, keep at -20°C.

mFreSR Media:

mFreSR media catalog number: **05855**

mFreSR media LoT number: **NA**

mFreSR media stock concentration: **NA**

mFreSR media working concentration: **NA**

mFreSR media Stock aliquot volume: **50 mL**



Diluting Matrigel

- 1 Thaw Matrigel stock [IM] 9-11 $\mu\text{g}/\mu\text{L}$ stored at $-80\text{ }^{\circ}\text{C}$ on ice
- 2 Make $400\text{ }\mu\text{L}$ aliquots of Matrigel to avoid multiple freeze-thaw cycles using previously cooled tips and tubes
- 3 Store $400\text{ }\mu\text{L}$ Matrigel aliquots at $-20\text{ }^{\circ}\text{C}$ until use
- 4 To prepare matrigel working solution: thaw $400\text{ }\mu\text{L}$ Matrigel aliquots On ice or at $4\text{ }^{\circ}\text{C}$
- 5 Dilute $400\text{ }\mu\text{L}$ Matrigel aliquot in 40 mL of MEM (Working Concentration 90-110 $\mu\text{g}/\text{mL}$)
- 6 Mix well before use. Store at $4\text{ }^{\circ}\text{C}$ protected from the light

Matrigel Coating

1h

- 7 Coat wells with 1 mL of Matrigel per well of 6 well plates or per 35 mm plate
- 8 Leave coated plates for 02:00:00 to 04:00:00 in the incubator at $37\text{ }^{\circ}\text{C}$
- 9 Store diluted matrigel at $4\text{ }^{\circ}\text{C}$

6h

Medium Preparation For Passaging

- 10 Prepare mTesr+ by adding 100 mL of mTesr+ 10x supplement into 400 mL of mTESR+ media and filter it.



- 11 Add  1 mL of [M] 50 Molarity (M) Primocin (final concentration 100ug/mL)
- 12 Take out  3 mL of MEM per well of 6 well plates or per 35 mm plate in a conical tube ( 1 mL /well to wash old medium and  2 mL /well to wash out acutase) and leave at room temperature to warm up
- 13 Take out  1 mL of acutase per well of 6 well plate or per 35mm plate in  15 mL conical tube to warm up at RT
- 14 Take out  3 mL of mTESR+ per well of 6 well plate or 35 mm plate ( 2 mL /well for plating and  1 mL /well for resuspension) to warm up at RT

Washing iPSCs

- 15 Put iPSC plate out of the incubator at  Room temperature
- 16 Aspirate old medium
- 17 Add  1 mL of MEM per well of 6 well plate or per 35 mm plate, rock the plate to wash
- 18 Aspirate Wash

Detaching Adherent Cells

4m

- 19 Add  1 mL of acutase per well of 6 well plate or per 35 mm plate
- 20 Place plate in the incubator at  37 °C for  00:05:00

5m

Centrifuging Cells



21 Get plate out of incubator and check that cells are detached, avoid pippeting too much.

Add  2 mL of MEM

22 Transfer detached cells with MEM ( 3 mL total volume) into 15 mL conical tube

23 Place conical tubes in centrifuge and spin at  1000 rpm for  00:05:00 at  Room temperature

5m

Preparing Seeding Medium

24 Add  1 μ L of CET per  1 mL of mTESR+ (1:1000 ratio)

25 Add  1 μ L of P per  1 mL of mTesr+ (1:1000 ratio)

26 Mix mTESR+ with CETP

27 Aspirate matrigel from coated plates

28 Add  2 mL of mTesr+ containing CETP to each well

Collecting Cells

29 Aspirate supernatant from  15 mL conical tube containing cells (should have pellet)

30 Resuspend pellet with  1 mL of mTESR+ very gently to avoid dissociating iPSCs into single cells

31 Add  50 μ L of cell solution per well in 6 well plate or per 35 mm plate (on top of the mTESR+ with CETP)



32 Gently mix the plate to evenly spread the cells

Incubation

5m

33 Ensure the presence of cells in each well by viewing the plate under the microscope (should see floating cells)

34 Incubate plate in the incubator at  37 °C  Overnight

5m

Changing iPSC medium

35 Aspirate old medium the next day to wash the CETP

36 Add  2 mL of fresh mTESR+ to each well of the 6 well plate or per 35mm plate (no wash with MEM first day after splitting)

Wash and iPSC Medium Change

37 Wash CETP from each well/plate with  1 mL of MEM

38 Aspirate wash

39 Add  2 mL of fresh mTESR+ to each well of the 6 well plate or per 35 mm plate

No need to keep washing plates after the first wash, you can directly change medium

iPSC Maintenance

40 Keep feeding cells everyday until they are ready for splitting again (cells are usually ready for another splitting in 2-3 days, when their confluency reach 70%)

Check cells under microscope to estimate confluency



iPSC Freezing

- 41 **For freezing purposes**, follow the same protocol shown above for cell maintenance for the following sections: Washing iPSCs, Detaching Afdherent cells, centirfuging cells (steps 10-23)
- 42 Instead of collecting cells have freezing vials set and barcoded with iPSC line number, passage number, gene mutation and date
- 43 Resuspend pellet with  1 mL of mFreSR very gently
- 44 Add  500 µL of cell solution to freezing vials

iPSC Banking

- 45 Store freezing vials at  -80 °C in a cryogenic freezing container to prevent ice crystals from forming within the cells in the freezing vials
- 46 After 24 hours, move freezing containers with vials from  -80 °C to a nitrogen tank