

Jan 08, 2025

Version 2

iPSCs Maintenance and Banking V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov1qd5ogr2/v2>

Mahmoud ElAchwah¹, Sol Diaz de Leon Guerrero¹, Zhiping Pang¹

¹Rutgers University



Mahmoud ElAchwah

Rutgers University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov1qd5ogr2/v2>

Protocol Citation: Mahmoud ElAchwah, Sol Diaz de Leon Guerrero, Zhiping Pang 2025. iPSCs Maintenance and Banking. protocols.io <https://dx.doi.org/10.17504/protocols.io.ewov1qd5ogr2/v2> Version created by **Sol Diaz de Leon Guerrero**

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 08, 2025



Last Modified: January 08, 2025

Protocol Integer ID: 117835

Keywords: hiPSC, cell culture, pluripotent stem cell, induced pluripotent stem cell, ipscs maintenance, cell, thorough description of the thawing

Abstract

This protocols offers a thorough description of the thawing, maintenance and banking of human induced pluripotent stem cells.



Materials

Matrigel:

Matrigel catalog number: **Cornig 354234**

Matrigel stock concentration: **9-11 µg/µL**

Matrigel working concentration: **90-110 µg/mL (1:100 dilution)**

Matrigel Stock aliquot volume: **400 µL**

Culture Antibiotic:

Primocin catalog number: **InvivoGen ant-pm-1**

Primocin stock concentration: **50 mg/mL**

Primocin working concentration: **100 µg/mL**

Primocin stock aliquot volume: **1 mL**

CEPT:

Chroman 1 catalog number: **MedChemExpress HY-15392**

Chroman 1 stock concentration: **50 µM**

Chroman 1 working concentration: **50nM**

Emricasan catalog number: **SelleckChem S7775**

Emricasan stock concentration: **5 mM**

Emricasan working concentration: **5 µM**

Polyamine supplement catalog number: **SigmaAldrich P8483**

Polyamine supplement working concentration: **1000x**

Polyamine supplement stock concentration: **1x**

Trans-ISRIB catalog number: **R&D systems 5284**

Trans-ISRIB stock concentration: **700 µM**

Trans-ISRIB working concentration: **700 nM**

Culture Media:

DMEM/F-12 catalog number: **cytiva SH30023.01**



MEM catalog number: **Gibco 51200-038**

Accutase catalog number: **Innovative cell technologies AT-104**

mTESR+ basal medium catalog number: **Stem Cell Technology 100-0274**

mTESR+ 10x supplement: **Stem Cell Technology 100-0275**

mTESR Plus kit: **Stem Cell Technology 100-0276**

mTESR+ supplemented media

- Add 100ml of mTESR 10x supplement into 400ml of mTESR+ media
- Filter and aliquot 45ml into conical tubes, keep at -20°C.

mFreSR media catalog number: **Stem Cell Technology 05855**



Matrigel preparation

- 1 Thaw Matrigel stock stored at $-80\text{ }^{\circ}\text{C}$ by placing on ice and leaving at $4\text{ }^{\circ}\text{C}$
 Overnight
- 2 Place a 10ml pipette, 1 ml pipette tips and 1.5ml eppendorf tubes to cool before preparing Matrigel aliquots
- 3 Place Matrigel stock on ice and mix with the previously cooled 10 ml pipette. Make $400\text{ }\mu\text{L}$ aliquots of Matrigel to avoid multiple freeze-thaw cycles using previously cooled tips and tubes.
- 4 Store $400\text{ }\mu\text{L}$ Matrigel aliquots at $-20\text{ }^{\circ}\text{C}$ until use
- 5 To prepare matrigel working solution thaw $400\text{ }\mu\text{L}$ Matrigel aliquots On ice or at $4\text{ }^{\circ}\text{C}$
- 6 Dilute $400\text{ }\mu\text{L}$ Matrigel aliquot in 40 mL of MEM (1:100 dilution)
- 7 Mix well before use. Store at $4\text{ }^{\circ}\text{C}$ protected from the light

Medium Preparation

- 8 Prepare mTESR+ media by adding 100 mL of mTESR+ 5x supplement into 400 mL of mTESR+ media and filter it.
- 9 Add 1 mL of Primocin (50mg/ml) to the mTESR+ media (final concentration 100ug/mL)

Matrigel Coating

5h



10 Before thawing or passaging cells, coat wells with  1 mL of Matrigel per well of 6 well plates or per 35 mm plate

11 Leave coated plates for  01:00:00 to  04:00:00 in the incubator at  37 °C

5h

12 Store diluted matrigel at  4 °C

iPSC thawing

5m

13 Take out  2 mL of mTESR+ per cell line and leave at room temperature to warm up

14 Take out  5 mL of DMEM/F-12 per cell line and leave at room temperature to warm up

15 Label  15 mL conical tubes per line to thaw

16 Take out iPSC vials in dry ice

17 Thaw iPSC vials using a  37 °C water bath until a small piece of ice is left in the vial

18 Using a 1 ml pipette mix the cells gently and transfer to labeled  15 mL conical tubes.

19 Add 5ml of DMEM/F-12 dropwise to the conical tube with the cells

20 Place conical tubes in centrifuge and spin at  300 x g for  00:05:00 at  Room temperature

5m

21 Supplement mTESR media with CETP

Add  1 µL of CET per  1 mL of mTESR+ (1:1000 ratio)

Add  1 µL of P per  1 mL of mTESR+ (1:1000 ratio)



- 22 Aspirate matrigel from coated plates
- 23 Add  1 mL of mTESR+ containing CETP to each well
- 24 Take out cells from the centrifuge. Aspirate supernatant from  15 mL conical tube containing cells, avoid disturbing the cell pellet
- 25 Resuspend pellet with  1 mL of mTESR+ very gently to avoid dissociating iPSCs into single cells. Transfer all the volume into the well of a 6 well plates or per 35 mm plate, label with the line.
- 26 Gently mix the plate to evenly spread the cells
- 27 Ensure the presence of cells in each well by viewing the plate under the microscope (should see floating cells)
- 28 Incubate plate in the incubator at  37 °C

iPSC passaging

10m

- 29 Take out  3 mL of MEM per well of 6 well plates or per 35 mm plate in a conical tube ( 1 mL /well to wash old medium and  2 mL /well to wash out accutase) and leave at room temperature to warm up
- 30 Take out  1 mL of accutase per well of 6 well plate or per 35mm plate in  15 mL conical tube to warm up at RT
- 31 Take out  3 mL of mTESR+ per well of 6 well plate or 35 mm plate ( 2 mL /well for plating and  1 mL /well for resuspension) to warm up at RT
- 32 Aspirate old medium



- 33 Add  1 mL of MEM per well of 6 well plate or per 35 mm plate, rock the plate to wash
- 34 Aspirate MEM
- 35 Add  1 mL of acutase per well of 6 well plate or per 35 mm plate
- 36 Place plate in the incubator at  37 °C for  00:05:00 5m
- 37 Get plate out of incubator and tap plate gently to check that cells are detached. Pipette cells once before removing from plate, avoid pipetting too much.
- 38 Transfer detached cells into a labeled 15 mL conical tube with  2 mL of MEM
- 39 Place conical tubes in centrifuge and spin at  300 x g for  00:05:00 at  Room temperature 5m
- 40 Supplement mTESR plating media with CETP
Add  1 µL of CET per  1 mL of mTESR+ (1:1000 ratio)
Add  1 µL of P per  1 mL of mTesr+ (1:1000 ratio)
- 41 Aspirate matrigel from coated plates
- 42 Add  2 mL of mTESR+ containing CETP to each well
- 43 Take out cells from the centrifuge. Aspirate supernatant from  15 mL conical tube containing cells, avoid disturbing the cell pellet
- 44 Resuspend pellet with  1 mL of mTESR+ very gently to avoid dissociating iPSCs into single cells



- 45 Add  70-100 μL of cell solution per well in 6 well plate or per 35 mm plate (on top of the mTESR+ with CETP)
- 46 Gently mix the plate to evenly spread the cells
- 47 Ensure the presence of cells in each well by viewing the plate under the microscope (should see floating cells)
- 48 Incubate plate in the incubator at  37 °C

iPSC Maintenance

- 49 After thawing or passaging cells aspirate old medium the next day to wash away the CETP
- 50 Add  2 mL of fresh mTESR+ to each well of the 6 well plate or per 35mm plate
- 51 Cell can be fed every other day (or daily depending on cell confluency and growth) until they are ready for splitting again when cells are 70-80% confluent. Cells are usually ready for splitting every 4 days.

Check cells under microscope to estimate confluency

iPSC Freezing

- 52 Have freezing vials set and labeled with the correct info (iPSC line number, passage number, gene mutation, date, etc).
- 53 Take out  3 mL of MEM per well of 6 well plates or per 35 mm plate in a conical tube ( 1 mL /well to wash old medium and  2 mL /well to wash out accutase) and leave at room temperature to warm up
- 54 Take out  1 mL of accutase per well of 6 well plate or per 35mm plate in  15 mL conical tube to warm up at RT
- 55 Follow the same protocol shown above for cell passaging for the steps 32-39.



- 56 Take out cells from the centrifuge. Aspirate supernatant from  15 mL conical tube containing cells, avoid disturbing the cell pellet
- 57 Resuspend pellet with  1 mL of mFreSR very gently
- 58 Add  500 μ L of cell solution per freezing vial.
- 59 Place freezing vials in cryogenic freezing container to prevent ice crystals from forming within the cells and store at  -80 °C  Overnight
- 60 Transfer vials to a liquid nitrogen tank for storage

Protocol references

Chen, Y., Tristan, C.A., Chen, L., Jovanovic, V.M., Malley, C., Chu, P.H., Ryu, S., Deng, T., Ormanoglu, P., Tao, D., et al. (2021). A versatile polypharmacology platform promotes cytoprotection and viability of human pluripotent and differentiated cells (Springer US).