



Oct 01, 2025

Irradiation and single-cell dissociation of hESCs and cortical spheroids

DOI

<https://dx.doi.org/10.17504/protocols.io.bp2l6xwbzlqe/v1>

Annika Martin¹, Hanqin Li^{1,2}, Dirk Hockemeyer^{1,2}

¹University of California, Berkeley, Berkeley, CA 94720, USA;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 2081, USA



Hanqin Li

UC Berkeley

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l6xwbzlqe/v1>

Protocol Citation: Annika Martin, Hanqin Li, Dirk Hockemeyer 2025. Irradiation and single-cell dissociation of hESCs and cortical spheroids. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bp2l6xwbzlqe/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 11, 2024

Last Modified: October 01, 2025

Protocol Integer ID: 93400

Keywords: ASAPCRN, cortical spheroids this protocol, cell dissociation of hesc, cortical spheroid, cortical spheroids for multi, procedure of irradiation, cell dissociation, irradiation, sequencing protocol, hesc, sequencing

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-000486

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-024409

Abstract

This protocol describes the procedure of irradiation and single-cell dissociation of hESCs and cortical spheroids for MULTI-Seq barcoding and sequencing.

Protocol Overview

A. Irradiation of hESC

B. hESCs derived cortical spheroids

Note

A list of reagents and relevant vendor information can be found in the table listed under the materials tab.

Attachments



[Irradiation of hESCs...](#)

62KB

Materials

Reagents

	Item	Vendor	Catalog Number
	10x HBSS (Ca and Mg Free)	Invitrogen	14185-052
	Sodium Pyruvate 100mM	Life Tech	11360070
	D-Glucose	Sigma	G8769-100ml
	HEPES pH 7.3	Invitrogen	15630-080
	Y-27632 – ROCK Inhibitor	Chemdea	CD0141
	Accutase	Thermo Fisher Scientific	SCR005
	Ethylenediaminetetraacetic acid disodium salt dihydrate (EDTA)	Sigma	E6635-100G
	2-Mercaptoethanol	Sigma	M3148
	L-Cystine solution	Sigma-Aldrich	C7352-100G
	Papain suspension	Worthington Biochemical	LS003126
	PBS- (without Ca and Mg)	Corning	MT21031CV
	PBS+ (with Ca and Mg)	Thermo Fisher	14040-182
	Trypsin Inhibitor	Sigma	T9253-5G
	Costar [®] 6-well Clear Flat Bottom Ultra-Low Attachment 6-well plate	Corning	3471



Irradiation of hESCs

10m

- 1 Culture hESCs with feeder-free system (dx.doi.org/10.17504/protocols.io.b4mcqu2w)
- 2 3 days before irradiation, passage cells and let them grow to approximately 50% confluence.
- 3 On the day of irradiation, change media to hESC media + 10uM Rock Inhibitor (RI).
- 4 Using a discrete cesium source, irradiate cells in plates at desired dosage. We used 0, 0.5, 2, 5 and 10 Gy.
- 5 24 hours post irradiation, dissociate cells by aspirating media and adding  1 mL of accutase per well of a 6-well plate. Return to incubator for  00:05:00 to 10 min. 5m
- 6 Resuspend cells using  10 mL of PBS- to dilute accutase. Spin down for  00:05:00 at  300 x g . 5m
- 7 Aspirate supernatant and resuspend cells in  1 mL PBS- per well.
- 8 Label cells with MULTIseq oligos as described in MULTI-Seq Barcoding and Library Preparation protocol (dx.doi.org/10.17504/protocols.io.kxygx3xzkg8j/v1) and proceed with sequencing.

hESCs Derived Cortical Spheroids

2h

9



Note

Dissociation Media (500mL):

48.9 mL 10x HBSS (Ca and Mg Free)

 440 mL H₂O

5 mL Sodium Pyruvate 100mM

1.1 mL D-Glucose

5 mL HEPES pH 7.3

Differentiate hESCs to Cortical spheroids as described in Cortical Spheroid Differentiation protocol (<https://doi.org/10.17504/protocols.io.5jyl8po57g2w/v1>).

- 10 At 25 days in differentiation, two organoids per condition were transferred with 1.5 mL of media into a 1.5mL Cryovial with screw top. Using a discrete cesium source, irradiate organoids in vials at desired dosage. We used 0, 0.5, 2, 5 and 10 Gy.
- 11 1. Post-irradiation, transfer both organoids per condition into one well of a 6-well low adherence plate with 2.5 mL additional media. Add Rock Inhibitor (RI) to final concentration of 10uM.
- 12 24 hours post irradiation, prepare 5 mL of Activated Papain Solution per irradiation condition:
 - 12.1 To each 5mL of Dissociation Media, add 12 µL of 0.5M EDTA, 23ul of 0.1uM β-mercaptoethanol, 5 µL of 5.5mM L-Cystine solution, and 172 µL of papain suspension.
 - 12.2 Transfer the solution to the 37 °C water bath for ~ 00:20:00 to 30 minutes to activate. The solution will gradually go from cloudy to transparent. 20m
 - 12.3 Sterilize by passing the solution through a 0.22 µm syringe filter.
- 13 Collect cortical spheroids and sediment. Aspirate the media and rinse once with PBS+.
- 14 Add 5 mL Activated Papain Solution to each well and return to the incubator for ~ 00:45:00 . 45m



- 15 While dissociating, prepare 15mL Trypsin Inhibitor solution per  5 mL of Activated Papain Solution.
- 15.1 For each 15mL Trypsin Inhibitor Solution: Add  15 mg Trypsin Inhibitor to  15 mL of usual organoid media. Mix thoroughly and sterilize by filtering through a 0.22 μ m syringe filter.
- 16 After  00:45:00 of dissociation, add  5 mL of Trypsin Inhibitor solution to each well and very gently triturate with a 5mL strippette to dissociate. Avoid excessive pipetting and bubbles. 45m
- 17 Strain mixture through a 70 μ m cell strainer into a 50mL conical tube.
- 18 Use remaining  10 mL Trypsin Inhibitor Solution to rinse the well to collect any cells left behind, then pass this through the same 70 μ m strainer to release any cells still stuck in the mesh.
- 19 Spin down cells at  300 x g for  00:05:00 and aspirate the media. 5m
- 20 Resuspend dissociated cells in  4 mL 0.2% BSA in PBS+ with RI per condition, strain through a 40 μ m mesh cell strainer and FACS sort for desired number of live cells and to remove debris.
- 21 Spin down live cells at  300 x g for  00:05:00 . Resuspend in 0.2% BSA in PBS+ with 10 μ M RI in appropriate volume according to single-cell sequencing protocol, count cells, and proceed to single-cell sequencing according to manufacturer protocol. 5m