

Feb 03, 2023

Isolation and Amplification of SARS-CoV-2 RNA from Nasopharyngeal Specimen

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g7y32kgwz/v1>



Frank Twum Aboagye^{1,2}, Maame Ekua Acquah³

¹Biomedical and Public Health Research Unit, CSIR - Water Research Institute, Accra, Ghana;

²Department of Medical Laboratory Technology, Faculty of Applied Sciences, Accra Technical University, Accra, Ghana.;

³West African Centre for Cell Biology of Infectious Pathogens, College of Basic and Applied Sciences, University of Ghana, Legon, Accra, Ghana



Frank Twum Aboagye

Biomedical and Public Health Research Unit, CSIR - Water Res...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7y32kgwz/v1>

Protocol Citation: Frank Twum Aboagye, Maame Ekua Acquah 2023. Isolation and Amplification of SARS-CoV-2 RNA from Nasopharyngeal Specimen. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.q26g7y32kgwz/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 27, 2023

Last Modified: September 09, 2023

Protocol Integer ID: 75979

Keywords: RNA, RNA Extraction, SARS-CoV-2, Zymo, COVID-19, RT-PCR, Virus, Nasopharyngeal specimen, Viral Transport Medium, rna viral kit, nasopharyngeal specimen, amplification of sar, fever, infection, allplex 2019ncov assay kit, rna, global pandemic, global pandemic by the world health organization, sar

Disclaimer

This is an optimised protocol for SARS-CoV-2 RNA isolation and amplification using Zymo Quick-RNA Viral Kit (200 prep) and Allplex 2019nCoV Assay Kit. The authors do not accept any liability for the collection and handling of both samples and reagents, results from the use of the protocol and its interpretation as well as any errors or omissions that may be made. The reader should make his/her own evaluation as to the appropriateness of the procedures described.

Abstract

COVID-19 caused by the SARS-CoV-2 was declared a global pandemic by the World Health Organization in March 2020. Classical symptoms associated with the infection include fever, cough, chills, headache, and muscle aches amongst others. To effectively diagnose the infection and contain its spread, efficient diagnostic tools are required. The current gold standard for the confirmation of SARS-CoV-2 infection is RT-PCR, this protocol outlines procedures for the isolation and amplification of SARS-CoV-2 RNA from Nasopharyngeal specimens using the Zymo Quick-RNA Viral Kit and Allplex 2019nCoV Assay Kit respectively.

Guidelines

Viral RNA Wash Buffer Preparation

Viral RNA Buffer Preparation

RNA Isolation Procedure (Zymo Quick-RNA Viral Kit)

RT-PCR Mix Preparation Procedure (Allplex 2019nCoV Assay Kit)

RT-PCR Reaction Mix Preparation Procedure

RT-PCR Amplification Protocol (Allplex 2019nCoV Assay Kit)



Materials

Reagent Preparation and RNA Extraction

Dithiothreitol (DTT)

Molecular grade ethanol (95-100%)

DNase/RNase-free water

Zymo Quick-RNA Viral Kit (Zymo Research, USA)

Cryo box

Automatic micropipette (10 μ L, 20 μ L, 200 μ L and 1000 μ L)

RT-PCR

Allplex 2019nCoV Assay Kit (Seegen Inc., Korea)

1.5 mL - 2.0 mL tube ice rack

0.2 mL PCR tube ice rack

Consumables

0.5 - 10 μ L filtered micropipette tips

2 - 20 μ L filtered micropipette tips

20 - 200 μ L filtered micropipette tips

100 - 1000 μ L filtered micropipette tips

1.5 mL microcentrifuge tubes

Viral Transport Medium (Shanghai Focusgen Biotechnology Co., Ltd, China)

96 well qPCR plate and plate cover

Equipment

Centrifuge

Vortex

Biosafety Cabinet II

Bio-Rad CFX 1000 series

Safety warnings

1. Perform all RNA isolation procedures in a Level II Biosafety Cabinet
2. Wear the appropriate PPEs before, during, and after the isolation and amplification of the SARS-CoV-2 RNA while still in the laboratory.
3. Dithiothreitol causes skin and eye irritation. Handle with care
4. Handle all reagents and specimens as potentially hazardous materials



Before start

1. DTT is not stable in solution. Only freshly-made DTT solutions should be used
2. Aliquot the needed volume of reagents from the stock for the procedure to prevent contamination of the large reagent stock.
3. Allow the amplification kit to thaw completely in a 4°C fridge, before preparing the PCR master mix. Avoid centrifuging the reagents to defrost them.
4. Decontaminate all workspace with 1% bleach followed by 70% alcohol



Reagent Preparation - Viral RNA Wash Buffer

- 1 To each bottle of  48 mL of RNA Wash Buffer (concentrate) add  192 mL of molecular grade ethanol (95-100%).
- 2 Invert the reagent bottle several times and label with the date of preparation

Reagent Preparation - Viral RNA Buffer

- 3 Weigh  0.75 g of dithiothreitol (DTT) into  15 mL of nuclease-free water
- 4 Allow the pellets to dissolve completely.
- 5 Add  12 mL of the freshly prepared DTT solution to the 100 mL Viral RNA Buffer (concentrate). Invert several times to mix completely.

Note

Invert the reagent bottle gently to avoid over-frothing the reagent. Allow the reconstituted reagent to sit at room temperature for 60 seconds before use.

RNA Extraction Protocol (Zymo Quick-RNA Viral Kit)

- 6 Vortex the specimen for  00:00:30 at  3000 rpm and allow it to sit on the bench for 1 minute 30s
- 7 Aliquot  200 µL of the specimen into a sterile 1.5 mL microcentrifuge tube
- 8 Add to the specimen  400 µL of the RNA Buffer and vortex at  3000 rpm, 00:00:30 30s



- 9 Transfer the solution in step 3 into a Zymo-Spin™ IC Column in a clean collection tube and centrifuge at 8000 rpm, 00:02:00 2m
- 10 Discard the flow-through liquid and transfer the Zymo-Spin™ IC Column into a new collection tube.
- 11 Transfer into the Zymo-Spin™ IC Column 500 µL of RNA Wash Buffer (refer to reagent preparation for reconstitution of RNA Wash Buffer).
- 12 Centrifuge the Zymo-Spin™ IC Column in the collection tube at 10000 rpm, 00:00:30 30s and discard the flow-through liquid. Repeat this step.
- 13 Add 500 µL of absolute ethanol to the Zymo-Spin™ IC Column and centrifuge at 10000 rpm, 00:01:00 1m and discard the flow-through liquid.
- 14 Perform a dry-spin step at 13000 rpm, 00:01:00 1m rpm to eliminate any trace of ethanol.
- 15 Transfer the Zymo-Spin™ IC Column into a new sterile 1.5 mL microcentrifuge tube and add to the spin column 25 µL of nuclease-free water.
- 16 Centrifuge the spin column for 1 minute at 8000 rpm. Cover the microcentrifuge tube while inspecting for the elute
- 17 Store the RNA at -20 °C pending downstream analysis.

Preparation of PCR Master Mix (Allplex™ 2019nCoV Assay Kit)

18

	A	B	C
	Reagent	X 1 (µL)	X 10 (µL)
	2019 nCoV MuDT* Oligo Mix (MOM)	5.0	50.0
	RNase Free Water (PCR Grade)	5.0	50.0
	Real-time One-step Buffer	5.0	50.0



A	B	C
Real-time One-step Enzyme	2.0	20.0
Total Volume	17.0	170.0

MuDT is the brand name of Seegene's oligo mixture

NB:

1. Add the One-step Enzyme as the last reagent and pipette up and down several times to wash the reagent out of the tip since it is slightly viscous.
2. Pulse vortex the master mix and centrifuge at 3000 rpm for 30 seconds.

Preparation of RT-PCR Reaction Mix

19

A	B	C
Reagent	X 1 (μL)	X 10 (μL)
2019 nCoV MuDT* Oligo Mix (MOM)	5.0	50.0
RNase Free Water (PCR Grade)	5.0	50.0
Real-time One-step Buffer	5.0	50.0
Real-time One-step Enzyme	2.0	20.0
Viral RNA	8.0	-
Total Reaction Volume	25.0	

MuDT is the brand name of Seegene's oligo mixture

NB:

1. Add 4 μL of the exogenous internal control (IC) to the extracted RNA immediately before pipetting the required volume for the reaction.
2. Centrifuge the qPCR plate or 8-strip qPCR tube for 30 seconds at 1000 rpm

RT-PCR Amplification Protocol (Allplex 2019nCoV Assay Kit)



20 **Cycling Conditions:** Reverse transcription at 50°C for 20 minutes, followed by an initial denaturation at 95°C for 15 minutes and 45 cycles of 94°C for 15 seconds and 58°C for 30 seconds (plate read: data acquisition)

Fluorophores: FAM (E gene), Cal Red 610 (RdRP gene), Quasar 670 (N gene), and HEX (Internal Control).

Cut-off Ct-value: 40