

Apr 10, 2023

Isolation of bacterial DNA with Gentra Puregene kit (modified protocol)

DOI

<https://dx.doi.org/10.17504/protocols.io.5qpvorwo7v4o/v1>

o.pogoutse¹, Megan Frederickson¹

¹University of Toronto



o.pogoutse

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpvorwo7v4o/v1>

Protocol Citation: o.pogoutse, Megan Frederickson 2023. Isolation of bacterial DNA with Gentra Puregene kit (modified protocol). protocols.io <https://dx.doi.org/10.17504/protocols.io.5qpvorwo7v4o/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 02, 2023



Last Modified: April 10, 2023

Protocol Integer ID: 77988

Keywords: isolation of bacterial dna, bacterial dna, qiagen protocol for purification, gentra puregene handbook, positive bacterial culture, genomic dna from gram, gentra puregene kit, genomic dna, using yeast, qiagen protocol, dna, purification, qiagen

Abstract

The Qiagen protocol for purification of genomic DNA from gram-positive bacterial cultures using Yeast/Bact. Kit (**Gentra Puregene Handbook - QIAGEN**) was modified to aid in the isolation of bacterial DNA from heterogeneous, low volume, low OD cultures



Cell lysis

- 1 transfer culture to 1.5ml microfuge tube
- 2 centrifuge for 10min at 12000rpm to pellet cells
- 3 decant or aspirate off supernatant
- 4 add 1 ml sterile PBS to pellet and vortex briefly to resuspend cells
- 5 centrifuge for 10 min at 12000rpm
- 6 decant of supernatant
- 7 repeat steps 4-6 one more time
- 8 resuspend pellet in 50ul PBS buffer
- 9 add 2ul of MetaPolzyme solution (5mg/ml) (Sigma#MAC4L)
- 10 incubate overnight at 35C
- 11 add 310ul lysis solution (Qiagen #158113) to each tube. Pipette up and down gently to resuspend cells
- 12 heat samples to 80°C for 5-10 minutes to complete the lysis
**stop step: lysed cells are stable at room temperature



RNase treatment

- 13 add 1.2ul of 10mg/ml RNase A per tube
- 14 invert tube gently to mix
- 15 incubate at 37°C for 30-60 minutes

Protein precipitation

- 16 cool sample to room temperature
- 17 add 105ul of protein precipitation solution (Qiagen #158123)
- 18 invert tube gently several times to mix
- 19 place on ice for 30 minutes
- 20 centrifuge at 12000rpm for 3 minutes.
(The precipitated protein should form a tight pellet).
- 21 transfer supernatant to a fresh sterile tube (by pipette)

DNA precipitation

- 22 add 1-2ul of glycogen (20mg/ml) (Thermofisher #R0561) to the supernatant and mix it by pipetting up and down
- 23 add 400ul of 100% isopropanol



- 24 invert tube gently 50 times to mix
(optionally: incubate at RT for 30 min to improve DNA yield)
- 25 pellet at 13000rpm for 3 minutes
- 26 carefully decant off isopropanol (the pellet might be loose)
- 27 add 500ul of 70% ethanol (cold), and invert tube a few times to wash the DNA pellet
- 28 centrifuge at 13000rpm for 3 minutes
- 29 pour off ethanol carefully – pellet may be loose
- 30 air dry pellet for 15-30 minutes
- 31 rehydrate DNA in ~35-50ul of 10mM Tris buffer (pH7.5-8.0)
(adjust for large or smaller pellet)