


Dec 05, 2024

Version 1

# Isolation of human primary neutrophils V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.dm6gpzq51lzp/v1>



Nataliia Pavliuchenko<sup>1</sup>, Tomáš dička<sup>1</sup>, Kateřina Paldusová<sup>2</sup>, Marek Cebecauer<sup>2</sup>

<sup>1</sup>Institute of Molecular Genetics, Czech Academy of Sciences, Prague;

<sup>2</sup>J. Heyrovský Institute of Physical Chemistry, Czech Academy of Sciences, Prague



Kateřina Paldusová

J. Heyrovský Institute of Physical Chemistry, Czech Academy ...

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.dm6gpzq51lzp/v1>

**Protocol Citation:** Nataliia Pavliuchenko, Tomáš dička, Kateřina Paldusová, Marek Cebecauer 2024. Isolation of human primary neutrophils. [protocols.io https://dx.doi.org/10.17504/protocols.io.dm6gpzq51lzp/v1](https://dx.doi.org/10.17504/protocols.io.dm6gpzq51lzp/v1)

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** August 16, 2024

**Last Modified:** December 05, 2024

**Protocol Integer ID:** 105737

**Keywords:** blood, isolation, neutrophil, neutrophils, primary, human, isolation of human primary neutrophil, human primary neutrophil, isolation

**Funders Acknowledgements:**

AMULET

Grant ID: CZ.02.01.01/00/22\_008/0004558

GAČR

Grant ID: 23-06676S

## Abstract

This protocol was designed for the isolation of human primary neutrophils from full blood or buffy coats. Its thorough optimisation was focused on minimisation of the non-specific activation and population purity.

## Image Attribution

neutrophil

## Materials

- RPMI-1640 + HEPES (Sigma-Aldrich, #R5886) supplemented with 0.2% heat-inactivated and filtered FBS (Sigma-Aldrich, #F7524)
- colour-free RPMI-1640 (Carl Roth, #9104.1) supplemented with 0.2% heat-inactivated and filtered FBS
- 1x PBS (Sigma-Aldrich, #P2272)
- 6% Dextran (Sigma-Aldrich, #31392), 150mM NaCl (Sigma-Aldrich, #S9888) in Milli-Q water, filtered
- 10x ACK - 1.5M NH<sub>4</sub>Cl (Sigma-Aldrich, #09718), 100mM KHCO<sub>3</sub> (Sigma-Aldrich, #60339), 1mM EDTA disodium salt dihydrate (Sigma-Aldrich, #E5134) in Milli-Q water, pH 7.2, filtered

## Protocol materials

 6% Dextran, 150mM NaCl

 1X PBS

 1X ACK buffer

 RPMI+, 0.2% FBS

 Ficoll

 colour-free RPMI, 0.2% FBS



## Safety warnings

- ! ☐ work carefully not to activate the neutrophils
- ☐ supplement the media with FBS right before use (FBS supplement optional or adjustable)
- ☐ prepare 1X ACK right before use

## Ethics statement

Prior ethics approval should be obtained before performing these experiments. This protocol requires prior approval by the user's Institutional Review Board (IRB) or equivalent ethics committee(s) for the work with human blood samples.

## Before start

- ☐ take all reagents out of the refrigerator, temper to room temperature
- ☐ supplement the media with FBS
- ☐ prepare 1X ACK



## Dextran separation

39m

- 1 transfer  blood sample from EDTA vacuette tube into 50ml Falcon tube 30s  

- 2 add  1 mL of  Room temperature  1X PBS per  1 mL of  blood sample 30s  
 
- 3 add  1 mL of  Room temperature  6% Dextran, 150mM NaCl per  1 mL of  blood sample 30s  
 
- 4 invert the tube 10 times 30s  

- 5 incubate at  Room temperature for  00:30:00 without inverting the tube anymore 30m  

- 6 transfer the upper phase into new 50ml Falcon tube, avoid transfer of any erythrocytes from the lower phase 1m  

- 7 spin at  290 rcf, Room temperature, 00:05:00 , swing-up rotor 5m  


### Equipment

centrifuge

NAME

5702

TYPE

Eppendorf

BRAND

022626001

SKU



8 using Pasteur pipette and 2 1ml-pipette tips, discard the supernatant

1m



## Erythrocyte lysis

11m 30s

9 gently resuspend pellet in  15 mL  Room temperature  1X ACK buffer



9.1 first, resuspend in  5 mL , then, add the remaining  10 mL



10 incubate at  Room temperature for  00:05:00 starting from the addition of ACK

5m



11 add  35 mL  Room temperature  1X PBS to neutralise the lysis

15s



12 invert the tube 10 times

15s



13 spin at  290 rcf, Room temperature, 00:05:00 , swing-up rotor

5m



14 using Pasteur pipette and 2 1ml-pipette tips, discard the supernatant, remove as much red debris as possible

1m



## Ficoll separation

37m

15 gently resuspend pellet in  10 mL  Room temperature  RPMI+, 0.2% FBS

1m



16 carefully load  10 mL  Room temperature  Ficoll at the bottom of the tube, make sure not to create any bubbles that would destroy the interface, load at the minimum speed

2m



17 spin at

30m

1260 rcf, Room temperature, 00:30:00 , swing-up rotor (S4×750), acceleration and brake ramp: 5

### Equipment

Centrifuge

NAME

5910 R


TYPE

Eppendorf

BRAND

5943000548

SKU



Phase layout after centrifugation

18 with a cut 1ml pipette tip or 5ml pipette tip, transfer PBMCs from the interface into new 15ml Falcon tube

2m

19 using Pasteur pipette and 2 1ml-pipette tips, remove the rest of the supernatant

1m



- 20 gently resuspend pellet in  5 mL (adjustable)  Room temperature
-  colour-free RPMI, 0.2% FBS (percentage of FBS adjustable)

1m



- 21 count the cells manually or using automated cell counter

