

Nov 12, 2019

Isolation of mononuclear cells using Septmate

DOI

dx.doi.org/10.17504/protocols.io.88chzsw

Ajit N Shah¹

¹UCSF



Ajit N Shah

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.88chzsw

External link: https://cdn.stemcell.com/media/files/pis/10000000612-PIS_01.pdf?_ga=2.182160867.53291116.1573512387-262712883.1573512387

Protocol Citation: Ajit N Shah 2019. Isolation of mononuclear cells using Septmate . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.88chzsw>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 11, 2019

Last Modified: November 12, 2019

Protocol Integer ID: 29668

Attachments



Sepmate_PBMC.pdf

565KB

Materials

MATERIALS

⊗ SepMate™-50 (IVD) 500 Tubes **STEMCELL Technologies Inc. Catalog #85460**

⊗ Dulbecco's Phosphate Buffered Saline with 2% Fetal Bovine Serum 500 mL **STEMCELL Technologies Inc. Catalog #7905**

⊗ Lymphoprep™ **STEMCELL Technologies Inc. Catalog #07801**

⊗ ACK Lysing Buffer **Thermo Fisher Catalog #A1049201**

- 1 Add density gradient medium to the SepMate™ tube by carefully pipetting it through the central hole of the SepMate™ insert. Refer to step 2 for required volumes. The top of the density gradient medium will be above the insert.


NOTE: Small bubbles may be present in the density gradient medium after pipetting. These bubbles will not affect performance.

2

SEPMATE™ TUBE	INITIAL SAMPLE (mL)	DENSITY GRADIENT MEDIUM (mL)
15	0.5 - 4.0	4.5
15	>4-5	3.5
50	4 - 17	15

Sample with Density Gradient Medium

2.1



Overview of Procedure

- 3 Dilute sample with an equal volume of PBS + 2% FBS. Mix gently. For example, dilute 5 mL of sample with 5 mL of PBS + 2% FBS.

- 4 Keeping the SepMate™ tube vertical, add the diluted sample by pipetting it down the side of the tube. The sample will mix with the density gradient medium above the insert.
NOTE: The sample can be poured down the side of the tube. Take care not to pour the diluted sample directly through the central hole.
- 5 Centrifuge at 1200 x g (see Notes) for 10 minutes at room temperature, with the brake on.
NOTE: For samples older than 24 hours, a centrifugation time of 20 minutes is recommended.
- 6 Pour off the top layer, which contains the enriched MNCs, into a new tube. Do not hold the SepMate™ tube in the inverted position for longer than 2 seconds.
NOTE: Some red blood cells (RBCs) may be present on the surface of the SepMate™ insert after centrifugation. These RBCs will not affect performance.

NOTE: To reduce platelet contamination in the enriched MNCs, pipette off some of the supernatant above the MNC layer before pouring.
- 7 Wash enriched MNCs with PBS + 2% FBS. Centrifuging at 300 xg for 8 minutes at RT.
- 8 Add ACK Lysis buffer and incubate for 3 minutes. Add 3-4 volumes of PBS + 2% FBS. Centrifuging at 400 xg for 5 minutes at RT.
- 9 Wash enriched MNCs with PBS + 2% FBS. Centrifuging at 400 xg for 5 minutes at RT.
- 10 Re-suspend in desired volume of media/buffer and perform cell count.
- 11 PMBC are ready for downstream process.