



Mar 15, 2024

Isolation of NeuN+ cells from brain tissue (for CUT and RUN)

DOI

<https://dx.doi.org/10.17504/protocols.io.4r3l27pejg1y/v1>

anita.adami¹

¹Laboratory of Molecular Neurogenetics, Department of Experimental Medical Science, Wallenberg Neuroscience Center and Lund Stem Cell Center, BMC A11, Lund University, 221 84 Lund, Sweden.



Raquel Garza

Lund University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l27pejg1y/v1>

Protocol Citation: anita.adami 2024. Isolation of NeuN+ cells from brain tissue (for CUT and RUN). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.4r3l27pejg1y/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 03, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 78067

Keywords: ASAPCRN, brain tissue, preparation for cut, cell, tissue

Funders Acknowledgements:

Aligning Science Across Parkinson's through the Michael J. Fox Foundation for Parkinson's Research

Grant ID: ASAP-000520

Swedish Research Council

Grant ID: 2018-02694

Swedish Brain Foundation

Grant ID: FO2019-0098

Cancerfonden

Grant ID: 190326

Barncancerfonden

Grant ID: PR2017-0053

NIHR Cambridge Biomedical Research Centre

Grant ID: NIHR203312

Swedish Society for Medical Research

Grant ID: S19-0100

National Institutes of Health

Grant ID: HG002385

Swedish Research Council

Grant ID: 2021-03494

Swedish Research Council

Grant ID: 2020-01660

Abstract

This protocol describes the steps to isolate NeuN+ cells from brain tissue in preparation for CUT and RUN



Isolation of NeuN+ cells

45m

- 1 Nuclei were isolated from frozen tissue as described [here](#).
- 2 Before FACSing, nuclei were incubated with Recombinant Alexa Fluor® 488 Anti-NeuN antibody [EPR12763] - Neuronal Marker (Abcam, Cat# ab190195, RRID:AB_2716282) at a concentration of 1:500 for ⌚ 00:30:00 🌡️ On ice as according to [Spalding, K.L. et al. \(2013\)](#). 30m
- 3 The nuclei were run through the FACS at 🌡️ 4 °C with a low flow rate using a 100 mm nozzle and 300.000 nuclei Alexa Fluor – 488 positive nuclei were sorted.
- 4 The sorted nuclei were pelleted at 1,300 x g for ⌚ 00:15:00 and resuspended in 1 mL of ice-cold nuclear wash buffer (20 mM HEPES, 150 mM NaCl, 0.5 mM spermidine, 1x cOmplete protease inhibitors, 0.1% BSA) and 10 µL per antibody treatment of ConA-coated magnetic beads (Epiccypher) added with gentle vortexing (Pipette tips for transferring nuclei were pre-coated with 1% BSA), to perform [CUT&RUN](#). 15m