

Oct 28, 2019

Isolation of single nuclei from postmortem fresh frozen human brain and immunostaining for NeuN

DOI

<https://dx.doi.org/10.17504/protocols.io.6tuhenw>

Single nuclei
isolation

Marcos Otero-Garcia¹, Inma Cobos¹

¹Stanford University

Neurodegeneration Method Development Community
Tech. support email: ndcn-help@chanzuckerberg.com



Inma Cobos

Stanford

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6tuhenw>

Protocol Citation: Marcos Otero-Garcia, Inma Cobos 2019. Isolation of single nuclei from postmortem fresh frozen human brain and immunostaining for NeuN. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.6tuhenw>



Manuscript citation:

Manuscript under review

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 26, 2019

Last Modified: October 28, 2019

Protocol Integer ID: 27220

Keywords: human brain tissue, postmortem, immunostaining, NeuN, single neuronal nuclei, isolation, FACS, Nuc-seq, single nuclei from fresh frozen postmortem brain, postmortem fresh frozen human brain, fresh frozen postmortem brain, fresh frozen human brain, fresh frozen mouse brain, isolation of all nuclei, isolation of nuclei, single nuclei, nucleus rna, subjects with neurodegenerative disorder, rna, nuclei, neurodegenerative disorder, seq, lab for the isolation

Abstract

- Protocol based on Krishnaswami *et al.*, Nat Protoc. 2016, 3:499-524
- Used routinely in our lab for the isolation of single nuclei from fresh frozen postmortem brains from subjects with neurodegenerative disorders and healthy controls for single-nucleus RNA-seq
- Works also for isolation of nuclei from fresh frozen mouse brains
- Used for isolation of all nuclei or antibody-enriched populations (i.e., NeuN⁺ nuclei)
- We always perform and strongly recommend FACS if used for single-nucleus RNA-seq

Attachments



[Single nuclei isolat...](#)

383KB

Materials

MATERIALS

- ✕ DNase I (RNase-free) - 1,000 units **New England Biolabs Catalog #M0303S**
- ✕ Bovine Serum Albumin (BSA) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A7906**
- ✕ Magnesium chloride hexahydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M2670**
- ✕ UltraPure™ DNase/RNase-Free Distilled Water **Thermo Fisher Scientific Catalog #10977023**
- ✕ DAPI **Invitrogen - Thermo Fisher Catalog #D3571**
- ✕ Sucrose **Fisher Scientific Catalog #S25590B**
- ✕ Triton™ X-100 **Fisher Scientific Catalog #AC215680010**
- ✕ Potassium chloride solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #60142**
- ✕ UltraPure™ 1M Tris-HCl pH 8.0 **Thermo Fisher Scientific Catalog #15568025**
- ✕ DL-Dithiothreitol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D9779-1G**
- ✕ cOmplete™ Protease Inhibitor Cocktail **Merck MilliporeSigma (Sigma-Aldrich) Catalog #11697498001**
- ✕ NxGen® RNase Inhibitor **Lucigen Catalog #30281-2**
- ✕ Anti-NeuN Antibody **Merck Millipore (EMD Millipore) Catalog #MAB377**
- ✕ Goat anti-Mouse IgG (H L) Alexa Fluor 647 **Thermo Fisher Scientific Catalog #A-21235**
- ✕ OptiPrep™ Density Gradient Medium **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D1556)**

STEP MATERIALS

- Material and tools needed: Forceps, spatula, blades, dounce tissue grinder, cell strainers, petri dishes

Equipment	
Dounce All-Glass Tissue Grinders	NAME
Tissue Grinder	TYPE
KIMBLE	BRAND
885300-0007	SKU
https://www.kimble-chase.com/advancedwebpage.aspx?cg=886&cd=4&SKUTYPE=202&SKUFLD=SKU&DM=1250&WEBID=6855	LIN K



Buffer preparation (All solutions should be RNase-free for single-soma RNAseq experiments):

Isolation Medium #1 (IM1), 45 ml (optional)

Prepare in a 50 ml Falcon tube and store at  4 °C up to 6 months.

	<i>Amo unt</i>	<i>Rea gent</i>	<i>Final conc entr atio n</i>
	750 0 µl	1.5 M Sucr ose	250 mM
	1125 µl	1 M KCl	25 mM
	225 µl	1M MgC l ₂	5 mM
	450 µl	1 M Tris (pH 8)	10 mM
	35.7 ml	RNa se- free H ₂ O	

Homogenization Buffer (3 ml per sample)

Prepare FRESH and keep ICED or at  4 °C . Discard after use.

	<i>Amo unt</i>	<i>Rea gent</i>	<i>Final conc entr atio n</i>
	292 5 µl	IM1	
	3 µl	DTT 1mM	1 µM
	30 µl	50x Prot ease Inhib itor	0.5 x



	15 μ l	RNaseIN 40U/ μ l	0.2 U/ μ l
	30 μ l	TritonX100 10 %v/ v	0.1 %

Iodixanol dilutions

Prepare in a 50 ml Falcon tubes and store at  4 °C up to 6 months. Accuracy with Iodixanol and sucrose concentrations is critical.

■ Iodixanol medium (IDM), 45 ml

	<i>Amount</i>	<i>Reagent</i>	<i>Final concentration</i>
	7500 μ l	1.5 M Sucrose	250 mM
	6750 μ l	1 M KCl	150 mM
	1350 μ l	1 M MgCl ₂	30 mM
	2700 μ l	1 M Tris (pH 8)	60 mM
	26.7 ml	Nuclease-free H ₂ O	

■ Iodixanol 50%(v/v), 20 ml

16.7 ml Iodixanol 60 %(v/v) + 3.3 ml IDM

■ Iodixanol 29%(v/v), 30 ml

14.5 ml Iodixanol 60 %(v/v) + 15.5 ml IDM

**Freezing Storage Buffer (FSB), 15ml** (optional)

Prepare in a 50 ml Falcon tube and store at  4 °C up to 6 months.

	<i>Amount</i>	<i>Reagent</i>	<i>Final concentration</i>
	1.665 ml	1.5M Sucrose	166.5 mM
	75 µl	1M MgCl ₂	5 mM
	150 µl	1M Tris (pH8)	10 mM
	13.1 ml	Nuclease-free H ₂ O	

Buffer for Immunostaining, 10ml

	<i>Amount</i>	<i>Reagent</i>	<i>Final concentration</i>
	9400 µl	RNAse-free PBS (pH7.4)	
	500 µl	RNAse-free 10 % BSA	0.5 %
	50 µl	1 M MgCl ₂	5 mM
	10 µl	DNAse I (2000)	2 U/ml



		U/ml)	
--	--	-----------	--

Buffer for Antibody Incubation

Add to 1 ml of Buffer for Immunostaining:

	<i>Amo unt</i>	<i>Rea gent</i>	<i>Final conc entr atio n</i>
	5 µl	RNA seIN 40 U/µl	0.2 U/µl

Collection medium for FACS (1 ml; 0.2 ml/vial)

	<i>Amo unt</i>	<i>Rea gent</i>	<i>Final conc entr atio n</i>
	950 µl	RNA se- free PBS (pH7 .4)	
	25 µl	RNA se IN	
	*	BSA 10% *Aft er colle cting !	1 %



Protocol materials

- ✕ cComplete™ Protease Inhibitor Cocktail **Merck MilliporeSigma (Sigma-Aldrich) Catalog #11697498001**
- ✕ Anti-NeuN Antibody **Merck Millipore (EMD Millipore) Catalog #MAB377**
- ✕ DL-Dithiothreitol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D9779-1G**
- ✕ DAPI **Invitrogen - Thermo Fisher Catalog #D3571**
- ✕ Bovine Serum Albumin (BSA) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A7906**
- ✕ Goat anti-Mouse IgG (H L) Alexa Fluor 647 **Thermo Fisher Scientific Catalog #A-21235**
- ✕ Magnesium chloride hexahydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M2670**
- ✕ UltraPure™ DNase/RNase-Free Distilled Water **Thermo Fisher Scientific Catalog #10977023**
- ✕ Potassium chloride solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #60142**
- ✕ DNase I (RNase-free) - 1,000 units **New England Biolabs Catalog #M0303S**
- ✕ OptiPrep™ Density Gradient Medium **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D1556)**
- ✕ Sucrose **Fisher Scientific Catalog #S25590B**
- ✕ NxGen® RNase Inhibitor **Lucigen Catalog #30281-2**
- ✕ Triton™ X-100 **Fisher Scientific Catalog #AC215680010**
- ✕ UltraPure™ 1M Tris-HCl pH 8.0 **Thermo Fisher Scientific Catalog #15568025**
- ✕ Anti-NeuN Antibody, clone A60 **Merck Millipore (EMD Millipore) Catalog #MAB377**



Safety warnings

- ⚠ Please see SDS (Safety Data Sheet) for hazards and safety warnings.
Uses fresh human brain tissue - Biosafety Level 2 lab work.

Before start

- All solutions and materials should be **RNase-free** and **kept iced or at**  4 °C at all times
- Glassware and metal tools are **sealed with aluminum foil** and **baked at**  220 °C **for**  06:00:00 .



1 Prepare **Homogenization Buffer** and cool it on ice.

	<i>Amount</i>	<i>Reagent</i>	<i>Final concentration</i>
	2925 μ l	IM1	
	3 μ l	DTT 1mM	1 μ M
	30 μ l	50x Protease Inhibitor	0.5 x
	15 μ l	RNaseIN 40U/ μ l	0.2 U/ μ l
	30 μ l	TritonX100 10 %v/v	0.1 %

Homogenization Buffer (3 ml per sample)

Prepare FRESH and keep ICED or at  4 °C . Discard after use.

1.1 Add  2925 μ L IM1 .

1.2 Add  3 μ L DTT 1mM .

1.3 Add  30 μ L 50x Protease Inhibitor .

1.4 Add  15 μ L RNaseIN 40U/ μ l .

1.5 Add  30 μ L TritonX100 10 %v/v .



- 2 Pre-cool the dounce tissue grinder on ice (**Kimble Kontes all glass tissue grinder**, 7 mL tubes, 0.02 - 0.056 mm clearance space between pestle and tubes)
- 3 Add  2.4 mL of **Homogenization buffer** to the dounce tissue grinder.
- 4 Collect the brain chunk (~  100 mg) and transfer to a Petri dish on ice.
- 5 Cut out into small pieces using a chilled scalpel or blade.
- 6 Transfer all pieces into the dounce tissue grinder.
- 7 Homogenize the tissue, on ice (**~30 firm strokes**).
Check on **hematocytometer** while homogenizing and adjust the number of strokes.

Note

Too few strokes – you will see cell clumps; too many strokes – you will get damaged nuclei.

- 8 Filter homogenate using **40 µm** Corning cell strainer to remove clumps.
→ Take sample for hemocytometer.
- 9 Transfer the homogenate into two precooled 1.5 ml Eppendorf tubes.
- 10 Centrifuge at  1000 x g for  00:08:00 at  4 °C .
- 11 Slowly aspirate the supernatant from the side of the tube.
Avoid disturbing the pellets, up to  50 µL of supernatant can be left.
- 12 Gently **resuspend** each pellet in  225 µL (final volume) of cold **Homogenization buffer** and pool both tubes (final volume  450 µL).



13 Add an **equal** volume of cold [M] 50 %v/v iodixanol .



Note

Critical! Be very exact with volumes. 450 µL suspension + 450 µL 50 % iod. Final iodixanol concentration is **25 %**.

14 Gently pipette mix.



15 Add an equal volume of **29 % iodixanol** (900 µL) into a 2 ml Eppendorf tube.

16 Slowly layer off the 25 % iodixanol suspension mix over the 25 % iodixanol, without mixing them.

17 Centrifuge at 13500 x g for 00:20:00 at 4 °C .



18 Prepare **buffer for immunostaining**.



	Amount	Reagent	Final concentration
	9400 µl	RNAse-free PBS (pH7.4)	
	500 µl	RNAse-free 10 % BSA	0.5 %
	50 µl	1 M MgCl ₂	5 mM
	10 µl	DNAse I (2000 U/ml)	2 U/ml



		U/ml)	
--	--	-----------	--

18.1 Add  9400 μ L RNase-free PBS (pH7.4) .

18.2 Add  500 μ L RNase-free 10 % BSA .

18.3 Add  50 μ L 1 M MgCl₂ .

18.4 Add  10 μ L DNase I (2000 U/ml) .

19 Prepare **buffer for Antibody incubation** (1.5 ml per sample) by adding  5 μ L of RNaseIN 40 U/ μ l to  1 mL of **buffer for Immunostaining** (final concentration 0.2U/ μ l)

20 After centrifugation, remove and discard the top myelin-rich debris layer (you can use a 1 ml pipette with the tip cut, or cotton swabs).

21 Then, remove and discard the aqueous supernatant, without disrupting the **nuclei pellet**. Avoid contaminating with the top layer.



22 Use a small amount of **buffer for Antibody Incubation** to resuspend the pellet and transfer the solution to a new tube.

23 Gently resuspend in  500 μ L of **buffer for Antibody Incubation**.

24 Incubate for  00:15:00 , at  4 $^{\circ}$ C or iced, for **blocking** nonspecific staining.



25 Take sample for hemocytometer.

26 Take a sample for **unstained** control. Take sample for **2AB-only** (nonspecific binding control).



- 27 Add primary antibody (**Ms-a-NeuN, 1:1000**) and incubate on a rotator in cold room ( 4 °C) for  00:40:00 .



Anti-NeuN Antibody, clone A60 **Merck Millipore (EMD Millipore) Catalog #MAB377**

Note

Use Eppendorfs coated with BSA for collection (to coat the Eppendorf tubes, fill them with 10 % BSA solution in PBS for  00:05:00 , rinse with PBS, and dry at  4 °C overnight).

Washing

- 28 Add  700 µL of **buffer for Immunostaining** and invert several times. 
- 29 Centrifuge at  500 x g for  00:05:00 at  4 °C . 
- 30 Carefully transfer all the supernatant (don't leave any supernatant! It's ok to drag some pellet) to a new Eppendorf tube and resuspend the pellet in  400 µL of **buffer for Antibody Incubation**. 
- 31 Centrifuge the supernatant at  500 x g for  00:05:00 at  4 °C to recover non-pelleted nuclei. 
- 32 Resuspend them in  200 µL of **buffer for Antibody Incubation** and pool for a final volume of 600 µL. 
- 33 Take samples for single color staining.
- 34 Add secondary antibody (**G-a-Ms Alexa Fluor 647, 1:1000**) and **Hoechst 33258** (5 µg/ml) to the same tube. (**Hoechst 1:1000 of 2,5 mg/ml stock**).
- 35 Check the quality of the sample and on hemocytometer.

**Note**

For 100 mg of tissue, should have ~ 2 Million nuclei in ~0.6 ml.

36 Transfer to 7 ml culture tubes.

37 Place on ice and bring them to the FACS facility for sorting.

Note

FACS sorting controls:

Non-stained

2AB only

Single stainings

38 Collect in BSA coated Eppendorfs containing **collection medium** (containing ~1/5 of the expected final volume after collection).

39 Add BSA after collection for final 1% BSA final concentration (tested for 10x Chromium v2 and v3).

Note

Different Single Cell assays tolerate different BSA concentrations, but lowering it may increase nuclei aggregation.