


Oct 25, 2019

🌐 Isolation of single nuclei from solid tissues

DOI

dx.doi.org/10.17504/protocols.io.ufketkw



Blue Lake¹, Kun Zhang¹

¹University of California, San Diego

Human Cell Atlas Metho...

KPMP

1 more workspace



Blue Lake

UCSD

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ufketkw>

External link: <http://genome-tech.ucsd.edu/ZhangLab/>

Protocol Citation: Blue Lake, Kun Zhang 2019. Isolation of single nuclei from solid tissues. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.ufketkw>

**Manuscript citation:**

Lake, B.B. et al. A single-nucleus RNA-sequencing pipeline to decipher the molecular anatomy and pathophysiology of human kidneys. *Nature Communications* 10, 2832 (2019).

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 09, 2018

Last Modified: October 25, 2019

Protocol Integer ID: 16588

Keywords: single nuclei, frozen tissues, genomic assays, isolation of single nuclei, solid tissues nuclei, isolation of nuclei, nucleus genomic assay, whole cell dissociation method, single cell, incomplete dissociation of solid tissue, adult human tissue atlas, human tissue atlas, nuclei, frozen tissue, discovery of molecular cell type, cell, solid tissue, rna, rna degradation, molecular cell type, tissue, overall cellular makeup, isolation

Abstract

Nuclei can be readily isolated from frozen tissues with a combination of chemical and physical treatments that can circumvent the non-uniform or incomplete dissociation of solid tissues into single cells. The isolation of nuclei can also circumvent RNA degradation or any introduction of technical artefacts (such as stress responses) that could be triggered during whole cell dissociation methods. Data generated from single-nucleus genomic assays permits discovery of molecular cell types that can be used to define the overall cellular makeup of a tissue or organ, and ultimately will inform upon adult human tissue atlases.

Materials

MATERIALS

⊗ DAPI **Thermo Fisher Scientific Catalog #D1306**

⊗ Dounce homogenizers **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D8938-1SET**

⊗ RNase Inhibitor **Enzymatics Catalog #Y9240L**

⊗ CellTrics Filters (30um) **Sysmex Catalog #04-004-2326**

STEP MATERIALS

⊗ RNase Zap **Merck MilliporeSigma (Sigma-Aldrich) Catalog #R2020-250ML**

⊗ RNAlater **Thermo Fisher Scientific Catalog #AM7020**

⊗ RNase Inhibitor **Enzymatics Catalog #Y9240L**

⊗ cOmplete™, Mini Protease Inhibitor Cocktail **Roche Catalog #11836153001**

⊗ DAPI **Invitrogen - Thermo Fisher Catalog #D3571**

⊗ RNAlater **Thermo Fisher Scientific Catalog #AM7020**



Protocol materials

⊗ RNase Inhibitor **Enzymatics Catalog #Y9240L**

⊗ CellTrics Filters (30um) **Sysmex Catalog #04-004-2326**

⊗ RNAlater **Thermo Fisher Scientific Catalog #AM7020**

⊗ RNase Inhibitor **Enzymatics Catalog #Y9240L**

⊗ cOmplete™, Mini Protease Inhibitor Cocktail **Roche Catalog #11836153001**

⊗ DAPI **Invitrogen - Thermo Fisher Catalog #D3571**

⊗ DAPI **Thermo Fisher Scientific Catalog #D1306**

⊗ Dounce homogenizers **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D8938-1SET**

⊗ RNase Zap **Merck MilliporeSigma (Sigma-Aldrich) Catalog #R2020-250ML**

⊗ RNAlater **Thermo Fisher Scientific Catalog #AM7020**

⊗ DAPI **Invitrogen - Thermo Fisher Catalog #D3571**

⊗ RNase Inhibitor **Enzymatics Catalog #Y9240L**

⊗ cOmplete™, Mini Protease Inhibitor Cocktail **Roche Catalog #11836153001**

⊗ RNase Zap **Merck MilliporeSigma (Sigma-Aldrich) Catalog #R2020-250ML**

⊗ RNAlater **Thermo Fisher Scientific Catalog #AM7020**

⊗ RNAlater **Thermo Fisher Scientific Catalog #AM7020**

Troubleshooting



Prepare Reagents and Tissue

- 1 Prepare NEB-complete (NEB containing 5 ug/ml DAPI and 0.04 U/ul RNase Inhibitor)

chill on ice

Final Concentration	Stock	Volume (25 ml)
20 mM Tris [pH 8]	1M	0.5ml
320 mM sucrose	1M	8ml
5 mM CaCl ₂	1M	125μl
3 mM MgAc ₂	0.5 M	75μl
0.1 mM EDTA	10% -	5μl
0.1% TritonX-100		250 μl
dH ₂ O		16ml


NEB Base Solution Composition


Note

For chromatin accessibility assays, include 1:100 dilution of cOmplete™ Protease Inhibitor Cocktail (stock one tablet in 0.5 ml H₂O)

 DAPI Invitrogen - Thermo Fisher Catalog #D3571

 RNase Inhibitor Enzymatics Catalog #Y9240L

 cOmplete™, Mini Protease Inhibitor Cocktail Roche Catalog #11836153001

- 2 Treat dounce with RNaseZap, rinse with sterile water (if possible: UV treat  00:15:00)

 RNase Zap Merck MilliporeSigma (Sigma-Aldrich) Catalog #R2020-250ML

- 3 Transfer vial containing tissue to ice.

Note

For solid tissues (e.g. adult human kidney), 40 µm cryosections can be used with the number of sections dependent on desired yield and the size and type of tissue. For kidney ~6 cubic mm will give ~150-200K nuclei.


- 4 For sections stored in a stabilizing solution (e.g. RNAlater), wash briefly with PBS and immediately proceed to Step 5 below

 RNAlater Thermo Fisher Scientific Catalog #AM7020

	Final Conc	Stock	Volume (50 ml)
	1x PBS	10x	5ml
	1 mM EGT A	0.1M	50µl
	dH ₂ O		45ml

PBSE Composition

Isolate Nuclei

- 5 Add  1 mL ice cold NEB buffer to tissue segments
- 6 Cut end off of a p1000 tip to increase bore size using a sterile scalpel, then pipette sections up and down to disperse and dissolve OCT, ~20 times
- 7 Using a regular p1000 tip, pipette ~10x to further dissociate tissues into manageable sizes.

Note

tissue needs to be passable through a p1000 tip easily before proceeding

Then transfer to dounce homogenizer

- 8 Gently dounce tissue on ice:

5 strokes with pestle A





~20 strokes with pestle B (minimize bubble formation)

**Note**

Increase number of pestle A strokes if the tissue appears too granular before proceeding with pestle B. Number of pestle B strokes used here is dependent on tissue toughness:
soft tissues use ~10-15 strokes
hard tissues use ~15-20 strokes

Note

Avoid making bubbles

- 9 Transfer solution to a  15 mL tube
- 10 Wash dounce with  1 mL NEB-complete buffer and add this into the same tube
- 11 Incubate on ice  00:10:00
- 12 Pass supernatant through 30 uM CellTrics filter to a new  15 mL conical tube

Equipment

new equipment

NAME

Sysmex

BRAND

04-004-2324

SKU


30 uM Celltrics Filter

SPECIFICATIONS





13 Bring up to  10 mL with PBSE

14 Pellet nuclei:  900 g

 00:10:00 at  4 °C

snRNA-Seq methods: nuclei can be stored in RNAlater

15 Remove supernatant and resuspend pellet in  100 µL -  1000 µL PBS + 0.1% RNAse Inhibitor

Note

Resuspension buffer and volume is dependent on downstream assays and nuclei concentration requirements. 1% BSA can be included here

16 **QA/QC:** Count nuclei (e.g. BioRad T20 Cell Counter)



Equipment

new equipment

NAME

Bio-Rad

BRAND

1450011

SKU

Cell Counting Slides for TC10™/TC20™ Cell Counter, Dual-Chamber

SPECIFICATIONS

- 17 **QA/QC:** Check nuclei integrity under fluorescent microscope using DAPI channel. Nuclei should appear distinct, have rounded borders and the majority occurring as singlets.

Note

High clumping rates would indicate damaged nuclei and would require re-filtering using 30-µm CellTrics filter or exclusion from downstream analyses.

At least **50,000 nuclei** are needed to proceed with snDrop-seq






At least **10,000 nuclei** are needed to proceed with 10X 3' RNA v3

Isolate Nuclei


- 18 To use nuclei directly for single nucleus assays, proceed to method

To use nuclei on a later date, proceed to Step 19

snRNA-Seq methods: nuclei can be stored in RNAlater



- 19 Add  900 µL RNAlater to  100 µL nuclei in PBS, incubate at  4 °C for  01:00:00 to  02:00:00



then transfer to  -20 °C for 1-2 months



RNAlater Thermo Fisher Scientific Catalog #AM7020

20 To remove RNAlater, centrifuge nuclei at 4000g,  00:10:00 at  4 °C .

Remove solution and resuspend in associated nuclei resuspension buffer (assay dependent)