


Dec 20, 2023

Isotopically labelled inorganic carbon delivered to algal cultures via bubbler bottle

DOI

dx.doi.org/10.17504/protocols.io.yxmvm3q8ol3p/v1



Usha F Lingappa¹, Sunnyjoy Dupuis¹, Xavier Mayali², Sabeeha S. Merchant¹

¹University of California, Berkeley; ²Lawrence Livermore National Laboratory

Merchant Lab UC Berkel...



Usha F Lingappa

University of California Berkeley

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvm3q8ol3p/v1>

Protocol Citation: Usha F Lingappa, Sunnyjoy Dupuis, Xavier Mayali, Sabeeha S. Merchant 2023. Isotopically labelled inorganic carbon delivered to algal cultures via bubbler bottle. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.yxmvm3q8ol3p/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 28, 2023

Last Modified: December 20, 2023

Protocol Integer ID: 91536

Keywords: stable isotope probing, algae, bicarbonate, carbonate chemistry, bubbling, chlamydomonas reinhardtii, photosynthesis, algal cultures via bubbler bottle, carbonate chemistry, solution as co2 gas, co2 gas, bubbling air, labelled inorganic carbon, atmospheric co2, algal culture, bubbler bottle, equilibrium with atmospheric co2, inorganic carbon, labelled gas, gas with air, bubbling preclude, label to culture

Funders Acknowledgements:

Sabeeha Merchant, Moore Foundation Symbiosis in Aquatic Systems Initiative Investigator Award

<https://doi.org/10.37807/GBMF9203>

Grant ID: GBMF9203

Xavier Mayali, US Department of Energy contract DE-AC52-07NA27344

Grant ID: SCW1039

Abstract

This protocol describes a method for delivering labelled inorganic carbon as $^{13}\text{CO}_2$ to algal cultures by bubbling air through a solution of $\text{H}^{13}\text{CO}_3^-$, and then into the culture. We developed this method to deliver label to cultures grown under continuous bubbling with air, without the use of $^{13}\text{CO}_2$ labelled gas and the necessary equipment to mix labelled gas with air at near-atmospheric levels. Bubbling precludes the more common approach of adding $\text{H}^{13}\text{CO}_3^-$ label directly to the media, because dissolved HCO_3^- is in equilibrium with atmospheric CO_2 . Thus, excess HCO_3^- added to a solution will leave the solution as CO_2 gas as it equilibrates. Bubbling rapidly accelerates this equilibration which is typically diffusion limited. The method described here takes advantage of this aspect of carbonate chemistry, and uses a solution of $\text{H}^{13}\text{CO}_3^-$ —which is less expensive and more convenient than $^{13}\text{CO}_2$ gas—to generate a flux of $^{13}\text{CO}_2$ that can be bubbled into a culture.



Materials

Aquarium pump

Bottle with vent and fill port assembly lid

Tubing and luer locks

Bubbling flask assembly (Erlenmeyer flask, foam plug, serological pipet, syringe filter)

Culture medium & inoculum

Na₂CO₃

H₂KPO₄

HK₂PO₄

Troubleshooting

The purpose and problem of bubbling cultures

- 1 Bubbling air into algal cultures stimulates photosynthetic growth by ameliorating diffusion limitation for CO_2 (fig. 1A-B). Stable isotope probing (SIP) experiments examining carbon fixation often involve ^{13}C label introduced directly into the culture medium as a $\text{H}^{13}\text{CO}_3^-$ salt. Unfortunately, because dissolved inorganic carbon (DIC) is in equilibrium with atmospheric CO_2 (eq. 1), this labelling method does not work for cultures that are bubbled in an open system. Bubbling accelerates equilibration, causing excess HCO_3^- to rapidly exit the solution so that the label is lost before it can be fixed into biomass (fig. 1C). This presents some inconvenience, as $\text{H}^{13}\text{CO}_3^-$ salts are cheaper and easier to work with than is $^{13}\text{CO}_2$ gas.





Figure 1. Bubbling enhances photosynthetic growth but precludes stable isotope probing with labelled bicarbonate. **A.** Diagrams of un-bubbled vs. bubbled culture formats. Culture data presented in this protocol are of the unicellular green alga *Chlamydomonas reinhardtii*, grown in minimal medium under continuous light at 28°C and shaken at 125 RPM. **B.** Growth curves of *C. reinhardtii* in un-bubbled vs. bubbled culture formats, with and without added bicarbonate. Bubbling greatly enhances autotrophic growth. At the concentrations shown, bicarbonate addition does not impact growth. **C.** ^{13}C enrichment of *C. reinhardtii* biomass grown with ^{13}C labelled bicarbonate added directly to the culture medium for 24 hours, measured by IRMS on lyophilized cell pellets. Unbubbled cultures exhibit substantial ^{13}C enrichment; ~90% of that signal is lost when bubbled.



Label delivery via bubbler bottle

- 2 We developed a SIP method that takes advantage of this liability of bubbling and carbonate chemistry (figure 2). Using an aquarium pump, we bubble air into a solution of $\text{H}^{13}\text{CO}_3^-$. The DIC in the solution exchanges with the air, releasing label in the form of $^{13}\text{CO}_2$ gas. That air is then bubbled into algal cultures. Using this method, we achieved substantial biomass ^{13}C enrichment in bubbling *C. reinhardtii* cultures (fig. 2B).




Figure 2. Bubbling air through a solution of $\text{H}^{13}\text{CO}_3^-$ is an effective strategy to label bubbling cultures. **A.** Diagram of our bubbler label delivery approach. Bubbler bottle contains 500 mL of 1.5 mM $\text{NaH}^{13}\text{CO}_3$, air flow is $\sim 1\text{L}/\text{min}$. **B.** ^{13}C enrichment of *C. reinhardtii* biomass labelled by this approach after 24 hours.

- 3 The rate of label release in this method can be tuned by buffering the bubbler solution at different pH values (fig. 4A). In lower pH solutions, carbonate equilibria (eq. 1) shift towards a higher fraction of the DIC pool speciated as dissolved CO_2 , which increases the rate at which that DIC exchanges into the bubbled air. We identified solution conditions suitable for experiments requiring rapid incorporation of ^{13}C into algal biomass over short timescales (minutes) and for steady-state release over longer timescales (days). With a bubbler solution buffered with phosphate at pH 7.5, we obtained ^{13}C enrichment detectable in algal biomass by IRMS in <15 minutes (fig. 4B, green). With a bubbler solution devoid of any additional buffer beyond the HCO_3^- itself, we obtained slower but longer-lived steady-state incorporation of the label (fig. 4B, orange).




Figure 3. Label release from the bubbler solution can be modulated by pH. **A.** Model describing the flux of CO₂ released from the bubbler solution when buffered at different pH values, presented as the fraction of initial DIC leaving the solution over time. At lower pH values, the DIC leaves the solution more rapidly. **B.** Data from four separate experiments illustrating the biomass ¹³C enrichment of *C. reinhardtii* cultures over time, with bubblers buffered with 10x excess H_xK_xPO₄ at pH 7.5 (green) vs no additional buffer (orange), examined over long (diamonds) or short (circles) timecourses. At pH 7.5, the equilibrium DIC concentration is >10x lower than the DIC added to the bottle as H¹³CO₃⁻, so most of the label leaves the solution rapidly. This results in rapid labelling of biomass but also rapid depletion of label from the system. Without a separate buffer, DIC leaving the system drives the pH up to the pH at which the DIC concentration in the solution is at equilibrium with the air, and then exchanges with the air at a steady-state rate. This results in a longer-lived signal of label in biomass.