



Aug 30, 2024

Version 2

JAX-Sen: Collection and shipment of specimen for single-nuclei RNA sequencing (snRNA-Seq) V.2

DOI

dx.doi.org/10.17504/protocols.io.3byl4982jgo5/v2

Laura Robinson¹, Susan Sheehan¹, Gaven Garland¹, Ron Korstanje¹

¹The Jackson Laboratory, Bar Harbor, ME, USA

Cellular Senescence Net...



Harshpreet Chandok

The Jackson Laboratory

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4982jgo5/v2>

Protocol Citation: Laura Robinson, Susan Sheehan, Gaven Garland, Ron Korstanje 2024. JAX-Sen: Collection and shipment of specimen for single-nuclei RNA sequencing (snRNA-Seq). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.3byl4982jgo5/v2> Version created by **Harshpreet Chandok**



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 30, 2024

Last Modified: August 30, 2024

Protocol Integer ID: 106767

Keywords: JAX-Sen, jackson laboratory for genomic medicine, nuclei rna, rna, sequencing, snrna, details on specimen collection, genomic medicine, jackson laboratory, specimen collection, shipment of specimen, sennet consortium, nuclei, sen project in the sennet consortium, seq, shipment to the robson laboratory, part of the jax, sen project, robson laboratory

Funders Acknowledgements:

National Institute on Aging (NIA) JAX-Sen Senescence Tissue Mapping Center

Grant ID: U54 AG079753

Abstract

These samples are part of the JAX-Sen project in the SenNet Consortium. Here we provide details on specimen collection and shipment to the Robson laboratory at The Jackson Laboratory for Genomic Medicine (JAX-GM) in Farmington, CT for its processing for single-nuclei RNA-sequencing (snRNA-seq).

Troubleshooting

Reagents and Materials

- 1
 - 2mL Cryovials/ screw-cap tubes
 - ice cold 1X PBS
 - Kimwipes
 - Petri dishes
 - Liquid Nitrogen
 - Dry ice
 - Tweezers (clean, sterile)

Quality Key Points:

- 2
 - The tissue specimen should be kept at 4 degrees Celsius and RNase-free until snap-freezing and maintained at -80 degrees C or on dry ice thereafter, until it is processed at JAX.
 - It is crucial to not store the tissue specimen at RT to avoid freeze-thaw cycles, any cell death, and tissue and/or RNA degradation.

Procedure:

- 3 **Timeline:** The daily shipping deadline at Jax BH is 12:00 noon. So, harvest the tissues before and closer to noon.

- 4 **Collection/Harvest:**

- 4.1 Animal was euthanized via cervical dislocation.
- 4.2 Animal was pinned to a necropsy tray which remained on wet ice throughout the harvest.
- 4.3 Animal was perfused with 20ml cold PBS.
- 4.4 Heart, pancreas, and kidney was collected and the specimen was collected and the specimen was rinsed with ice-cold 1x PBS in a Petri dish to remove blood and other debris.
- 4.5 Pat the specimen dry using clean Kimwipes and freeze as quickly as possible after harvest.



4.6 To flash freeze (snap freeze), place the tissue in a cryo-tube and submerge the tube in liquid nitrogen or immerse deeply in dry ice.

4.7 Wait at least 2-3 minutes for the tissue to completely freeze, and transfer the tube containing the tissue to -80 degrees C.

4.8 Keep at -80 degrees Celsius thereafter, until shipping.

5 **Shipment:**

5.1 Place sample tubes in a plastic box (cardboard boxes insulate the samples from the cold ice) or in double Ziplock bags after checking that they are completely sealed.

5.2 Ship the sample box on dry ice (O/N shipping) to:
The Jackson Laboratory for Genomic Medicine,
Farmington, CT, 06032