



Dec 13, 2024

JAX-Sen: Mouse Heart and Kidney dissociation for bulk RNA sequencing

DOI

<https://dx.doi.org/10.17504/protocols.io.5jyl8d24rg2w/v1>

Patrick Fleming¹, Jessica Garofalo¹, Ramalakshmi Ramasamy¹, Juliana Alcoforado Diniz¹, Paul Robson^{1,2}

¹The Jackson Laboratory for Genomic Medicine, Farmington, CT, USA;

²Department of Genetics and Genome Sciences, University of Connecticut School of Medicine, Farmington, CT, USA

Cellular Senescence Net...



Harshpreet Chandok

The Jackson Laboratory

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl8d24rg2w/v1>

Protocol Citation: Patrick Fleming, Jessica Garofalo, Ramalakshmi Ramasamy, Juliana Alcoforado Diniz, Paul Robson 2024. JAX-Sen: Mouse Heart and Kidney dissociation for bulk RNA sequencing. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.5jyl8d24rg2w/v1>



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 13, 2024

Last Modified: December 13, 2024

Protocol Integer ID: 115152

Keywords: JAX-Sen, kidney dissociation for bulk rna, bulk rna, rna, sequencing, senescence, dissociation of the heart, kidney dissociation, mouse heart, sen project in the sennet consortium, sennet consortium

Funders Acknowledgements:

National Institute on Aging (NIA) JAX-Sen Senescence Tissue Mapping Center Grant

Grant ID: U54 AG079753

Abstract

These samples are part of the JAX-Sen project in the SENNET Consortium. We aim to study and characterize senescence in the C57Bl/6 mouse heart and kidney. This protocol describes the dissociation of the heart and kidney before library preparation and sequencing for Bulk RNA-seq.

Reagents and Materials

- 1
 - BeadbugTM 6 Homogenizer (recommended)
 - BeadBugTM prefilled tubes, 2.0 mL capacity (Sigma Aldrich, Z763772-50EA) (2 tubes per sample)
 - RNeasy Plus Mini Kit (Qiagen, 74134) (50 reactions)
 - Qias shredder tubes (Qiagen, 79654) (50 tubes total)
 - 2-mercaptoethanol (Gibco, 21985023)
 - 2 mL tubes (1 tube per sample)
 - 6cm dishes
 - Scalpel/tweezers

Procedure

- 2 Thaw fresh frozen sample on ice for 30 minutes.
- 2.1 Prepare two Beadbug tubes per sample by adding 600 uL RLT buffer (supplied) to each tube. Ensure that 2-mercaptoethanol has been added to RLT buffer before use as per RNeasy Plus Mini Kit instructions.
- 2.2 On a dish, mince individual organs to ~2 mm pieces using a scalpel
- 2.3 Add half of the minced organ to each prepared beadbug tube
- 2.4 Using the Beadbug tissue homogenizer, run the sample through two cycles of 90 seconds set at 3000 speed
- 2.5 Remove the lysate from the beadbug tube and collect in a 2 mL tube. Avoid pieces of tissue. Lysate will likely contain many bubbles, collect as much sample as possible. Briefly spinning the beadbug tube using a minifuge may aid in lysate collection
- 2.6 Add 600 uL RLT buffer with 2-mercaptoethanol to the beadbug tube containing the residual of each sample.
- 2.7 Using the Beadbug tissue homogenizer, run the sample through two cycles of 90 seconds set at 3000 speed.



- 2.8 Remove the residual lysate from the beadbug tube and combine the lysate in the same 2 mL tube as initial tissue lysate. Avoid pieces of tissue. Lysate will likely contain many bubbles, collect as much sample as possible. Briefly spinning the beadbug tube using a minifuge may aid in lysate collection.
- 2.9 At this step in the protocol, whole samples (~1200 uL) have been pooled together in a single 2 mL tube.
- 2.10 For kidney, freeze a half volume (~600 uL) of tissue lysate at -80°C. For heart, freeze three quarters a volume (~900 uL) of tissue lysate at -80°C.
- 2.11 Begin to follow the RNeasy Plus Mini Kit Protocol as directed by manufacturer. For kidney, proceed with ~600 uL (about half) of sample lysate. For heart, proceed with ~300 uL (about a quarter) of sample lysate.
- 2.12 Add lysate to a Qias shredder tube.
- 2.13 Centrifuge 30 seconds x 8000g.
- 2.14 Transfer flowthrough to a gDNA binding column.
- 2.15 Centrifuge 30 seconds x 8000g.
- 2.16 Discard gDNA column, save flowthrough.
- 2.17 Add 600 uL of 70% ethanol to the sample and mix by pipetting.
- 2.18 Transfer no more than 700 uL of sample to a RNeasy spin column (supplied).
- 2.19 Centrifuge 30 seconds x 8000g.
- 2.20 Repeat steps 19-20 until all sample volume has passed through the RNeasy spin column.



- 2.21 Add 600 uL buffer RW1 (supplied) to RNeasy spin column. Ensure to add pure ethanol to buffer RW1 as directed in manufacturer instructions.
- 2.22 Centrifuge 30 seconds x 8000g, discard flowthrough.
- 2.23 Add 500 uL buffer RPE (supplied) to RNeasy spin column.
- 2.24 Centrifuge 30 seconds x 8000g, discard flowthrough.
- 2.25 Add 500 uL buffer RPE to RNeasy spin column.
- 2.26 Centrifuge 2 minutes x 8000g, discard flowthrough.
- 2.27 Transfer RNeasy spin column to a new 2 mL tube (supplied).
- 2.28 Centrifuge 1 minute at full speed.
- 2.29 Transfer spin column to a final 1.5 mL collection tube (supplied).
- 2.30 Add 30 uL of RNase-free water (elution buffer, supplied).
- 2.31 Centrifuge 1 minute at 8000g.
- 2.32 Measure RNA concentration and record.