

Aug 31, 2024

JAX-Sen: Mouse heart dissociation for single-cell RNA sequencing

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vzj55xlx1/v1>

Ramalakshmi Ramasamy¹, Juliana Alcoforado Diniz¹, Patrick Fleming¹, Jessica Garofalo¹, Paul Robson^{1,2}

¹The Jackson Laboratory for Genomic Medicine, Farmington, CT, USA;

²Department of Genetics and Genome Sciences, University of Connecticut School of Medicine, Farmington, CT, USA

Cellular Senescence Net...



Harshpreet Chandok

The Jackson Laboratory

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzj55xlx1/v1>

Protocol Citation: Ramalakshmi Ramasamy, Juliana Alcoforado Diniz, Patrick Fleming, Jessica Garofalo, Paul Robson 2024. JAX-Sen: Mouse heart dissociation for single-cell RNA sequencing. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.rm7vzj55xlx1/v1>



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 03, 2024

Last Modified: August 31, 2024

Protocol Integer ID: 102774

Keywords: JAX-Sen, mouse heart dissociation, cell dissociation of the heart, cell rna, using scrnasequencing, mouse heart, cell dissociation, senescence, sequencing, heart, sen project in the sennet consortium, rna

Funders Acknowledgements:

National Institute on Aging (NIA) JAX-Sen Senescence Tissue Mapping Center Grant

Grant ID: ID: U54 AG079753

Abstract

These samples are part of the JAX-Sen project in the SenNet Consortium. We aim to study and characterize senescence in the C57Bl/6 mouse heart. We compare young (6 months old) and old (24 months old) mouse hearts using scRNAsequencing. This protocol describes the single-cell dissociation of the hearts before library preparation and sequencing.

Reagents and Materials

- 1
 - 5 ml Protein lobind tubes
 - 1.5 or 2 ml Protein lobind tubes
 - Micro scissors/ scalpel
 - Ice-cold PBS
 - Wide-bore pipette tips
 - 0.25% Trypsin/EDTA (Gibco, 25200056). - 2 ml/sample
 - 20 mg/mL collagenase A and B (Sigma, 10103578001, 11088807001) – 2 ml/sample
 - DMEM + 10% FBS
 - 100 µm cell strainer (Corning, 431752)
 - HBSS (Ca²⁺/Mg²⁺ free) (Gibco, 14170120) – 2 ml/sample
 - ACK lysis buffer – 2 ml/sample
 - SS_buffer – 5 ml/sample
 - 40 um cell strainer
 - 70 uM strainer

SS_buffer

	A	B	C	D	E	F	G	H	I
		Stock conc.	Workin g conc.	50 ml	100 ml	250 ml	500 ml	20 ml	
	EDTA	0.5 M	2 mM	0.2	0.4	1	2	0.08	ml
	BSA	10%	2%	10	20	50	100	4	ml
	FBS	100%	15%	7.5	15	37.5	75	3	ml
	DPBS	1X	make up	32.2	64.4	161	322	12.88	ml

Procedure

- 2 Collect shipment from the warehouse.
- 2.1 Fill up ice buckets and prep them with reagents:
Bucket 1(small): PBS, HBSS, DMEM+FBS, SS_buffer



- Bucket 2 (large): For mincing tissue
- 2.2 Confirm if the sample is cold. Take pictures.
 - 2.3 Place tubes on ice and mince the tissue (<1mm) with sterile micro scissors.
 - 2.4 After mincing, transfer the sample into labeled 5 ml tubes.
 - 2.5 Centrifuge at **500 x g, 4°C, 5 min** to get rid of the supernatant.
 - 2.6 Resuspend the pellet in **ice-cold PBS (Ca+ and Mg+ free)** and pipette up and down gently with wide-bore pipette tip.
 - 2.7 Centrifuge at **500 x g, 4°C, 5 min** to get rid of the supernatant.
 - 2.8 Resuspend the pellet in **ice-cold PBS (Ca+ and Mg+ free)** and pipette up and down gently with wide-bore pipette tip.
 - 2.9 Centrifuge at **500 x g, 4°C, 5 min** to get rid of the supernatant.
 - 2.10 Add 2 ml of **0.25% Trypsin/EDTA** to the pellet and enzymatically digest at **37°C for 10 min** with slow rocking. ***Check at 5 mins for dissociation, under the microscope.***
 - 2.11 Stop the digestion by adding an equal volume (1ml) of **DMEM/10% FBS**.
 - 2.12 Centrifuge at **500 x g, 4°C, 5 min** to get rid of the supernatant.
 - 2.13 Resuspend the pellet 2 ml of **20 mg/mL collagenase A and B** mixture and incubate the samples at 37 °C until the tissues were into single cells.
 - 2.14 Vigorously pipette samples with 1 ml wide-bore pipette tips or if needed, regular tips.



- 2.15 Filter cells through **MACS SmartStrainers (70um)** and wash 1 ml **HBSS (Ca²⁺/Mg²⁺ free)** through the strainer.
- 2.16 Centrifuge at **500 x g, 4°C, 5 min** to get rid of the supernatant.
- 2.17 Resuspend the cell pellet in 2 mL **ACK Lysing buffer**, and incubate on **ice for 3 min** to exclude RBCs.
- 2.18 Add 5 mL **SS_buffer** to stop ACK Lysing buffer reaction, then collect cells by centrifugation **500 x g, 4°C, 5 min**.
- 2.19 Resuspend the pellet in **1ml of HBSS (Ca²⁺/Mg²⁺ free)** (Gibco, 14170120).
- 2.20 Filter the cell suspension through **Flowmi 40um cell strainers** tube to remove debris and non-dissociated tissue fragments.
- 2.21 Collect cells by centrifugation **500 x g, 4°C, 5 min** and proceed to Flex protocol.