

Aug 31, 2024

JAX-Sen: Mouse pancreas dissociation for single-cell RNA sequencing

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov19zzylr2/v1>

Ramalakshmi Ramasamy¹, Juliana Alcoforado Diniz¹, Patrick Fleming¹, Jessica Garofalo¹, Dylan Baker¹, Paul Robson^{1,2}

¹The Jackson Laboratory for Genomic Medicine, Farmington, CT, USA;

²Department of Genetics and Genome Sciences, University of Connecticut School of Medicine, Farmington, CT, USA

Cellular Senescence Net...



Harshpreet Chandok

The Jackson Laboratory

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov19zzylr2/v1>

Protocol Citation: Ramalakshmi Ramasamy, Juliana Alcoforado Diniz, Patrick Fleming, Jessica Garofalo, Dylan Baker, Paul Robson 2024. JAX-Sen: Mouse pancreas dissociation for single-cell RNA sequencing. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.ewov19zzylr2/v1>



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 03, 2024

Last Modified: August 31, 2024

Protocol Integer ID: 102769

Keywords: JAX-Sen, cell dissociation of the pancrea, mouse pancrea, cell rna, using scrnasequencing, pancrea, cell dissociation, sequencing, senescence, sen project in the sennet consortium, rna

Funders Acknowledgements:

National Institute on Aging (NIA) JAX-Sen Senescence Tissue Mapping Center Grant

Grant ID: ID: U54 AG079753

Abstract

These samples are part of the JAX-Sen project in the SenNeT Consortium. We aim to study and characterize senescence in the C57Bl/6 mouse pancreas. We compare young (6 months old) and old (24 months old) mouse pancreas using scRNAsequencing. This protocol describes the single-cell dissociation of the pancreas before library preparation and sequencing.

Reagents and Materials

- 1
 - 15 ml conical tubes – 4 per sample
 - 1.5 or 2 ml Protein lobind tubes
 - Micro scissors/ scalpel
 - Ice-cold PBS
 - Histopaque-1077 – Sigma H8889
 - Wide-bore pipette tips
 - Collagenase solution
 - Wash buffer
 - 100 µm cell strainer (Corning, 431752)
 - ACK lysis buffer – 2 ml/sample
 - SS_buffer – 5 ml/sample
 - 40 um cell strainer
 - 70 uM strainer
 - 10 ml syringe and 18G needle if available

2 **SS_buffer:**

	A	B	C	D	E	F	G	H	I
		Stock conc.	Workin g conc.	50 ml	100 ml	250 ml	500 ml	20 ml	
	EDTA	0.5 M	2 mM	0.2	0.4	1	2	0.08	ml
	BSA	10%	2%	10	20	50	100	4	ml
	FBS	100%	15%	7.5	15	37.5	75	3	ml
	DPBS	1X	make up	32.2	64.4	161	322	12.88	ml

2.1 **Collagenase solution:**

For each mouse, dissolve 4mg collagenase powder in 10mL serum-free Medium 199 (or MEM). Prepare fresh daily. Fill 5mL syringe with 3mL collagenase P (Roche Applied Science 11 249 002 001) and 15mL tube with 5mL collagenase; label tubes and store on ice until ready to use. Note: each lot of collagenase must be tested experimentally to determine optimal working concentration. Use 0.4mg/mL as starting point.




- **Wash buffer:**
DMEM containing 20% FBS.
- **Stopping media:** DMEM + 10% FBS + 1:100 dilution Glutamate CMRL media

Procedure

- 2.2 Collect shipment from the warehouse.
- 2.3 Fill up ice buckets and prep them with reagents:
Bucket 1(small): PBS, HBSS, DMEM+FBS, SS_buffer
Bucket 2 (large): For mincing tissue
- 2.4 Confirm if the sample is cold. Take pictures.
- 2.5 Place tubes on ice and mince the tissue (<1mm) with sterile micro scissors.
- 2.6 After mincing, transfer the sample into labeled 15 ml tubes containing **5 ml of 0.4mg/ml collagenase solution**. Incubate at 37°C for 10 mins.
- 2.7 Remove and shake tubes vigorously until pancreatic tissue is completely dissociated with no large chunks visible, 10-20 seconds. Alternatively, triturate using 10 mL pipet until well-dispersed.
- 2.8 **Immediately** dilute digest with cold wash buffer to 15mL total volume. Invert several times to mix. **Keep all the tubes cold.**
- 2.9 Centrifuge at **300 x g, 4°C, 1 min** to pellet pancreatic tissue.
- 2.10 Decant supernatant into waste beaker. Add **12mL wash buffer** and resuspend pellet by shaking or with 10mL pipette.
- 2.11 Repeat wash (steps #2.9-2.10) **twice**. After the final spin, decant supernatant, then invert open tubes onto a paper towel and tap gently to remove all residual buffer.
- 2.12 Add **5mL Histopaque-1077** and gently resuspend pellet using 10mL pipette.

- 2.13 Hold tube at 45° angle and *slowly* layer **5mL wash buffer** on top of gradient layer (**10mL syringe** is helpful for slowly dispensing buffer). Do not allow buffer to mix with Histopaque! Very carefully transfer tubes to centrifuge.
- 2.14 **Turn off centrifuge brake.** Centrifuge at 300xg/4°C for 25 minutes.
- 2.15 *Islets should be visible floating at interface between Histopaque and wash buffer, or slightly below.* Remove islets from interface using sterile transfer pipette and transfer to new 15mL conical tube. Centrifuge at 300xg/4°C for 25 minutes. *Other cells (mostly acinar): in the pellet.*



Split the 2 layers and treat them as separate samples going forward.

Proceed to islets protocol and acinar+ducts protocols that we use for human pancreatic cells.

Acinars and Ducts:

- 2.16 Wash pellets from histopaq gradient **3x** with 10 ml cold PBS
- 2.17 Add 3 ml TrypLE E to the pellet and incubate at 37C (waterbath), mixing with pipette approximately every 5 mins. Check at 10 mins and every 5 mins following for single cells. (Ducts should take 10-15 mins, and acinar should take slightly longer ~20 mins. But we have a mixture here, so check under the microscope and stop reaction accordingly)
- 2.18 Add **3X stopping media (9 ml)**.



- 2.19 Centrifuge at 1300 rpm/4°C for 2 minutes. Discard supernatant.
- 2.20 Wash again with **3 ml stopping media** and pellet cells. Resuspend in 1 ml of stop media.
- 2.21 Filter through 40 uM Flowmi filter.
- 2.22 Count cells and proceed to Flex protocol.

Islets:

- 2.23 Centrifuge the islets suspension at 130g/RT for 5 minutes.
- 2.24 Add 1.5 ml accutase, mix with pipette and incubate at 37C. Mix every 2 mins, checking at 6 mins.
- 2.25 Add CMRL media to approximately 15 ml total volume.
- 2.26 Centrifuge at 500xg/4°C for 5 mins.
- 2.27 Resuspend in 1 ml of CMRI media.
- 2.28 Filter through 40 uM filter.
- 2.29 Count cells and proceed to Flex protocol.