

Apr 28, 2025

# JAX-Sen: Mouse placenta single-cell dissociation for single-cell RNA sequencing

DOI

<https://dx.doi.org/10.17504/protocols.io.x54v9korml3e/v1>

Ramalakshmi Ramasamy<sup>1</sup>, Juliana Alcoforado Diniz<sup>1</sup>, Jessica Garofalo<sup>1</sup>, Patrick Fleming<sup>1</sup>, Paul Robson<sup>1,2</sup>

<sup>1</sup>The Jackson Laboratory for Genomic Medicine, Farmington, CT, USA;

<sup>2</sup>Department of Genetics and Genome Sciences, University of Connecticut School of Medicine, Farmington, CT, USA

Cellular Senescence Net...



Harshpreet Chandok

The Jackson Laboratory

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.x54v9korml3e/v1>

**Protocol Citation:** Ramalakshmi Ramasamy, Juliana Alcoforado Diniz, Jessica Garofalo, Patrick Fleming, Paul Robson 2025. JAX-Sen: Mouse placenta single-cell dissociation for single-cell RNA sequencing. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.x54v9korml3e/v1>



**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** April 26, 2025

**Last Modified:** April 28, 2025

**Protocol Integer ID:** 174722

**Keywords:** cell dissociation of the placenta tissue, cell rna, placenta tissue, cell dissociation, senescence, sequencing, cell, single cell rna

**Funders Acknowledgements:**

National Institute on Aging (NIA) JAX-Sen Senescence Tissue Mapping Center Grant

Grant ID: U54 AG079753

## Abstract

These samples are part of the JAX-Sen project in the SenNet Consortium. We aim to study and characterize senescence in the C57Bl/6 mouse heart. This protocol describes the single-cell dissociation of the placenta tissue before library preparation and sequencing.

## Reagents and Materials:

- 1
  1. 6 cm Petri dishes
  2. Ice-cold MACS buffer
  3. 5 ml Protein lobind tubes
  4. Micro scissors/ scalpel
  5. Ice-cold PBS
  6. Wide-bore pipette tips
  7. TrypLE express (10 volumes)
  8. Dispase (0.4 g/ml)
  9. Collagenase
  10. DMEM + 10% FBS
  11. ACK lysis buffer
  12. CS\_buffer
  13. 40 um cell strainer (Flowmi)
  14. 70 uM strainer

## Procedure:

- 2 Add sample labels and collect shipment from the warehouse.

## Procedure:

- 3 Fill up ice buckets and prep them with reagents:  
Bucket 1: PBS, MACS buffer, DMEM, DMEM+FBS, SS\_buffer  
Bucket 2: For mincing tissue, either in 6 cm dishes or in tubes – protocol A  
Bucket 3: For mincing tissue, either in 6 cm dishes or in tubes – protocol B  
Bucket 4: extra ice
- 4 Confirm if the sample is cold.  
Notes:  
*Optional: Place the sample on a dish with **MACS buffer** and take pictures.*
- 5 Place tubes on ice and mince the tissue (<1mm) with sterile micro scissors (or) mince using a scalpel on 6 cm dish w/MACS buffer on ice.
- 6 After mincing, transfer the sample into labeled 5 ml/ 10 ml tubes.
- 7 Centrifuge at 800 x g, 4°C, **5 min** to get rid of the supernatant.



- 8 Resuspend the pellet in **ice-cold PBS** and pipette up and down gently with wide bore pipette tip.
- 9 Centrifuge at 800 x g, 4°C, **5 min** to get rid of the supernatant.
- 10 Resuspend the pellet in **ice-cold PBS** and pipette up and down gently with wide bore pipette tip.
- 11 Centrifuge at 800 x g, 4°C, **5 min** to get rid of the supernatant.
- 12 Add 10 volumes (~500 uL) of **[1937.5 uL TrypLE + 50 uL Dispase + 12.5 uL RNase inhibitor + 250 uL DNaseI]** to the pellet and enzymatically digest at **4°C for 35 min** with slow rocking.  
***Check every 5-10 mins minutes for dissociation, under the microscope.***
- 13 Stop the digestion by adding an equal volume (~500 uL) of **DMEM/10% FBS**.
- 14 Filter cells through MACS SmartStrainers (70um)
- 15 Wash 1.5 ml DMEM through the strainer. -----à this is the filtrate cell suspension  
Filtrate: Leave on ice.  
Residue:
  1. Collect the residue with wide bore pipette tips and transfer them to a new 2 ml conical tube.
  2. Add collagenase Type IV (0.5 mg/mL) and DNase I (0.2 mg/mL) to the residue and enzymatically digest for 10 min at 37 °C.  
[500 ul DNaseI + 1.25 mg collagenase + 12.5 uL RNase inhibitor + 1862.5 uL DPBS (make up to 2.5 ml)]
  3. Filter cells through MACS SmartStrainers (70um)
  4. Wash 1.5 ml DMEM through the strainer. —→ this is the residue cell suspension
- 16 Pool the filtrate and residue cell suspensions and proceed to step 16



- 17 Centrifuge at 800 x g, 4°C, **5 min** to get rid of the supernatant.
- 18 Resuspend the cell pellet in 2 mL ACK Lysing buffer, and incubate on ice for **3 min** to exclude RBCs.
- 19 Add **5 mL CS\_buffer** to stop ACK Lysing buffer reaction, then collect cells by centrifugation 800 x g, 4°C, **5 min**.
- 20 Resuspend the pellet in **1ml of CS\_buffer**.
- 21 Filter the cell suspension through 40um cell strainers tube to remove debris and non-dissociated tissue fragments.
- 22 Collect cells by centrifugation 800 x g, 4°C, **5 min** and proceed to 10X Flex protocol/ cell sorting.
- 23 CS\_buffer (Cell suspension buffer)

|  | A    | B           | C             | D     | E      | F      | G      | H     | I  |
|--|------|-------------|---------------|-------|--------|--------|--------|-------|----|
|  |      | Stock conc. | Working conc. | 50 ml | 100 ml | 250 ml | 500 ml | 20 ml |    |
|  | EDTA | 0.5 M       | 2mM           | 0.2   | 0.4    | 1      | 2      | 0.08  | ml |
|  | BSA  | 10%         | 2%            | 10    | 20     | 50     | 100    | 4     | ml |
|  | FBS  | 100%        | 15%           | 7.5   | 15     | 37.5   | 75     | 3     | ml |
|  | DPBS | 1X          | make up       | 32.2  | 64.4   | 161    | 322    | 12.88 | ml |