

Jan 09, 2026

 JCVSeq

DOI

<https://dx.doi.org/10.17504/protocols.io.4r3l21w14g1y/v1>

Ellie Bourgikos<sup>1</sup>, Mallery I Breban<sup>1</sup>, Nicole M Feriancek<sup>1</sup>, Verity Hill<sup>1</sup>, Chantal Vogels<sup>1</sup>, Nathan Grubaugh<sup>1</sup>

<sup>1</sup>Department of Epidemiology of Microbial Diseases, Yale School of Public Health



Ellie Bourgikos

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l21w14g1y/v1>

**Protocol Citation:** Ellie Bourgikos, Mallery I Breban, Nicole M Feriancek, Verity Hill, Chantal Vogels, Nathan Grubaugh 2026. JCVSeq. protocols.io <https://dx.doi.org/10.17504/protocols.io.4r3l21w14g1y/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** November 27, 2025

**Last Modified:** January 09, 2026

**Protocol Integer ID:** 233667

**Keywords:** segments of the jcv genome, jcv genome, jcv diversity across lineage, specific small segment utr primer, segment specific small segment utr primer, large utr primer, jcv diversity, genome coverage, multiplexed pcr primer, separate primer set, developed jcv segment, primer design, multiplexed approach to primer design, sequencing protocol, jcv segment, specific primer, jcvseq in order, jcvseq, lineage, primer, bp amplicon, final primer, pcr

## Abstract

In order to capture all three segments of the JCV genome, we employed a multiplexed approach to primer design by combining previously developed JCV segment-specific primers and an internal set generated by PrimalScheme. To maximize genome coverage, we developed our own highly multiplexed PCR primers for the medium and large segments using PrimalScheme. Based on publicly available sequences, we chose to make two separate primer sets for the highly divergent A and B lineages, each forming 800 bp amplicons. Because the A and B primers share identical binding sites, they can be mixed together by pool during the sequencing protocol to cover JCV diversity across lineages. The final primer set consists of segment specific small segment UTR primers, segment specific medium and large UTR primers, and internal A and B lineage primers for the medium and large segments.

## Guidelines

It is recommended that steps performed up to amplicon generation be performed at a different workstation and with different equipment than steps post-amplicon generation.

## Materials

### Equipment

- 96-well format Thermocycler (two needed for 49+ sample runs)
- Qubit
- Bioanalyzer
- Magnetic rack fit for a 96-well PCR plate
- Magnetic rack fit for 1.5 / 2-mL microcentrifuge tubes
- Cold block fit for a 96-well PCR plate
- Pipettes (assorted sizes)

### Consumables

- 8-strip PCR tubes; may be substituted with PCR plates with heat-sealing film
- 1.5 / 2mL / 5 mL microcentrifuge tubes
- Filtered pipette tips (assorted sizes)
- Reservoirs
- Waste containers

### Reagents

- [Illumina COVIDSeq Test Kit](#)
- [Illumina Nextera XT DNA Library Preparation kit](#)
- [New England Biolabs NEBNext ARTIC SARS-CoV-2 Library Prep Kit \(Illumina\)](#)
- [Oxford Nanopore Technologies Rapid barcoding kit](#)

### Bioinformatics pipeline

The Grubaugh Lab's virus agnostic analysis pipeline can be adapted to perform consensus assembly and generate quality control plots for JCVSeq. The pipeline can be found here:

[https://github.com/grubaughlab/DENV\\_pipeline](https://github.com/grubaughlab/DENV_pipeline). To run for this primer scheme and assemble Jamestown Canyon virus genomes, replace reference and BED files with the following files.

### Reference genomes

 JCVSb.fasta 0B

 JCVSa.fasta 0B

 JCVMb.fasta 4KB

 JCVMa.fasta 4KB

 JCVLb.fasta 6KB

 JCVLa.fasta 6KB

### BED files

 JCVSb.bed 0B

 JCVSa.bed 0B

 JCVMb.bed 0B

 JCVMa.bed 0B

 JCVLb.bed 1KB

 JCVLa.bed 1KB

### Troubleshooting

## qPCR

- 1 Briefly vortex and centrifuge all reagents.
- 2 Prepare 10  $\mu\text{M}$  working stocks of the primers and probes below, by adding 10  $\mu\text{L}$  of 100  $\mu\text{M}$  stock to 90  $\mu\text{L}$  of nuclease-free water.

Primer/Probe Name	Sequence
CAES_JCV154_F	TAATGCAG CAAAAGCC AAAG
CAES_JCV307_R	AAGCCGAT GGATGGTA AGAT
CAES_JCV181_P	CGCTCGTA AACCGGAG CGGA

qPCR primers developed by the Connecticut Agricultural Experiment Station (CAES) as part of unpublished JCV research

- 3 On ice, prepare a master mix containing the following:

Reagent Name	Volume per Sample
Luna Master Mix	10 $\mu\text{L}$
RT	1 $\mu\text{L}$
CAES_JCV154_F	1 $\mu\text{L}$
CAES_JCV307_R	1 $\mu\text{L}$
CAES_JCV181_P	0.5 $\mu\text{L}$
Nuclease Free Water	1.5 $\mu\text{L}$

- 4 Place the 96-well PCR plate on a cold block and add 15  $\mu\text{L}$  of PCR master mix to each designated well.



- 5 Add 5  $\mu$ L of RNA to each designated well. Mix by pipetting, avoiding bubbles.
- 6 Add 5  $\mu$ L of positive control and non-template control to the final two wells. Mix by pipetting, avoiding bubbles.
- 6.1 Our recommended positive control is a gene block targeting the small segment. Its sequence is:

Name	Sequence
JCVS_gBlock	CAGATGCA GGGTTTGT GGCATTAT GGCTGACC ATGGGGAA TCTGTCAG TCTGTCAG CCGTTAGG ATCTTCTT CCTTAATG CAGCAAAA GCCAAGGC TGCTCTCG CTCGTAAA CCGGAGCG GAAAGCTA CTCCTAAA TTTGGAGA GTGGCAGG TGGA AATT GTCAATAAT CATT TTTCC TGGA AACA GGAACAAC CCAATTGG TAACAACG ATCTTACC ATCCATCG GCTTTCAG GATATCTAG CTAGATGG GTTCTTGA GCATTTTA CTACAGAT GATGATGA GTCCCAGA GAGA ACTC ATAAGGAG CACCATCA TTAATCCA ATTGCAGA GTCCAATG GCATT CAT TGGAACAA TGGCCCAG AGATTTATC TTTCATTCT



Name	Sequence
	TTCCAGGA ACAGAAAT GTTTTTGG AAATTTTC AAATTCTAT CCCTTGAC CATTGGAA TTTACAGA GTCAAGCA TGGTATGAT GGACCCTC AGTATCTGA AGAAGGCT CTCAGACA GCGC

7 Seal with a transparent plastic PCR seal. Centrifuge in the plate spinner.

8 Set the thermocycler to read FAM.

9 Run the following thermocycler conditions:

Temperature	Time	Cycles
55 °C	10 minutes	
95 °C	1 minute	44 cycles
95 °C	10 seconds	(included in above 44)
60 °C	30 seconds	(included in above 44)
	<b>READ PLATE</b>	

## Primer Pooling and Preparation

10 Combine the following primers at equal volume into a 20  $\mu$ M pool. This will be AB1, containing pool 1 of primers covering diversity for both major lineages A and B.



Primer Name	Sequence
NYS_JCVSM F1	AGTGTACTA CCAAGTAT AGAAAACG TTCA
NYS_JCVSM 991R1	AGTAGTGT GCTCCACT GAATACATT TAA
NYS_JCVMF 1	AGTAGTGTA CTACCAAG TATAGAAAA CGTT
Yale_JCVM8 37R1B	GGGTGATA AACTAACC CGCACA
Yale_JCVM8 37R1A	GGGTGATA CACCAAAC CACATAGT
Yale_JCVM1 257F3B	GTGATATGT ATCATGAA AAAGCCGG T
Yale_JCVM1 257F3A	GTGACATG TACCATGA AAAAGCTG G
Yale_JCVM2 057R3B	TGTCTTTT CGGGCCAT AGCTTC
Yale_JCVM2 057R3A	TATCTTTGC GAGCCATA GCTTCTG
Yale_JCVM2 611F5B	TGGGAACT GTGATGTT CAAGAAAA TG
Yale_JCVM2 611F5A	TGGAAATT GCAATGTT CAAGAAAA TGATT
Yale_JCVM3 440R5B	TGGGAGTT ATAGCATTT AAGTTTGT GC
Yale_JCVM3 439R5A	AGGAGTTA CAGCATTT



Primer Name	Sequence
	AAATTAGT GCAGTA
Yale_JCVM3 680F7B	GACAACGA TTACCAAG CTTGCAA
Yale_JCVM3 680F7A	GACAATGA TTATCAAG CTTGCAA TTTCT
NYS_JCVM4 510R2	AGTAGTGT GCTACCAA GTATATCTA AATGA
NYS_JCVLF1	AGTAGTGTA CTCCTATTT ACAAAAC TACAAATAC
Yale_JCVL8 73R1B	GGGCTGTT TGTAATCTC CAGTGG
Yale_JCVL8 73R1A	GGGCTGTT TGTAGTTC CCTGT
Yale_JCVL11 77F3	AGCTGGCG ACAAGTTAT GAATAAGA
Yale_JCVL2 014R3B	GTGCAGGC TCAGTTAG AGATAACA
Yale_JCVL2 014R3A	GTGCAGGC TCAGTTAA GGACAA
Yale_JCVL2 315F5B	CCAAAGGT CTGCATGA GAAGCA
Yale_JCVL2 315F5A	CTAAGGGT TTACATGAG AAACATCA TGTT
Yale_JCVL31 21R5B	CCCTATTT TTTTGCCT AGTAGTTT CAACT
Yale_JCVL31 21R5A	CCCTATTC TTTTGTCTT



Primer Name	Sequence
	GTAGTCTC AAC
Yale_JCVL3 638F7B	CAATAGTT CAAGATAA GGCCCCTG A
Yale_JCVL3 638F7A	CGATAGTT CAAGACAA GGCTCCA
Yale_JCVL4 432R7B	CATTCATG CCTCCTGG TGATGATA
Yale_JCVL4 432R7A	CATTCATTC CTCCTGGT GATGACA
Yale_JCVL4 975F9B	GCAGACCC AACAGAGA TGTCAA
Yale_JCVL4 975F9A	GCTGATCC AACAGAGA TGTC AAGA
Yale_JCVL5 804R9B	GCCATATTT TCGAATTT TAGACCAT GC
Yale_JCVL5 804R9A	GCCATATTT TCAAATTT CAAGCCAT GT
Yale_JCVL6 087F11B	AAAGAAAG CACATTTTA GCAAAATG GTATC
Yale_JCVL6 087F11A	AAAGAAGG CGCATTTTC AGCAAAA
NYS_JCVL6 957R2	AGTGTGCT CCTATTTAC AAATATATA CTATAAGC

Primers developed using PrimalScheme, in collaboration with New York State Wadsworth Center

- 11 Combine the following primers at equal volume into a 20  $\mu$ M pool. This will be AB2, containing pool 1 of primers covering diversity for both major lineages A and B.

Primer Name	Sequence
Yale_JCVM5 82F2B	ACCAACAT ATGTCATGT GTACGGT
Yale_JCVM5 82F2A	ATCAGCATA TGTCATGC GTACGAT
Yale_JCVM1 367R2B	TGTGTACTA TCAATCCA GTAACATC TTCA
Yale_JCVM1 367R2A	TATGAACC ATTAACCC AGTAATGTC TTCA
Yale_JCVM1 910F4B	TGGGATTT TGCAAATG AAATGAAG ACA
Yale_JCVM1 910F4A	TGGGATTT CGCAAATG AAATGAAA ACT
Yale_JCVM2 747R4B	TACTGCAT CTGGGACT AAGGCA
Yale_JCVM2 747R4A	TGTTGCAC CTAGGGCT TAGACA
Yale_JCVM3 303F6B	TCTTCGGG TCCTGTCA AGACATTA
Yale_JCVM3 299F6A	TGTGTATTT GGATCCTG TCAAGATAT CATA
Yale_JCVM4 134R6B	TGGTCTAC TGGTGCAA GCTCTAG
Yale_JCVM4 134R6A	TGGTCAAC TGGTGCAA GTTCTAA



Primer Name	Sequence
Yale_JCVL7 46F2B	CATATGAAT CAGAAAGA TGGAACAC TAACC
Yale_JCVL7 46F2A	CATACGAA TCAGAGAG GTGGAACA
Yale_JCVL15 77R2B	CATTGCCA GAATTTGG TCTCAAAA ATAG
Yale_JCVL15 77R2A	CACTGCCA GAACCTAG TCTCAAA
Yale_JCVL16 80F4B	TGCATTAGT TTACCCTT CAGCTGAT
Yale_JCVL16 80F4A	TGCACTAG TCTACCCA TCTGCA
Yale_JCVL2 453R4B	AGGATTTG TAAATTAAC TGTTTGCT TGGT
Yale_JCVL2 453R4A	AGGATCTG CAAATTGA CAGTTTGT TT
Yale_JCVL2 940F6B	TGTAGGAG AATATGAAG CCAAAATG TGC
Yale_JCVL2 940F6A	TGTGGGAG AATACGAG GCTAAAAT G
Yale_JCVL3 744R6	TGTCTTTTT CATATTGG CTTGACAT CC
Yale_JCVL4 301F8B	TGGGTGAG ACAAGTGA TATGAGGG
Yale_JCVL4 301F8A	TGGGTGAA ACAAGCGA TATGAGAG



Primer Name	Sequence
Yale_JCVL51 27R8B	TCTTGTCAT TTCTTTCA ATTCAAAC ACAAT
Yale_JCVL51 26R8A	CTCGTCAT CTCCTTCA GTTCAAAA AC
Yale_JCVL5 681F10B	TTGGTGAA GACAACAA GCTAACTT ATTC
Yale_JCVL5 681F10A	TCGGTGAA GACAATAA GCTAACCT AC
Yale_JCVL6 499R10B	TAGCCTCT AAGCAACC TAAATTTTC A
Yale_JCVL6 499R10A	TGGCTTCT AAACAACC TAAATTTTC TGG

Primers developed using PrimalScheme, in collaboration with New York State Wadsworth Center

## cDNA & Amplification

- 12 Combine the reagents from the table below to make PCR Master Mix. Multiply each volume by the number of samples (reagent overage already included) and create separate master mixes for each primer pool (AB1 and AB2)

Reagent Name	Volume per Sample
IPM	15 $\mu$ l
FSM	3.2 $\mu$ l
20 $\mu$ M primer pool	1.2 $\mu$ l
Nuclease Free Water	3.6 $\mu$ l

	Reagent Name	Volume per Sample
	RVT	1.0 $\mu$ l

- 13 Add 20  $\mu$ l of PCR Master Mix to each well of a 96-well plate.
- 14 Add 5  $\mu$ l of extracted RNA to each well of the 96-well plate (for a total volume of 25  $\mu$ l per well)
- 15 Seal and shake at 1600 rpm for 1 minute.
- 16 Centrifuge at 280  $\times$  g for 1 minute.
- 17 Place on thermocycler and run the PCR program to generate cDNA and amplify. The thermocycler conditions can be found below

	Temperature	Time	Cycles
	25 °C	5 minutes	
	50 °C	10 minutes	
	80 °C	5 minutes	
	98 °C	3 minutes	
	98 °C	15 seconds	35 cycles
	63 °C	5 minutes	(included in above 35)
	4 °C	hold	

Hybrid "one-step" cDNA generation and amplification thermocycler conditions

- 17.1 **This is a safe stopping point, marking the end of day 1 of the protocol. Store all samples at 4 °C**

## Tagment Amplicons and Cleanup

18 Using the Qubit, quantify the last column of samples (including both controls) from both amplicon sets.

19 Prepare the tagmentation master mix:

	Reagent Name	Volume per Sample
	TB1	12 $\mu$ l
	EBLTS	4 $\mu$ l
	Nuclease Free Water	20 $\mu$ l

20 Add 30  $\mu$ l of the tagmentation master mix to 96 new PCR tubes.

20.1 Keep these new tubes on a cold block.

21 Combine 10  $\mu$ l of the generation amplicons (pool 1 and 2) into the new tubes.

21.1 Open tubes strip-by-strip and only keep tubes you are immediately using open.

22 Mix for 1 minute at 1400 rpm on the plate mixer (set to 4 °C) and spin down.

23 Load into the thermocycler and run the following program.

	Temperature	Time
	55 °C	3 minutes
	10 °C	hold

24 Once the program has finished, **immediately** remove tubes, spin down, and keep on a cold block.



- 25 Add 10  $\mu$ l ST2 to each tube with a multichannel pipette and mix for 1 minute on the plate mixer.
- 26 Incubate at RT for 5 minutes and spin down.
- 27 Place on magnetic stand and wait until liquid is clear (a few minutes).
- 28 Remove and discard all supernatant from each well with a multichannel pipette.
- 29 Remove from the magnetic stand and add 100  $\mu$ l TWB to each tube.
- 29.1 Mix for 1 minute at 1400 rpm on the plate mixer and spin down. Be careful not to introduce bubbles.
- 30 Repeat steps 27-29, but do not remove the second wash.

## Amplify Tagmented Ampicons

- 31 Prepare the following amplification master mix:

	Reagent Name	Volume per Sample
	EPM	24
	Nuclease Free Water	24

- 32 Place the tubes with tagmented ampicons on the magnetic stand and wait until the liquid is clear (a few minutes).
- 33 Working in sets of 3 strips at a time, remove all TWB from each tube.
- 33.1 Remove any residual TWB from tubes with lower volume pipette.

- 34 Remove the tubes from the magnetic stand and add 40  $\mu$ l of master mix to each tube.
- 35 Clean the foil of the index adapter plate with DNA-Away and 70% ethanol.
- 36 Add 10  $\mu$ l index adapters to each tube using a multichannel pipette.
- 36.1 Make note of which index adapter set was used.
- 37 Mix for 1 minute at 1400 rpm on the plate mixer and spin down.
- 38 Load into the thermocycler and run the following program:

	Temperature	Time	Cycles
	72 °C	3 minutes	
	98 °C	3 minutes	
	98 °C	20 seconds	7 cycles
	60 °C	30 seconds	(included in above 7)
	72 °C	1 minute	(included in above 7)
	72 °C	3 minutes	
	10 °C	hold	

## Pool and Clean Up

- 39 Spin down tubes, place on magnetic stand, and wait until liquid is clear.
- 40 Transfer 5  $\mu$ l of each library to new tubes.

	Strips Pooled	Volume per Sample
	1-3	15 $\mu$ l
	4-5	10 $\mu$ l
	6-12	5 $\mu$ l

40.1 Change tips in between tubes; vortex and spin down when finished.

41 Transfer 55  $\mu$ l (total volume - 5  $\mu$ l) from each tube to a new 1.5 mL tube

42 Add 0.9 x total pooled volume  $\mu$ l to the tube and vortex.

43 Incubate at RT for 5 minutes

44 Spin down and place on magnetic stand and wait until liquid is clear.

45 Remove and discard supernatant; do not remove tube from magnetic rack.

46 Add 1000  $\mu$ l of freshly prepared 80% ethanol and incubate for 30 seconds.

47 Repeat steps 45-46.

48 Remove the ethanol wash.

49 Add RSB and vortex to mix. Consult the following table for how much RSB to add:

	Strips Pooled	Volume per Sample
	1-2	30 $\mu$ l
	3-4	40 $\mu$ l



	Strips Pooled	Volume per Sample
	5-12	55 $\mu$ l

50 Incubate at RT for 2 minutes.

51 Transfer total volume - 5  $\mu$ l of the pooled libraries to a new 1.5 mL tube.

51.1 **This is the final pooled library tube.**

52 Quantify the library using the Qubit.

53 Analyze the fragment distribution using the bioanalyzer.

## Protocol references

JCV ecology and evolution. forthcoming.

Bioinformatics pipeline: [https://github.com/grubaughlab/DENV\\_pipeline](https://github.com/grubaughlab/DENV_pipeline)