



Dec 09, 2025

Jess western blot automated system

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g7nq11lwz/v1>

Ellen Hertz¹

¹National Institutes of Health



Ellen Hertz

National Institutes of Health

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7nq11lwz/v1>

Protocol Citation: Ellen Hertz 2025. Jess western blot automated system. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.q26g7nq11lwz/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 09, 2025

Last Modified: December 09, 2025

Protocol Integer ID: 234607

Keywords: ASAPCRN, jess western blot, samples for jess, system brief description of preparation, western blot, detailed instruction, automated system, sample, reagent pack, preparation, jess

Funders Acknowledgements:

ASAP

Grant ID: ASAP-000458

Abstract

Brief description of preparation of samples for Jess western blot automated system. Detailed instructions are included with the reagent pack.



Preparation of samples for Jess

- 1 There is a detailed instruction included in the reagent pack from Bio Techne. The protocol was applied without changes.
- 2 Open EZ Standard pack. Reconstitute DTT (40 uL dH₂O), Fluorescent 5X Master Mix (20 uL 10X Sample buffer and 20 uL DTT) and Biotinylated Ladder (20 ul dH₂O)
- 3 Dilute 10x sample buffer to 0.1x sample buffer in dH₂O
- 4 Combine one-part fluorescent master mix with four parts of sample (final concentration of 1 µg/µl)
- 5 Boil samples and biotinylated ladder at 95 °C for 5 min
- 6 Load Separation Capillary Cartridges with 3 ul of sample per well and antibodies
- 7 Spin plate at 1000 x g for 5 min
- 8 Load Separation Capillary Cartridges with wash buffer, luminol-peroxide mix and RePlex reagent (if running replex protocol)