



Oct 29, 2025

JMN-MSMP Mouse Muscle Fibrosis scRNA-seq

DOI

<https://dx.doi.org/10.17504/protocols.io.kqdg31n4zl25/v1>

Lexi Rindone¹

¹Johns Hopkins University



Lexi Rindone

Johns Hopkins University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kqdg31n4zl25/v1>

Protocol Citation: Lexi Rindone 2025. JMN-MSMP Mouse Muscle Fibrosis scRNA-seq. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.kqdg31n4zl25/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 28, 2025



Last Modified: October 29, 2025

Protocol Integer ID: 231020

Keywords: seq dataset of murine muscle fibrosi, murine muscle fibrosi, muscle fibrosi, enriched scrna, profiling of senotype, polycaprolactone, seq dataset, volumetric muscle loss injury, jmn

Funders Acknowledgements:

NIH

Grant ID: 1U54AG079779

Abstract

This protocol was used to generate a CD45+ cell-enriched scRNA-seq dataset of murine muscle fibrosis induced via a volumetric muscle loss injury and synthetic biomaterial (polycaprolactone) implant. This tissue is enriched in p16+ and p21+ cells, enabling the profiling of senotypes involved in muscle fibrosis.

Tissue Harvest and Cell Isolation

- 1 Euthanize mice according to approved IACUC protocol.
- 2 Harvest each quadricep from each animal. Take extra care to avoid cutting the femoral artery, as excess bleeding will contaminate the tissue with immune cells from the blood.
- 3 Place both quads in a petri dish containing chilled RPMI-1640 with L-Glutamine + 25 mM HEPES (Base Media). Keep petri dishes on ice.
- 4 Dice quads into <3 mm pieces using a razor blade.
- 5 Add 1:1 v/v of pre-warmed solution of Base Media + 1 mg/mL Liberase TL + 0.4 mg/mL DNase I (2X Digestion Media) to yield a 1X digestion solution.
- 6 Incubate quads at 37°C with gentle agitation for 45 min.
- 7 Add 2:1 v/v of chilled Base Media + 1% BSA to neutralize the Digestion Media.
- 8 Strain tissues through a 100 µm cell strainer into a 50 mL tube. Use the end of a syringe plunger to force cells through the strainer. Rinse the strainer with 1X DBPS twice, mashing the tissues in between each wash.
- 9 Centrifuge cells at 400 g for 5 min.
- 10 Resuspend cells in 5 mL of ACK Lysis Buffer at room temperature for 3 min.
- 11 Add >20 mL of 1X DBPS and strain cells through a 40 µm strainer.
- 12 Centrifuge cells at 400 g for 5 min.



Dead Cell Magnetic Bead Removal

- 13 Resuspend cells in Dead Cell Removal MicroBeads (Miltenyi # 130090101; 10^7 cells/100 μ L) for 15 minutes at room temperature.
- 14 Prepare LS column with 3 mL of 1X Binding Buffer.
- 15 Add 3 mL of 1X Binding Buffer to cells and add to LS column to enable flow-through of live cells.
- 16 Collect the live cells from the flow-through solution.
- 17 Wash column with 4×3 mL of 1X Binding Buffer.

CD45+ Cell Enrichment

- 18 Resuspend cells in 80 μ L of Wash Buffer (1X DBPS + 1% w/v BSA + 1 mM EDTA) per 10^7 cells.
- 19 Add 20 μ L of CD45 Microbeads (Miltenyi # 130045801) per 10^7 cells (yielding 100 μ L of solution per 10^7 cells) and incubate for 15 minutes at 4°C.
- 20 Wash cells by adding 2 mL of Wash Buffer per 10^7 cells.
- 21 Spin down cells at 300 g for 10 minutes.
- 22 Resuspend cells in 500 μ L of Wash Buffer for 10^8 cells.
- 23 Prime column 3×3 mL of Wash Buffer.
- 24 Add cell suspension to column.



- 25 Wash column 3×3 mL of Wash Buffer.
- 26 Collect flow-through cell solution. This solution contains the CD45⁻ cell fraction.
- 27 Remove column from magnet, add 5 mL of Wash Buffer to column, and plunge over a new tube. Collect flow-through cell solution. This is the CD45⁺ cell fraction.
- 28 Count cells in the CD45⁺ and CD45⁻ cell fractions.
- 29 Combine cell solutions to yield a 1:1 mixture of CD45⁺ and CD45⁻ cells.
- 30 Submit samples to sequencing core to perform 10X 3' Chromium scRNA-seq. (We use the Johns Hopkins Single Cell & Transcriptomics Core.)

Single Cell Sequencing

- 31 Prepare libraries using 10X 3' Chromium v3.1 chemistry protocols and reagents.
- 32 Sequence libraries using NovaSeq platform (Illumina) with S2 100 flow cell at a target depth of 100,000 unique reads/cell.