

Sep 27, 2023

Version 2

# JMN-MSMP Muscle Single Cell Isolation V.2

DO

dx.doi.org/10.17504/protocols.io.j8nlko8xxv5r/v2

ccherry<sup>1</sup>

<sup>1</sup>SenNet JMN-MSMP

Cellular Senescence Net...



ccherry

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account





DOI: https://dx.doi.org/10.17504/protocols.io.j8nlko8xxv5r/v2

Protocol Citation: ccherry 2023. JMN-MSMP Muscle Single Cell Isolation. protocols.io <a href="https://dx.doi.org/10.17504/protocols.io.j8nlko8xxv5r/v2">https://dx.doi.org/10.17504/protocols.io.j8nlko8xxv5r/v2</a> Version created by <a href="mailto:ccherry">ccherry</a>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: September 27, 2023

Last Modified: September 27, 2023

**Protocol Integer ID:** 88425

Keywords: SenNet JHM-MSMP, msmp muscle single cell isolation single cell isolation, single cell, jmn, cell, isolation

#### **Abstract**

Single cell isolation from muscle

## Troubleshooting



## Single cell isolation

- 1 Tissue was finely diced in RPMI 1640 (Gibco) with 1.6 Wunsh U/mL Liberase TL (Roche Diagnostics) and 0.2 mgl/mL deoxyribonuclease I (Roche Diagnostics).
- 2 Finely diced tissue was incubated at 37C for 45 minutes.
- 3 Digested tissues were ground subsequently through 100um and 70um strainers (Thermo Fisher Scientific).
- 4 Cells were centrifuged at 300g for 5 minutes and washed twice with 1x PBS