

Sep 13, 2019

Version 2

K/2 Ian / K-ET V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.7bxhipn>



Ian Probert¹

¹CNRS-Sorbonne Université, Station Biologique, Place G. Tessier 29680 Roscoff FRANCE

Roscoff Culture Collection



Daniel Vaultot

Station Biologique, Roscoff, France , Nanyang Technologica...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.7bxhipn>

Protocol Citation: Ian Probert 2019. K/2 Ian / K-ET. protocols.io <https://dx.doi.org/10.17504/protocols.io.7bxhipn>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 13, 2019

Last Modified: September 13, 2019

Protocol Integer ID: 27735

Keywords: oceanic ultraphytoplankton, media for the culture, medium

Abstract

Modified from K medium by Ian Probert

Reference

Keller, M.D., Selvin, R.C., Claus, W. & Guillard, R.R.L. 1987. Media for the culture of oceanic ultraphytoplankton. J. Phycol. 23:633–8.

Before start

Please refer to our general recommendations to grow cultures :

<https://www.protocols.io/private/A48906DC1374AD6281495CB86A8F092F>



- 1 ■ To 994 ml of seawater (pH 8.2, adjusted with NaOH) add:

	Quantity	Compound	Stock solution (sterile)	Final conc. in K medium
	0.5ml	NaNO ₃	48.9542 g/litre H ₂ O	288 µM
	0.5ml	NH ₄ Cl *	0.535g/litre H ₂ O	5µM
	0.5ml	KH ₂ PO ₄	4.8992g/litre H ₂ O	18µM
	0.5ml	FeEDTA solution	(see recipe below)	(see below)
	0.5ml	Trace metal solution	(see recipe below)	(see below)
	1.0ml	f/2 vitamin solution	(see recipe below)	(see below)

* optional

FeEDTA solution

- 2 ■ To 950ml distilled H₂O add:

	Quantity	Compound	Stock solution	Final conc. in K
--	----------	----------	----------------	------------------

				medium
	4.3g	(Na) FeE DTA	-	5.85 μ M

- Make up to 1 litre with milliQ H₂O, sterilize (filter 0.22 μ m) and store in fridge.

Trace metal solution

- 3
- To 950ml distilled H₂O add:

	Quantity	Compound	Stock solution	Final conc. in K medium
	37.2 2g	Na ₂ EDT A.2H 2O	-	50 μ M
	1.0m l	CuS O _{4.5} H ₂ O	2.49 7g/litre H ₂ O	0.00 5 μ M
	1.0m l	Na ₂ MoO 4.2H 2O	7.25 85g/ litre H ₂ O	0.01 5 μ M
	1.0m l	ZnS O _{4.7} H ₂ O	23.0 g/litre H ₂ O	0.00 4 μ M
	1.0m l	CoS O _{4.7} H ₂ O	14.0 55g/ litre H ₂ O	0.02 5 μ M
	1.0m l	MnC l _{2.4} H ₂ O	178.1 1g/litre H ₂ O	0.45 μ M
	1.0m l	H ₂ S eO ₃	1.29 g/litre H ₂ O	0.00 5 μ M



	1.0ml	NiCl ₂ ·6H ₂ O	1.49 g/litre H ₂ O	0.00314 μM
--	-------	--------------------------------------	-------------------------------	------------

- Make up to 1 litre with milliQ H₂O, sterilize (filter 0.22 μm) and store in fridge.

f/2 Vitamin solution

- 4
- To 950ml distilled H₂O add

	Quantity	Compound	Stock solution	Final conc. in K medium
	1.0ml	Vit. B ₁₂ (cyanocobalamin)	0.5g /litre H ₂ O	0.37 nM
	1.0ml	Biotin	5.0mg /litre H ₂ O	2.0nM
	100.0mg	Thiamine HCl	-	0.3 μM

- Make up to 1 litre with milliQ H₂O, filter sterilize into plastic vials and store in freezer.

Sterilization of medium

- 5
- Optional: Heat to 80°C for 2 hours and leave to cool – this should kill most organisms but should not chemically modify the medium too much
 - Filter sterilize through 0.22 μm filters (e.g. Millipore Steritop units) into sterile (autoclaved) polycarbonate bottles.

For K-ET

- 6
- Add 10-30 ml marine soil extract (ET) as detailed in the protocol below



Protocol



NAME

Soil extract for algal media

CREATED BY

Daniel Vaultot[Preview](#)

- 6.1 Sample 10 g of dry soil from a place that does not contain any pesticide nor pollutant and where you can be sure to be able to come back in order to always use the same soil. The soil must be dry because the boiling step will be longer if the soil is not dry.
- 6.2 Add the 10 g soil to 400ml of milliQ water.
- 6.3 Boil during 1 hour.
- 6.4 Filter on 0.2 μm or first through 0.8/0.4 μm and then through 0.2 μm in order not to clog the filter.
- 6.5 Aliquot and freeze.
- 6.6 Add 30 à 50ml of soil extract per 1 L of medium.