

Jun 12, 2019 Version 1

K/2 Ian / K-ET V.1

DOI

dx.doi.org/10.17504/protocols.io.s53eg8n



Ian Probert¹

¹CNRS-Sorbonne Université, Station Biologique, Place G. Tessier 29680 Roscoff FRANCE

Roscoff Culture Collection



Daniel Vaultot

Station Biologique, Roscoff, France , Nanyang Technologica...

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.s53eg8n

Protocol Citation: Ian Probert 2019. K/2 Ian / K-ET. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.s53eg8n>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 30, 2018

Last Modified: June 12, 2019

Protocol Integer ID: 15259



Abstract

Modified from K medium by Ian Probert

Reference

Keller, M.D., Selvin, R.C., Claus, W. & Guillard, R.R.L. 1987. Media for the culture of oceanic ultraphytoplankton. J. Phycol. 23:633–8.

Before start

Please refer to our general recommendations to grow cultures :

<https://www.protocols.io/private/A48906DC1374AD6281495CB86A8F092F>

- 1
- To 994 ml of seawater (pH 8.2, adjusted with NaOH) add:

Quantity	Compound	Stock solution (sterile)	Final conc. in K medium
0.5ml	NaNO ₃	48.95 g/litre H ₂ O	288μM
0.5ml	NH ₄ Cl*	0.535 g/litre H ₂ O	5μM
0.5ml	KH ₂ PO ₄	4.899 g/litre H ₂ O	18μM
0.5ml	FeEDTA solution	(see recipe below)	(see below)
0.5ml	Trace metal solution	(see recipe below)	(see below)
1.0ml	f/2 vitamin solution	(see recipe below)	(see below)

* optional

FeEDTA solution

2

- To 950ml distilled H₂O add:

Quantity	Compound	Stock solution	Final conc. in K medium
4.3g	(Na)FeEDTA	-	5.85μM

- Make up to 1 litre with milliQ H₂O, sterilize (filter 0.22μm) and store in fridge.

Trace metal solution

- 3 ■ To 950ml distilled H₂O add:

Quantity	Compound	Stock solution	Final conc. in K medium
37.22 g	Na ₂ EDTA.2 H ₂ O	-	50 μM
1.0ml	CuSO ₄ .5H ₂ O	2.497 g/litre H ₂ O	0.005 μM
1.0ml	Na ₂ MoO ₄ .2 H ₂ O	7.258 5g/litre H ₂ O	0.015 μM
1.0ml	ZnSO ₄ .7H ₂ O	23.0g /litre H ₂ O	0.004 μM
1.0ml	CoSO ₄ .7H ₂ O	14.05 5g/litre H ₂ O	0.025 μM
1.0ml	MnCl ₂ .4H ₂ O	178.11 g/litre H ₂ O	0.45 μM
1.0ml	H ₂ SeO ₃	1.29g /litre H ₂ O	0.005 μM
1.0ml	NiCl ₂ .6H ₂ O	1.49 g/litre H ₂ O	0.003 14 μM

- Make up to 1 litre with milliQ H₂O, sterilize (filter 0.22 μm) and store in fridge.

f/2 Vitamin solution

- 4 ■ To 950ml distilled H₂O add

Quantity	Compound	Stock solution	Final conc. in K medium
1.0ml	Vit. B ₁₂ (cyanocob)	0.5g/litre H ₂ O	0.37 nM



		alamin)		
1.0ml	Biotin	5.0mg/litre H ₂ O	2.0nM	
100.0mg	Thiamine HCl	-	0.3μM	

- Make up to 1 litre with milliQ H₂O, filter sterilize into plastic vials and store in freezer.

Sterilization of medium

- 5
 - Optional: *Heat to 80°C for 2 hours and leave to cool – this should kill most organisms but should not chemically modify the medium too much*
 - Filter sterilize through 0.22μm filters (e.g. Millipore Steritop units) into sterile (autoclaved) polycarbonate bottles.

For K-ET

- 6
 - Add 10–30 ml marine soil extract (ET)