

Jan 16, 2026

## KAPP-Sen TMC: 10x Xenium Prime with Segmentation Kit

DOI

<https://dx.doi.org/10.17504/protocols.io.n92ld1mr9l5b/v1>

Santhosh Sivajothi<sup>1</sup>, Emily Soja<sup>1</sup>, Allison McNeilly<sup>1</sup>, Grace Frohock<sup>1</sup>, Shruti Bhargava<sup>1</sup>, William F Flynn<sup>1</sup>, Elise T Courtois<sup>1</sup>

<sup>1</sup>Single Cell Biology Lab, The Jackson Laboratory, USA

Cellular Senescence Net...

KAPP-Sen TMC



Sergii Domanskyi

The Jackson Laboratory for Genomic Medicine

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.n92ld1mr9l5b/v1>

**Protocol Citation:** Santhosh Sivajothi, Emily Soja, Allison McNeilly, Grace Frohock, Shruti Bhargava, William F Flynn, Elise T Courtois 2026. KAPP-Sen TMC: 10x Xenium Prime with Segmentation Kit . **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.n92ld1mr9l5b/v1>

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**We use this protocol and it's working**

**Created:** January 15, 2026

**Last Modified:** January 16, 2026

**Protocol Integer ID:** 238659

**Keywords:** SenNet KAPP-Sen TMC, situ spatial transcriptomic analysis of tissue section, 10x genomics xenium platform, 10x xenium v1 with segmentation kit, integrated 10x genomics workflow, resolution spatial gene expression profiling, 10x genomics workflow, detailed molecular characterization of tissue architecture, spatial gene expression profiling, spatial transcriptomic analysis, 10x xenium v1, tissue architecture, tissue section, segmentation kit, tissue processing, sen tmc, enabling detailed molecular characterization, xenium segmentation kit for accurate cell segmentation, 10x xenium prime with segmentation kit, xenium segmentation kit, broad transcriptomic coverage, 10x xenium prime, 5k standard human gene panel, accurate cell segmentation, gene add, gene expression

**Funders Acknowledgements:**

**National Institute on Aging (NIA) KAPP-Sen Tissue Mapping Collaborative**

**Grant ID:** U54 AG075941

## Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://dx.doi.org/10.17504/protocols.io.n92ld1mr9l5b/v1) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://dx.doi.org/10.17504/protocols.io.n92ld1mr9l5b/v1), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

## Abstract

This protocol describes the preparation, staining, and in situ spatial transcriptomic analysis of tissue sections using an integrated 10x Genomics workflow. Standardized procedures were followed for tissue processing, slide preparation, and staining to ensure consistency and quality. The 10x Genomics Xenium platform was employed for high-resolution spatial gene expression profiling, enabling detailed molecular characterization of tissue architecture. The analysis utilized the Xenium Segmentation Kit for accurate cell segmentation, the 5k Standard Human Gene Panel for broad transcriptomic coverage, and an additional 100-gene add-on panel to enhance detection of specific targets of interest.



## Section preparation

- 1 Block handling, QC, and sectioning were performed according to **CG000758 Rev E**.
- 2 After sectioning, slides were placed in a slide drying rack and baked at 42°C for 3 hours, then placed in a desiccator at room temperature overnight.

## Xenium In Situ

- 3 The assay was performed following Deparaffinization & Crosslinking according to **CG000758 Rev E**, section 4, immediately followed by probe priming, polishing, hybridization, ligation, and amplification according to **CG000760 Rev C**, then cell segmentation staining and quenching according to the same protocol.
- 4 The Xenium Analyzer was utilized according to **CG000584 Rev J**.

