

Sep 28, 2023

KAPP-Sen TMC: Fixation of Cells and Nuclei for Chromium Fixed RNA Profiling

DOI

<https://dx.doi.org/10.17504/protocols.io.x54v9py5zg3e/v1>

Juliana Alcoforado Diniz¹, Dylan Baker¹, Jessica Garofalo¹, Paul Robson^{1,2,3}

¹The Jackson Laboratory for Genomic Medicine, Farmington, CT, USA;

²Department of Genetics and Genome Sciences, University of Connecticut School of Medicine, Farmington, CT, USA;

³Institute for Systems Genomics, University of Connecticut, Farmington, CT, USA

Cellular Senescence Net...

KAPP-Sen TMC



Sergii Domanskyi

The Jackson Laboratory for Genomic Medicine

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.x54v9py5zg3e/v1>

Protocol Citation: Juliana Alcoforado Diniz, Dylan Baker, Jessica Garofalo, Paul Robson 2023. KAPP-Sen TMC: Fixation of Cells and Nuclei for Chromium Fixed RNA Profiling. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.x54v9py5zg3e/v1>



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 29, 2023

Last Modified: September 28, 2023

Protocol Integer ID: 85646

Keywords: SenNet KAPP-Sen TMC, chromium fixed rna profiling cell, nuclei for chromium fixed rna profiling cell, 10x genomics protocol cg000478, rna, fixation of cell, sen tmc, seq, fixation

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **[protocols.io](#)** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **[protocols.io](#)**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

Cells are fixated prior to scRNA-seq according to 10X Genomics protocol CG000478.

Fixation of Cells & Nuclei

- 1 Centrifuge sample at 350 rcf for 5 min at 4°C.
- 2 Remove the supernatant without disturbing the pellet.
- 3 Add 1 ml Fixation Buffer to the sample pellet and pipette mix 5x.
- 4 Incubate for 1 h at room temperature (20°C) or for 16-24 h at 4°C for long term storage.
- 5 Centrifuge at 850 rcf for 5 min at room temperature.
- 6 Remove the supernatant without disturbing the pellet.
- 7 Add 1 ml chilled Quenching Buffer to the sample pellet and pipette mix 5x and keep on ice.
- 8 Determine cell concentration of the fixed sample using AO/PI (acridine orange/propidium iodide) Cell Viability Kit for Luna-FL automated cell counter.
- 9 Thaw Enhancer (**10x Genomics PN-2000482**) for 10 min at 65°C. Vortex and centrifuge briefly. Keep warm and verify no precipitate before use.
- 10 Add 1 volume pre-warmed Enhancer to fixed sample in Quenching Buffer. Pipette mix.
- 11 Store sample at 4°C for up to 1 week.
- 12 For long term storage (up to 6 months) also add 50% glycerol for a final concentration of 10% with Quenching buffer and store at -80°C.