

Sep 06, 2024

KAPP-Sen TMC: Placenta tissue blocks paraffin embedding

DOI

<https://dx.doi.org/10.17504/protocols.io.261ge51bjg47/v1>

Ramalakshmi Ramasamy¹, Anne Marchini², Tim Adams², Elaine Bechtel², Lesley Bechtold²

¹The Jackson Laboratory for Genomic Medicine, Farmington, CT, USA;

²The Jackson Laboratory, Bar Harbor, ME, USA

Cellular Senescence Net...

KAPP-Sen TMC



Sergii Domanskyi

The Jackson Laboratory for Genomic Medicine

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261ge51bjg47/v1>

Protocol Citation: Ramalakshmi Ramasamy, Anne Marchini, Tim Adams, Elaine Bechtel, Lesley Bechtold 2024. KAPP-Sen TMC: Placenta tissue blocks paraffin embedding. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.261ge51bjg47/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 06, 2024

Last Modified: September 06, 2024

Protocol Integer ID: 107030

Keywords: KAPP-Sen TMC, embedding placenta formalin, placenta tissue, placenta formalin, fixed tissue blocks paraffin, tissue blocks paraffin, jax bar harbor histology service, tissue, paraffin, sen tmc

Funders Acknowledgements:

National Institute on Aging (NIA) KAPP-Sen Tissue Mapping Collaborative

Grant ID: U54 AG075941

Abstract

Placenta formalin-fixed tissue blocks paraffin embedding at JAX Bar Harbor Histology Services.



- 1 Procedure addition between sample collection at Mayo Clinic and JAX Bar Harbor Histology Services processing:
 - 1.1 Samples were received from Mayo Clinic in 10% NBF on ice.
 - 1.2 Formalin was replaced with 70% Ethanol and left overnight at room temperature.
 - 1.3 Samples were shipped at room temperature to JAX Bar Harbor Histology Services for tissue processing and paraffin embedding (FFPE).
- 2 At JAX Bar Harbor Histology Services samples were embedded in paraffin (FFPE) according to: <https://dx.doi.org/10.17504/protocols.io.rm7vzj9z2lx1/v1>