

Aug 30, 2024

KAPP-Sen TMC: Whole human kidney preparation

DOI

dx.doi.org/10.17504/protocols.io.eq2lyw4kwvx9/v1

Juliana Alcoforado Diniz¹, Ramalakshmi Ramasamy¹, Alejandro Peña², Qian Wu³, Jillian L. Woodworth², Vladislav Bugay², Manpreet Semwal², Paul Robson^{1,4,5}

¹The Jackson Laboratory for Genomic Medicine, Farmington, CT, USA;

²University of Texas Health Science Center at San Antonio, San Antonio, TX, USA;

³Department of Pathology and Laboratory Medicine, University of Connecticut School of Medicine, Farmington, CT, USA;

⁴Department of Genetics and Genome Sciences, University of Connecticut School of Medicine, Farmington, CT, USA;

⁵Institute for Systems Genomics, University of Connecticut, Farmington, CT, USA

Cellular Senescence Net...

KAPP-Sen TMC



Sergii Domanskyi

The Jackson Laboratory for Genomic Medicine

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2lyw4kwvx9/v1>



Protocol Citation: Juliana Alcoforado Diniz, Ramalakshmi Ramasamy, Alejandro Peña, Qian Wu, Jillian L. Woodworth, Vladislav Bugay, Manpreet Semwal, Paul Robson 2024. KAPP-Sen TMC: Whole human kidney preparation. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.eq2lyw4kwvx9/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 27, 2024

Last Modified: August 30, 2024

Protocol Integer ID: 106495

Keywords: KAPP-Sen TMC, whole human kidney preparation this protocol, whole human kidney preparation, whole human kidney, sen tmc, dead donor, tmc, kapp, formalin preservation, preparation

Funders Acknowledgements:

National Institute on Aging (NIA) KAPP-Sen Tissue Mapping Collaborative

Grant ID: U54 AG075941

Abstract

This protocol describes removal, shipping, dissection, formalin preservation of a whole human kidney from a brain-dead donor.

Troubleshooting



- 1 Remove the whole kidney from a brain-dead donor. Flush the kidney with UW solution and perfuse with 1 liter of Waters IGL Pulsatile Perfusion Solution (PPS) for shipping.
- 2 Use overnight (ON) shipping, cold with ice.
- 3 Once received, bring the sample still on ice to the pathology department at UConn Health for dissection. The whole process was performed at room temperature (RT) and non-sterile.
- 4 Cut the whole organ in half along the middle plane (coronal section) and remove tissue pieces with a scalpel. The samples are roughly 1 cm x 1 cm and 0.5 cm thick and placed directly in labelled containers with 10 ml of 10% Neutral buffered formalin (NBF).
- 5 The samples are kept in RT and shipped the same day to JAX Histology Services in Bar Harbor. The samples are submitted to paraffinization (embedding of tissue blocks in paraffin) the same day they are received.