Jan 20, 2026

Kappaphycus alvarezii DNA Extraction by CTAB (2X) and precipitation in ammonium acetate and isopropanol

DOI

<https://dx.doi.org/10.17504/protocols.io.4r3l2z2rxl1y/v1>

Yasmin Cunha^{1,2}

¹BioBureau Biotecnologia; ²UFRJ

Yasmin Cunha

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS

DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l2z2rxl1y/v1>

Protocol Citation: Yasmin Cunha 2026. Kappaphycus alvarezii DNA Extraction by CTAB (2X) and precipitation in ammonium acetate and isopropanol . **protocols.io** <https://dx.doi.org/10.17504/protocols.io.4r3l2z2rxl1y/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 19, 2026

Last Modified: January 20, 2026

Protocol Integer ID: 238919

Keywords: Biologia molecular, DNA, Extração, tecido, kappaphycus alvarezii dna extraction by ctab, quality genomic dna extraction from kappaphycus alvarezii, kappaphycus alvarezii dna extraction, quality genomic dna extraction, kappaphycus alvarezii, red seaweed, ctab protocol, optimized ctab protocol, ctab

Abstract

Our optimized CTAB protocol offers a reliable solution for high-quality genomic DNA extraction from *Kappaphycus alvarezii*, a red seaweed, samples.

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Before start

Considerations for bench and solution prep:

Make sure the bench and materials have been thoroughly cleaned and sterilized. Use diluted chlorine (at 0.5%), followed by an ethanol 70% cleanup, and after it has dried, spray LGC Cleaner, let it soak for 5 minutes, and then cleanup with ethanol and a piece of paper.

Always prepare fresh solutions and aliquots, and use molecular biology-grade water for dilutions.

CTAB (2X)

25m

- 1 This buffer consists of Tris-HCl (at a final concentration of 0.1 M), NaCl (1 M), EDTA (0.2 M), and CTAB (at 2%), without PVP and β -mercaptoethanol (to be added immediately prior to use).
- 1.1 To prepare a 100 ml solution, add 10 ml of Tris-HCl (at 1 M, pH 8.0), 46.66 ml of NaCl solution (at 3 M) and 4 ml of EDTA (at 0.5 M).
- 1.2 Add 2 g of CTAB, and adjust the volume up to 100 ml with distilled water.
- 1.3 Autoclave the solution for 20 minutes, at 1 ATM and 121°C.
- 1.4 Before use of the solution, for 15 ml (one sample) add 0.3 g of PVP (for a 2% final concentration) and 30 μ l of β -mercaptoethanol (for 0.2%).

5m

20m

TE high salt

25m

- 2 This buffer consists of Tris-HCl (at a final concentration of 0.01 M), NaCl (1 M), and EDTA (1 mM).
- 2.1 To prepare a 50 ml solution, add 0.5 ml of Tris-HCl (at 1 M, pH 8.0), 16.66 ml of NaCl solution (at 3 M) and 100 μ l of EDTA (at 0.5 M). Adjust the volume up to 50 ml with distilled water.
- 2.2 Autoclave the solution for 20 minutes, at 1 ATM and 121°C.

5m

20m

TE high minus NaCl

- 3 This buffer consists of Tris-HCl (at a final concentration of 0.01 M), and EDTA (1 mM).
- 3.1 To prepare a 50 ml solution, add 0.5 ml of Tris-HCl (at 1 M, pH 8.0), 16.66 ml of NaCl solution (at 3 M) and 100 μ l of EDTA (at 0.5 M). Adjust the volume up to 50 ml with distilled water.
- 3.2 Autoclave the solution for 20 minutes, at 1 ATM and 121°C.

Tissue prep

10m

- 4 Use properly sterilized material, as described above.
- 5 If the samples are not being extracted immediately after sampling, place the algae pieces in aluminum foil, wrap and properly identify the sample. Immerse the sample in liquid nitrogen, and store it at -80 oC until processing.
- 6 Place small algae pieces in a crucible, and pour liquid nitrogen into it until the pieces are fully covered. Macerate the pieces into a fine powder, constantly replacing the liquid nitrogen, not letting it evaporate completely.

10m

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6h 20m

- 7 Transfer 5 g of the algae powder into a 50 ml falcon, with 15 ml of CTAB (with 2% PVP and 0.2% b-mercaptoethanol). Add 15 ml of chloroform:isoamyl alcohol (24:1). Gently invert the tubes to homogenize the content.
- 8 Centrifuge the falcons at maximum speed (5000 rpm) for 20 minutes, at room temperature.
- 9 After centrifugation, two phases will form. Transfer the aqueous phase into a new falcon and add 1 V (about 15 ml) of isopropanol.
- 10 Keep samples at -20 oC overnight or at -80 oC for 1 hour.
- 11 Centrifuge the falcons at maximum speed (5000 rpm) for 20 minutes, at room temperature. Discard the supernatant.
- 12 Add 5 ml of ethanol 80% and 0.08 V (400 ul) of ammonium acetate 7.5 M.
- 13 Centrifuge the falcons at maximum speed (5000 rpm) for 20 minutes, at room temperature. Discard the supernatant.
- 14 Add 5 ml of ethanol (100%). Incubate the tubes for 20 min at room temperature.

5m

20m

1h

20m

20m

- 15 Centrifuge the falcons at maximum speed (5000 rpm) for 20 minutes, at room temperature. Discard the supernatant. 
- 16 Resuspend the pellet in 5 ml in "TE high salt", incubate at 60 oC por 1-2 hours. 2h  
- 17 Add 5 ml of CTAB 2X (without PVP and b-mercaptoethanol) and 5 ml of chloroform:isoamyl alcohol (24:1).
- 18 Centrifuge the falcons at maximum speed (5000 rpm) for 20 minutes, at room temperature. Discard the supernatant. 20m 
- 19 After centrifugation, two phases will form. Transfer the aqueous phase into a new falcon and add 1 V (about 15 ml) of isopropanol.
- 20 Keep samples at -20 oC overnight or at -80 oC for 1 hour. 1h   
- 21 Centrifuge the falcons at maximum speed (5000 rpm) for 20 minutes, at room temperature. Discard the supernatant.
- 22 Add 5 ml of ethanol 80% and 0.08 V (400 ul) of ammonium acetate 7.5 M.
- 23 Centrifuge the falcons at maximum speed (5000 rpm) for 20 minutes, at room temperature. Discard the supernatant. 20m 
- 24 Add 5 ml of ethanol (100%). Centrifuge the falcons at maximum speed (5000 rpm) for 20 minutes, at room temperature. Discard the supernatant and let the pellets dry at room temperature. 20m 
- 25 Resuspend the pellets in 250 ul of "TE minus NaCl".
- 26 For RNA removal, add 0.25 ul of RNase A and incubate at 37 oC for 15 min 15m  

- 27 Assess sample quality and DNA concentration by NanoDrop and Qubit, and integrity by eletrophoresis by agarose gel (0.8%).

Acknowledgements

This procedure was originally created by *Milica Markovic_*.

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