

Feb 24, 2025

Kinetic assay with waveRAPID grating-coupled interferometry waveCreoptix system

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlk87xwl5r/v1>

Eda Capkin^{1,2,3}

¹Diamond Light Source; ²Research Complex at Harwell; ³ASAP Discovery Consortium

ASAP Discovery

OpenBind Consortium



Eda Capkin

--

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlk87xwl5r/v1>

Protocol Citation: Eda Capkin 2025. Kinetic assay with waveRAPID grating-coupled interferometry waveCreoptix system. protocols.io <https://dx.doi.org/10.17504/protocols.io.j8nlk87xwl5r/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 22, 2024

Last Modified: February 24, 2025

Protocol Integer ID: 110525

Keywords: Molecular interactions, Label-free detection, Grated-coupled interferometry, Automated sampling, coupled interferometry wavecreoptix system, interferometry wavecreoptix system this protocol, binding assay, creoptix waverapid grating, molecular interaction measurement, throughput molecular interaction measurement, waverapid grating, interferometry system, coupled interferometry system

Funders Acknowledgements:

National Institutes of Health/National Institute Of Allergy and Infectious Diseases (NIH/NIAID)

Grant ID: Grant ID: U19AI171399

Disclaimer

The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Institutes of Health.

Abstract

This protocol outlines the setup of kinetic binding assays using the Creoptix waveRAPID grating-coupled interferometry system. The method covers automated sample handling in 96- or 384-well formats, including DMSO corrections, control samples, and buffer blanks, enabling high-throughput molecular interaction measurements.

Materials

10x HBS-P:  HBS-P Buffer 10x **Cytiva Catalog #BR100671**

DMSO:  Pierce™ Dimethylsulfoxide (DMSO), Sequencing Grade **Thermo Scientific Catalog #10127403**

Echo 384 well plate:

 Echo® Qualified 384-Well Polypropylene Microplate 2.0 **Beckman Coulter Catalog #001-14555**

Creoptix 384 well plate:  240 µl 384 Deepwell Plate, Square Wells **StarLab Catalog #E2384-2400**

Creoptix 96 well plate:

 PlateOne 96 well Conical Bottom (V) Plates, max. Working Volume 320 µl **StarLab Catalog #S1833-9600**

Safety warnings

- ❗ Do not set the Flow Cell Temperature below 10 and/or 15°C. Lower temperature causes to the condensation in the chip.

Before start

Once the concentration value is specified in the method, it can not be changed in the evaluation part.

RAPID kinetics assay

5m

- 1 Variables that are for RAPID kinetic assay steps indicated the table.

A	B
Variable	Abbreviation for variable
Running buffer+DMSO %	[X]
Sample concentration	[Y]
Association and dissociation duration	[Z]
Flow rate	[F]

These specifics will differ for each protein target.

- 2 Build the method in the **Measure tab**/ RAPID kinetics.

- 2.1 In the Wizard choose the correct sample container type for the experiment: 48 vials, 96 or 384 well plate.
- 2.2 Enter the details of samples including sample name, concentration, and molecular weight.
- 2.3 Enter the configuration details for the method. Specify flow rate **[F]** and association and dissociation duration **[Z]**.
- 2.4 In the Reagent Rack section, choose 10 startup cycles.
- 2.5 Select Blanks: 1 blank after every 5 samples.
- 2.6 Choose DMSO correction (0.5%DMSO added to running buffer) as DMSO calibration solution at the beginning, end + every 10 th sample.



2.7 Enter a control sample and indicate its molecular weight and repeats depending on the sample size.
(i.e. For 96 samples, choose 1 sample injection after 50 th sample.)

2.8 Select DMSO wash: 1 DMSO wash (**[M]** 50 % volume DMSO) after every 5 samples.

2.9 Click create series.

3 Check the required buffer volume from cycles tab within the configuration. Prepare running buffer including DMSO **[X]** about 200 ml more than calculated buffer volume.

3.1 Degas the running buffer with DMSO% for 2 minutes.

2m

4 **Sample preparation**

4.1 Transfer samples from **[M]** 100 millimolar (mM) stocks 1536LDV to 96 or 384 well plate and add running buffer +DMSO% **[X]** to make final μ M sample concentration **[Y]**.

5 **Blank, DMSO calibration solution, and Control preparation**

5.1 Go to the Autosampler section within the configuration and prepare other solution (running buffer as blank and start-ups, DMSO calibration solution and control sample).

5.2 Place the samples according to the series in autosampler.

6 Run 2x Priming and Chip Prime steps in **Device**/Chip Prime

7 Go to the **Series**/Run Measurement.



Protocol references

Kartal Ö, Andres F, Lai MP, Nehme R, Cottier K. waveRAPID-A Robust Assay for High-Throughput Kinetic Screens with the Creoptix WAVEsystem. *SLAS Discov.* 2021 Sep;26(8):995-1003. doi: 10.1177/24725552211013827. Epub 2021 May 28. PMID: 34049465.