

Jul 06, 2023

Kompetitive Allele Specific PCR (KASP) with BioRad Software

DOI

<https://dx.doi.org/10.17504/protocols.io.dm6gpj9njgzp/v1>

Olivia E Todd¹, Eric L Patterson²

¹USDA-ARS Fort Collins, CO; ²Colorado State University, Fort Collins, CO



Olivia Todd

USDA-ARS

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gpj9njgzp/v1>

Protocol Citation: Olivia E Todd, Eric L Patterson 2023. Kompetitive Allele Specific PCR (KASP) with BioRad Software. protocols.io <https://dx.doi.org/10.17504/protocols.io.dm6gpj9njgzp/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 19, 2022



Last Modified: November 07, 2023

Protocol Integer ID: 74213

Keywords: KASP, Primer design, Hex/Fam, PCR, kompetitive allele specific pcr, lgc genomics kasp master mix, using lgc genomics kasp master mix, biorad software, biorad analyzation software, primer design with hex, basic kasp protocol, pcr, kasp, primer design, primer, allele

Abstract

A short guide to primer design with HEX/FAM tags and a basic KASP protocol using LGC Genomics KASP master mix and BioRad analyzation software.

Attachments



[590-1238.docx](#)

38KB

Materials

Materials

- KASP master mix
- primer mix

Primer Design

- 1 Design regular primers over the mutation you wish to test.

Note

The primer sequence must include this mutation to be specific.

- 2 Copy and paste the FAM and HEX tags to your wild type and mutant alleles. Pay attention to which tag is which.

FAM tag: GAAGGTGACCAAGTTCATGCT

HEX tag: GAAGGTCGGAGTCAACGGATT

Example: IAA16 GG to RR mutation Primers in *Bassia scoparia*

- 3

Note

Order these sequences through your primer manager.

	A	B
	Bs_IAA16_S_FP(He x)	GAAGGTCGGAGTCAACGGATT TGTTCTTCAGGACACAAGTTGTAGG
	Bs_IAA16_R_FP(Fa m)	GAAGGTGACCAAGTTCATGCT TGTTCTTCAGGACACAAGTTGTAAA
	Bs_IAA16_RP	AGTTTGATCATCGGACGTCTTCTT

Remove 2x KASP master mix from freezer, place  On ice (LGC Genomics, Beverly, MA, USA).

- This volume of  432 µL contains KASP enzyme, buffer, cofactors, dNTPs.

- 4 Thaw primers; one reverse primer, two forward primers, each specific to a SNP and corresponding fluorophore.
- 5 Make primer mix:



	A	B
	EACH forward primer	18 µl
	Reverse primer	45 µl
	H2O	69 µl
	Total primer mix	150

6 Add  12 µL of primer mix to the 2x KASP tube, this now the KASP master mix.



7 Add  4 µL of KASP master mix to each sample well of a 96 well plate.



8 Add  4 µL of template DNA @  5 ng/µl -  20 ng/µl .



9 Include 3 NTCs, as well as appropriate controls.

Bio-Rad machine and software

10 Set up thermocycling conditions according to the protocol:




	A	B	C
Step		Temperature	Length
1		94 C	15 minutes
2		94 C	20 seconds
3		61 C	1 min
4		Go to step 2	10x
5		94 C	20 seconds
6		55 C	1 min
7		30 C	10 seconds
8		Go to 5	35x (take read each cycle)

- 11 Make sure that each well of the plate has both FAM and HEX fluorophores selected.
- 12 Click Start Run.

Analysis

- 13 On the “Allelic Discrimination” tab the relative fluorescence units (RFU) for FAM and HEX are displayed, these are used to make a call on the SNP(s) present in each sample.
- 14 This can be compared to the controls for identification of genotype or species.
Example output of allelic discrimination:



Y axis: wild type, X axis: mutant