

May 24, 2022

Labeled microtubules for single-molecule imaging

DOI

dx.doi.org/10.17504/protocols.io.bp2l6bdedgqe/v1

John Salogiannis¹

¹Department of Cellular and Molecular Medicine, University of California, San Diego, La Jolla, CA 92093



Mariusz Matyszewski

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bp2l6bdedgqe/v1

Protocol Citation: John Salogiannis 2022. Labeled microtubules for single-molecule imaging. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bp2l6bdedgqe/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 20, 2021



Last Modified: June 01, 2024

Protocol Integer ID: 46496

Keywords: single-molecule, microtubule, imaging, ASAPCRN, labeled microtubule, stabilized microtubule, molecule imaging assay, molecule imaging this protocol, molecule imaging, biotin slide, labeled taxol, imaging

Abstract

This protocol is for making labeled taxol-stabilized microtubules to be used for single-molecule imaging assays by adhering to biotin slides.

Materials

Required Buffers:

5x BRB80:

- [M] 400 millimolar (mM) PIPES, pH 6.8 with KOH
- [M] 10 millimolar (mM) MgCl₂
- [M] 5 millimolar (mM) EGTA

2x Polymerization Mix:

- [M] 2 x BRB80
- [M] 2 millimolar (mM) DTT
- [M] 2 millimolar (mM) GTP (Add last)
- [M] 2 millimolar (mM) MgCl₂ (Add second to last)
- [M] 20 % DMSO
- Mix well between adding each ingredient

Safety warnings

! For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet).



Before start
















Please take notice of the buffer preparation in the Materials section



Create taxol-stabilized labeled microtubules (can be reused for multiple weeks):

50m

- 1 In a prechilled 1.5 mL Eppendorf tube, make a  10 μL mixture of  10 mg/mL tubulin .

The mixture should be 80% unlabeled, 10% biotin-tubulin, and 10% 405-tubulin. Mix by gently flicking.
- 2 Let it sit  On ice for  00:10:00 .
- 3 Add equal volume of  2 x polymerization buffer ( 10 μL). Mix by gently flicking.
- 4 Incubate at  37 $^{\circ}\text{C}$ for  00:30:00 . Make a  1 x BRB80 +  1 millimolar (mM) DTT +  20 micromolar (μM) Taxol stock and incubate at  37 $^{\circ}\text{C}$ at the same time.
- 5 Add equal volume of prewarmed  1 x BRB80 + DTT + Taxol ( 20 μL). Mix by gently flicking.
- 6 Incubate at  37 $^{\circ}\text{C}$ for at least  00:10:00 (solution will be stable for hours at this point).
- 7 Store in the dark at  Room temperature . Should be usable for several weeks, but more aggregates will appear over time.



10m



30m



10m

