


Dec 12, 2019

## Lactophenol-aniline blue solution

DOI

[dx.doi.org/10.17504/protocols.io.baf4ibqw](https://dx.doi.org/10.17504/protocols.io.baf4ibqw)



Andrea Sweigart<sup>1</sup>

<sup>1</sup>University of Georgia

Mimulus



Makenzie Whitener

UGA

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.baf4ibqw](https://dx.doi.org/10.17504/protocols.io.baf4ibqw)

**Protocol Citation:** Andrea Sweigart 2019. Lactophenol-aniline blue solution. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.baf4ibqw>

**Manuscript citation:**

Kearns, CA and Inouye, DW. Techniques for Pollination Biologists. University Press of Colorado. 1993. pp. 111–111.

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** December 12, 2019

**Last Modified:** December 13, 2019

**Protocol Integer ID:** 30940

**Keywords:** pollen, staining, lactophenol, blue

## Abstract

Lactophenol-aniline blue recipe

## Guidelines

Wear gloves and prepare in the hood.

## Materials

### MATERIALS

-  Water
-  Phenol
-  Glycerine/Glycerol
-  Lactic Acid (85% v/v)
-  1% aqueous aniline blue

### STEP MATERIALS

-  Water
-  Lactic Acid (85% v/v)
-  Phenol
-  Glycerine/Glycerol
-  1% aqueous aniline blue

## Protocol materials

☒ Lactic Acid (85% v/v)

☒ Glycerine/Glycerol

☒ 1% aqueous aniline blue

☒ Phenol

☒ Glycerine/Glycerol

☒ 1% aqueous aniline blue

☒ Water

☒ Phenol

☒ Water

☒ Lactic Acid (85% v/v)

☒ Glycerine/Glycerol

☒ Water

☒ Lactic Acid (85% v/v)

☒ Phenol

☒ 1% aqueous aniline blue

## Safety warnings

❗ Please refer to the MSDS for further safety details

## Make Lactophenol

1

To make  100 mL of stock Lactophenol:

Combine:

 20 mL

 Phenol

 40 mL

 Glycerine/Glycerol

 20 mL

 Water

 20 mL

 Lactic Acid (85% v/v)

## Lactophenol-aniline blue recipe

2

To prepare  100 mL working solution, combine:

 100 mL stock Lactophenol

 1 mL to  5 mL

 1% aqueous aniline blue