



Apr 18, 2023

## Lakes ABPS Protocol - Optimized protocol for the extraction of fish DNA from freshwater sediments (Thomson-Laing et al., 2022)

DOI

[dx.doi.org/10.17504/protocols.io.yxmvm24d6g3p/v1](https://dx.doi.org/10.17504/protocols.io.yxmvm24d6g3p/v1)

Georgia Thomson-Laing<sup>1</sup>

<sup>1</sup>Cawthron Institute, Nelson, New Zealand



Grayson Huston

University of Maine

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvm24d6g3p/v1>

**Protocol Citation:** Georgia Thomson-Laing 2023. Lakes ABPS Protocol - Optimized protocol for the extraction of fish DNA from freshwater sediments (Thomson-Laing et al., 2022). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.yxmvm24d6g3p/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**We use this protocol and it's working**

**Created:** April 17, 2023

**Last Modified:** April 19, 2023

**Protocol Integer ID:** 80636

**Keywords:** lakes abps protocol, lake sediment samples by an alkaline lysis method, fish dna from freshwater sediment, fish seddna from surface sediment, comparison of multiple seddna extraction method, lake sediment, range of lake sediment, multiple seddna extraction method, freshwater sediment, detecting fish seddna, protocol for the extraction, fish dna, sediment sample, surface sediment, abps protocol, optimized abps protocol, alkaline buffer, highest concentrations of target genet, extraction, alkaline lysis method, lake, sediment, multiple powersoil extraction, dna, abp

## Abstract

DNA was extracted from lake sediment samples by an alkaline lysis method with ethanol precipitation adapted from method described by Kuwae et al. (2020); Sakata et al. (2020b); Sakata et al. (2020a).

In a comparison of multiple sedDNA extraction methods, the ABPS (Alkaline buffer - power soil) protocol yielded the highest concentrations of target genets across a range of lake sediments. This protocol was further optimized (65C incubation temperature, pooling of multiple PowerSoil extractions) to overcome technical challenges related to co-precipitation of organic content in lake-surface sediments.

The optimized ABPS protocol is called the "Lakes ABPS protocol"










This protocol has proven to be successful at detecting fish sedDNA from surface sediments in multiple systems for multiple species.

## Troubleshooting



## Alkaline extraction

1h 51m









- 1 **INTO** a 50 mL tube add:
  -  10 g of sediment sample
  -  6 mL sodium hydroxide (0.33M)
  -  3 mL Tris-EDTA (pH 8)
- 2 **VORTEX** for  00:01:00
  - INCUBATE** at  65 °C for  00:50:00
- 3 **ALLOW** samples to cool to  Room temperature
  - CENTRIFUGE** at  15000 x g for  01:00:00

51m

1h

## Ethanol precipitation

2h

- 4 **TRANSFER**  7.5 mL of supernatant to a new 50 mL tube
  - ADD**  7.5 mL of Tris HCl (1M, pH 6.7) to neutralize
- 5 **ADD**  1.5 mL sodium acetate (3M, pH 5.2)
  - ADD**  30 mL of molecular grade 100% ethanol
  - INCUBATE** samples at  -20 °C for <  01:00:00
- 6 **CENTRIFUGE** samples at  10000 x g for  01:00:00
  - DISCARD** supernatant
  - RETAIN** precipitated pellet

1h

1h

## DNeasy PowerSoil extraction

- 7 **EXTRACT** the total pellet using multiple DNeasy PowerSoil DNA Isolation Kit extractions following the manufacturer's instructions



#### Note

🧪 0.25–0.5 g of pellet per extraction

**POOL** resultant DNA elutes

**DNA is now ready for downstream applications**

## Protocol references

Thomson-Laing, G., Howarth, J. D., Vandergoes, M. J., & Wood, S. A. (2022). Optimised protocol for the extraction of fish DNA from freshwater sediments. *Freshwater Biology*, 67, 1584–1603. <https://doi.org/10.1111/fwb.13962>