

Oct 21, 2019

## LAMP assay

DOI

[dx.doi.org/10.17504/protocols.io.8jahuie](https://dx.doi.org/10.17504/protocols.io.8jahuie)

Niels Appelman<sup>1</sup>

<sup>1</sup>Wageningen University

iGEM Wageningen 2019



Niels Appelman

Wageningen University

---

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.8jahuie](https://dx.doi.org/10.17504/protocols.io.8jahuie)

**Protocol Citation:** Niels Appelman 2019. LAMP assay. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.8jahuie>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 21, 2019

**Last Modified:** October 21, 2019

**Protocol Integer ID:** 28994

## Abstract

This is a protocol used to perform LAMP assays. It was adapted to

<https://international.neb.com/protocols/2016/08/15/warmstart-colorimetric-lamp-2x-master-mix-typical-lamp-protocol-m1800>

- 1 Prepare a 10X primermix according to the following scheme:

Prime r	Conc entrat ion
FIP	16 uM
BIP	16 uM
F3	2 uM
B3	2 uM
LOOP F	8 uM
LOOP B	8 uM

- 2 Add the following components to a microcentrifuge tube

Comp onent	Volu me adde d (25 uL reacti on)	Volu me adde d (10 uL reacti on)
Mast ermix	12.5 uL	5 uL
Lamp prime rs (10X)	2.5 uL	1 uL
Temp late	x uL	x uL
MilliQ	10 – x uL	4 – x uL

Add a no-template control to ensure template-specificity

- 3 Incubate at 65 C for 30 minutes.
- 4 Record the color of all samples. Yellowing of a sample indicates a positive result.



- 5 In case of doubt regarding the result of a sample, run it on an agarose gel. The presence of many bands, often not individually recognizable, indicates a positive result.