

Dec 11, 2020

Library Adapter Preparation for Dual-Index Double Stranded DNA Illumina Sequencing

 In 2 collections

DOI

<https://dx.doi.org/10.17504/protocols.io.bem5jc86>

Franziska Aron¹, Guido Brandt²

¹Friedrich-Schiller Universität Jena;

²Department of Archaeogenetics, Max Planck Institute for the Science of Human History


WarinnerGroup

MPI EVA Archaeogenetics



Franziska Aron

Friedrich-Schiller Universität Jena, Max Planck Institute fo...



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bem5jc86>

External link: <https://doi.org/10.1371/journal.pone.0241883>



Protocol Citation: Franziska Aron, Guido Brandt 2020. Library Adapter Preparation for Dual-Index Double Stranded DNA Illumina Sequencing. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bem5jc86>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 06, 2020

Last Modified: December 11, 2020

Protocol Integer ID: 35229

Keywords: ancient DNA, palaeogenetics, archaeogenetics, DNA, Illumina, aDNA, nucleic acids, paleogenetics, archeogenetics, library preparation, genomic DNA, genomics, ancient dna library preparation for illumina, ancient dna libraries for illumina next, stranded ancient dna library preparation, sequencing adapter preparation protocol, stranded ancient dna library, based ancient dna library, ancient dna library, library adapter preparation, library adapter preparation for dual, index double stranded dna, adapter preparation protocol, stranded dna, sequencing, adapter mix, production of an adapter mix, dna, library preparation, library, illumina next

Abstract

Adapter preparation protocol for double-stranded ancient DNA libraries for Illumina Next-Generation-Sequencing, based on Meyer and Kircher *et al.* (2010) *Cold Spring Harb. Protoc.* (doi: [10.1101/pdb.prot5448](https://doi.org/10.1101/pdb.prot5448)).

The protocol describes the production of an adapter mix for 2000 libraries.

This protocol is used in conjunction with Meyer and Kircher-based **ancient DNA libraries** and is described accordingly (see library construction like: Non-UDG treated double-stranded ancient DNA library preparation for Illumina sequencing).

Image Attribution

Matthäus Rest



Guidelines

Working in an Ancient DNA Laboratory

- All steps of the protocol should take place in a clean room facility specifically designed for ancient DNA.
- The researcher performing lab work should wear correspondingly suitable lab-wear, such as:
 - full-body suit with hood (e.g., Tyvek)
 - hairnet
 - face mask
 - two pairs of clean gloves
 - clean shoes
 - protective glasses
- Sample processing should be carried out in separated work benches with integrated UV irradiation (e.g. Dead Air PCR work bench)
- Surfaces and equipment should be regularly decontaminated with e.g. bleach solution or Thermofisher's DNA AWAY (or similar) and irradiated with UV.
- All home-made buffers should be prepared in a separate dedicated PCR-free ultra-clean room and UV-irradiated for 30 min.

Please see the following for more detailed guidance:

Llamas, B. et al., 2017. From the field to the laboratory: Controlling DNA contamination in human ancient DNA research in the high-throughput sequencing era. *STAR: Science & Technology of Archaeological Research*, 3(1), pp.1–14. Available at: <https://doi.org/10.1080/20548923.2016.1258824>.

All steps of the protocol should take place in a clean room facility specifically designed for ancient DNA.

Protocol Specific Guidelines

This protocol requires the use of a two rooms - a buffer preparation room and a library preparation room.

Materials

MATERIALS

✕ 2 ml LoBind Tubes **Eppendorf Catalog #0030108078**

✕ 0.2 ml PCR Tube strips **Eppendorf Catalog #0030124359**

✕ 1.5 ml LoBind tubes **Eppendorf Catalog ##0030108051**

✕ EDTA (0.5 M) pH 8.0 **Life Technologies Catalog #AM9261**

✕ 5 M Sodium chloride (NaCl) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S5150-1L**

✕ Tris-HCL **Life Technologies Catalog #15568025**

✕ Water HPLC Plus **Merck MilliporeSigma (Sigma-Aldrich) Catalog #34877-2.5L-M**

Lab Equipment

PCR Thermocycler (e.g. Eppendorf Mastercycler Nexus)

UV irradiation box or cross linker (e.g. Vilber Lourmat Bio-Link BLX-254)

Additional Reagents

Oligos (e.g. SigmaAldrich Custom DNA Oligos)

Oligo ID	Sequence (5'-3')	Concent ration
IS1_adapter.P5.F	A*C*A*C*TCTTTCCCTACACGACGCTCTTCCG*A *T*C*T	500 µM
IS2_adapter.P7.F	G*T*G*A*CTGGAGTTCAGACGTGTGCTCTTCCG* A*T*C*T	500 µM
IS3_adapter.P5+P7 .R	A*G*A*T*CGGAA*G*A*G*C	500 µM

* indicates a PTO bond

Safety warnings

! Reagents

Household bleach solution (2-6%) diluted to a working concentration of 0.2-0.5 % NaClO in total

- H290 May be corrosive to metals.
- H314 Causes severe skin burns and eye damage.
- H411 Toxic to aquatic life with long lasting effects.
- EUH206 Warning! Do not use together with other products. May release dangerous gases (chlorine). Remove from surface after recommended incubation time with water-soaked tissue.



DNA AWAY

- H314 Causes severe skin burns and eye damage.



Note: Both bleach solutions and DNA AWAY are used for decontamination. DNA AWAY is less corrosive than bleach and should be preferred for decontamination of sensitive equipments such as surfaces of electric devices.

EDTA

- H373 May cause damage to organs through prolonged or repeated exposure.



Sodium Chloride

- H290 May be corrosive to metal
- H314 Causes severe skin burns and eye damage
- H400 Very toxic to aquatic life



Equipment

UV radiation

- UV radiation can damage eyes and can be carcinogenic in contact with skin. Do not look directly at unshielded UV radiation. Do not expose unprotected skin to UV radiation.
- UV emitters generate ozone during operation. Use only in ventilated rooms.





Before start

Planning

This protocol takes around 3 hours including the cleaning process of the workspace afterwards.

Preparation of Reagents

Only the oligo hybridisation buffer can be prepared within buffer preparation room, with a DNA-free hood. All other steps are performed in the the library preparation room.

HPLC-Water should be decontaminated with a **30 min UV irradiation** before use.

Equipment

Make sure all necessary equipment is available (see Materials)

Abbreviations

EDTA = Ethylendiaminetetraacetic acid

HPLC = High Performance Liquid Chromatography (-Grade Water)

NaCl = Sodium chloride

Tris-HCl = Tris hydrochloride

UV = Ultraviolet (radiation)



Hybridization buffer preparation (Buffer Preparation Room)

- 1 Prepare oligo hybridization buffer (10× concentration,  1 mL for 50 reactions)

	Reagent	Stock concentration [M]	Final Concentration [M]	1× Volume [μl]
	NaCl	5	0.5	100
	Tris-HCl	1	0.01	10
	EDTA	0.5	0.001	2
	UV HPLC-water			888
	Total			1000

Note

The oligo hybridization buffer can be stored at room temperature for up to one year.

- 2 Irradiate the buffer with UV for  00:30:00 without the lid or with an open lid.

Adapter Preparation (Library Preparation Room)

- 3 **Prepare P5 adapter** ( 100 μL per reaction)

Use one tube of a 0.2 ml PCR strip to set up the P5 adapter mix.

	Reagent	Stock concentration	Final concentration	1× Volume [μl]
	Oligo Hybridization Buffer	10 ×	1 ×	10
	IS1_adapter.P5	500 μM	200 μM	40
	IS3_adapter.P5+P7	500 μM	200 μM	40



	UV HPLC-water		10
	Total		100

Gently pipette up and down to mix. Split the  100 μL reaction into two tubes of the 0.2ml PCR strip with  50 μL each .

4 **Prepare P7 adapter** (100 μL per reaction)

Use one tube of a 0.2 ml PCR strip to set up the P7 adapter mix.

	Reagent	Stock concentration	Final concentration	1 $\times$ Volume [μl]
	Oligo Hybridization Buffer	10 $\times$	1 $\times$	10
	IS2_adapter.P7	500 μM	200 μM	40
	IS3_adapter.P5+P7	500 μM	200 μM	40
	UV HPLC-water			10
	Total			100

Split the  100 μL reaction into two tubes of the 0.2ml PCR strip with  50 μL each .

5 Prepare ready-to-use adapter mix.

- 5.1 Combine one P5 adapter reaction with one P7 adapter reaction in one tube, mix thoroughly by flicking the tubes with a finger, spin down briefly.



Repeat with the second P5 and P7 adapter reaction.

**Note**

In the end you should have two 0.2 ml tubes with a total of $\text{100 } \mu\text{L}$ volume, each containing $\text{50 } \mu\text{L}$ of the P5 adapter and $\text{50 } \mu\text{L}$ of the P7 adapter.

5.2 Incubate both reactions in a thermocycler with a heated lid at $\text{95 } ^\circ\text{C}$ for

00:00:10 ,

followed by a cool down ramp from $\text{95 } ^\circ\text{C}$ to $\text{12 } ^\circ\text{C}$ at a rate of

$\text{0.1 } ^\circ\text{C per sec}$.

**Combine and Dilute**

6 Combine both reactions into a 2 ml tube to obtain a ready-to-use double-stranded library adapter mix with $\text{100 micromolar (}\mu\text{M) adapter (each)}$.



7 Add $\text{1800 } \mu\text{L}$ of UV HPLC-water to dilute the double-stranded library adapter mix to $\text{10 micromolar (}\mu\text{M)}$ and aliquot the dilution in $10 \times 1.5 \text{ ml LoBind tubes}$, each containing $\text{200 } \mu\text{L}$. Briefly vortex and spin down before freezing.

8 Store the adapter mix at $\text{-20 } ^\circ\text{C}$

Note

Each aliquot contains sufficient double-stranded library adapter mix for 200 reactions (overall 2000 libraries).

The adapter mix aliquot can be thawed and re-frozen with no detriment to the reagent quality.