


Apr 16, 2020

Lightsheet Tissue Intake - Photodocumentation and Tracking

DOI

dx.doi.org/10.17504/protocols.io.bef6jbre



Seth Currlin¹, Marda Jorgensen¹

¹University of Florida

Optical Clearing of Tissue

Human BioMolecular Atl...



Seth Currlin

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bef6jbre

Protocol Citation: Seth Currlin, Marda Jorgensen 2020. Lightsheet Tissue Intake - Photodocumentation and Tracking. [protocols.io](https://dx.doi.org/10.17504/protocols.io.bef6jbre) <https://dx.doi.org/10.17504/protocols.io.bef6jbre>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 31, 2020

Last Modified: April 16, 2020

Protocol Integer ID: 35038

Keywords: lightsheet microscopy, tissue processing, tissue tracking



Abstract

In order to track tissue position and orientation each sample is photographed when received for lightsheet microscopy (Figure 1, A). Large samples ($>1\text{ cm}^3$), such as the spleen in Figure 1, are cut into 2 mm sections (Figure 1, B) to facilitate tissue clearing, staining, and three-dimensional imaging. These tissues are batch processed and withdrawn from the lipid clearing pipeline as needed (Figure 1, C); at this time a novel identifier is assigned to the sample.

Guidelines

Photodocumentation, in our pipeline, occurs following tissue-hydrogel polymerization and multiple washes in PBS. This serves to remove trace fixatives for safer handling.

Materials

MATERIALS

✂ 50 ml conical tubes

✂ Forceps (tweezers), 12.5cm, Blunt End **Bio Basic Inc. Catalog #FC003.SIZE.1**

✂ Razor blades **Fisher Scientific Catalog #12-640**

✂ Hexagonal Polystyrene Weighing Dishes **Thermo Fisher Catalog #02202103**

Safety warnings

- ! Wash tissues of any fixatives, or other dangerous chemicals, prior to manipulating for photodocumentation.



- 1 Photograph whole tissue; capture from various angles if possible. A dissection ruler or grid should accompany each photograph.
- 2 Using a razor blade and forceps slice the tissue into sections approximately 2 mm wide. Try to slice in a continuous motion, while applying a gentle downward force.
- 3 Set aside tissue section and repeat Step 2 as needed. Maintain or keep track of slice order and orientation.
- 4 Photograph the sequential 2 mm slices, again using a dissection ruler or grid.
- 5 Retain all slices within one container for subsequent processing (lipid removal).