

Jun 12, 2025

## Lipidomic analysis in HeLa cells

DOI

[dx.doi.org/10.17504/protocols.io.eq2ly69yegx9/v1](https://dx.doi.org/10.17504/protocols.io.eq2ly69yegx9/v1)

Rosanne Wouters<sup>1,2</sup>, Igor Beletchi<sup>1</sup>, Chris Van den Haute<sup>3,4,2</sup>, Veerle Baekelandt<sup>3,2</sup>, Shaun Martin<sup>1</sup>,  
jan eggermont<sup>1</sup>, Peter Vangheluwe<sup>1,2</sup>

<sup>1</sup>Laboratory of Cellular Transport Systems, Department of Cellular and Molecular Medicine, KU Leuven;

<sup>2</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network;

<sup>3</sup>Research Group for Neurobiology and Gene Therapy, Department of Neurosciences, KU Leuven;

<sup>4</sup>Leuven Viral Vector Core, KU Leuven



Rania Abou El Asrar

KU Leuven

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2ly69yegx9/v1>

**Protocol Citation:** Rosanne Wouters, Igor Beletchi, Chris Van den Haute, Veerle Baekelandt, Shaun Martin, jan eggermont, Peter Vangheluwe 2025. Lipidomic analysis in HeLa cells . **protocols.io** <https://dx.doi.org/10.17504/protocols.io.eq2ly69yegx9/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** January 09, 2025

**Last Modified:** June 12, 2025

**Protocol Integer ID:** 117995

**Keywords:** lipidomic analysis in hela cell, lipidomic analysis, protocol for the lipidomic analysis, hela cell, cell

**Funders Acknowledgements:**

**Aligning Science Across Parkinson's**

Grant ID: ASAP-000458

**C1 KU Leuven grant**

Grant ID: C15/15/073

**Fonds voor Wetenschappelijk Onderzoek (FWO) Flanders**

Grant ID: S006617N SBO Neuro-Traffic

**Michael J. Fox Foundation**

Grant ID: MJFF-008610

## Abstract

A protocol for the lipidomic analysis in HeLa cells is described.

## Troubleshooting



## Lipidomic analysis in HeLa cells

- 1 Collect cells and homogenize cell pellets in 0.1 ml water.
- 2 To the homogenate, add lipid standards PC25:0, PC43:6, PI25:0, PI31:1, PI43:6, PS25:0, PS31:1, PS37:4 (Avanti Polar Lipids), on the amount of [protein or DNA] in the original sample.
- 3 Lipid extraction: 700  $\mu$ l of sample was mixed with 800  $\mu$ l 1 N HCl:CH<sub>3</sub>OH 1:8 (v/v), 900  $\mu$ l CHCl<sub>3</sub> and 200  $\mu$ g/ml of the antioxidant 2,6-di-tert-butyl-4-methylphenol (BHT; Sigma Aldrich).
- 4 Evaporate the organic fraction using a Savant Speedvac spd111v (Thermo Fisher Scientific) at room temperature.
- 5 Store the remaining lipid pellet at -20 °C under argon.
- 6 Just before mass spectrometry, reconstitute lipid pellets in running solution (CH<sub>3</sub>OH:CHCl<sub>3</sub>:NH<sub>4</sub>OH; 90:10:1.25; v/v/v)
- 7 Analyze **phospholipid** species by **electrospray ionization tandem mass spectrometry** (ESI-MS/MS) on a hybrid triple quadrupole/linear **ion trap** mass spectrometer (4000 QTRAP system; Applied Biosystems SCIEX) equipped with a TriVersa NanoMate (Advion Biosciences) robotic nanosource for automated sample injection and spraying.
- 8 Perform **phospholipid** profiling by (positive or negative) precursor ion or neutral loss scanning at a collision energy of 50 eV/45 eV, 35 eV, -35 eV, and -60 eV for precursor 184 (phosphatidylcholine (PC)), neutral loss 87 (phosphatidylserine (PS)) and precursor 241 (phosphatidylinositol (PI)), respectively.
- 9 Perform phospholipid quantification by **multiple reactions monitoring** (MRM), the transitions being based on the neutral losses or the typical product ions as described in step 8.
- 10 Relative quantification of phospholipids is obtained as ratio of analytes to their internal standard. The lipid data is normalized to nucleic acid concentration.