

Aug 26, 2023

Liposome binding

DOI

<https://dx.doi.org/10.17504/protocols.io.yxmvmndqng3p/v1>

Xinbo Wang^{1,2}, Pietro De Camilli^{1,2}

¹1. Departments of Neuroscience and of Cell Biology, Howard Hughes Medical Institute, Program in Cellular Neuroscience, Neurodegeneration and Repair, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

²2. Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815



xinbo.wang

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvmndqng3p/v1>

Protocol Citation: Xinbo Wang, Pietro De Camilli 2023. Liposome binding. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.yxmvmndqng3p/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 19, 2022

Last Modified: May 31, 2024

Protocol Integer ID: 68884

Keywords: Liposome binding, Co-sedimentation, Fluorescence microscopy, ASAPCRN, liposome, binding experiment, confocal fluorescence microscopy, microscopy, protocol details methods for Irrk2, Irrk2

Abstract

This protocol details methods for LRRK2-liposome binding experiments analyzed by co-sedimentation or confocal fluorescence microscopy, respectively.

Attachments



[iuudbv9p.docx](#)

19KB



Co-sedimentation analysis

- 1 Prepare samples in Beckman microfuge tubes with [M] 300 nanomolar (nM) LRRK2 or RCKW, [M] 20 micromolar (μ M) PS liposomes in the absence or presence of [M] 1 millimolar (mM) GMPPNP. 
- 2 Incubate samples at  37 °C for  00:30:00 . 

- 3 Centrifuge samples at  49000 rpm ( 100000 x g) for  00:20:00 in a Beckman Optima TLX ultracentrifuge. 

- 4 Collect Supernatants and Pellets.

Note
Note: Pellets were resuspended with the same volume of protein buffer as the supernatant.
- 5 Analyze samples by SDS-PAGE and coomassie blue staining.

Confocal fluorescence microscopy analysis

30m

- 6 For PS liposome binding. Prepare samples in PCR tubes with [M] 20 micromolar (μ M) Rhod-PE labeled PS liposomes and different concentrations of LRRK2 as indicated in the main text. 
- 6.1 For GC/PS nanotube binding. Prepare samples in PCR tubes with [M] 20 micromolar (μ M) Cy5-PE labeled GC/PS nanotubes and [M] 100 nanomolar (nM) GFP-LRRK2. 
- 7 Immediately deposit  6 μ L -  10 μ L samples of step 6 on  35 mm glass bottom dishes and incubate at  37 °C for  00:30:00 . 




Note

Note: Drop some buffer in the dish to prevent samples from drying out due to evaporation during incubation.

- 8 After incubation, capture the images with a spinning disk confocal (SDC) microscopy at 🌡️ Room temperature on a Nikon Ti-E inverted microscope using the Improvision UltraVIEW VoX system (Perkin-Elmer).

