



Aug 25, 2023

Liposome preparation

DOI

<https://dx.doi.org/10.17504/protocols.io.6qpvr612ovmk/v1>

Xinbo Wang^{1,2}, Pietro De Camilli^{1,2}

¹1. Departments of Neuroscience and of Cell Biology, Howard Hughes Medical Institute, Program in Cellular Neuroscience, Neurodegeneration and Repair, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

²2. Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815



xinbo.wang

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvr612ovmk/v1>

Protocol Citation: Xinbo Wang, Pietro De Camilli 2023. Liposome preparation. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.6qpvr612ovmk/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 19, 2022

Last Modified: May 31, 2024

Protocol Integer ID: 68882

Keywords: Liposome, LRRK2, Lipid nanotube, ASAPCRN, ps lipid nanotube, liposome, lrrk2 binding

Abstract

This protocol details methods for the preparation of 100% PS liposome and GC/PS lipid nanotubes used for LRRK2 binding, tubulation assays.

Attachments



[iuucbv9p.docx](#)

18KB

Materials

⊗ Brain PS **Avanti Polar Lipids, Inc. Catalog # 840032**

⊗ Galactosylceramide (GC) **Avanti Polar Lipids, Inc. Catalog #860546P**

⊗ Rhod-PE **Avanti Polar Lipids, Inc. Catalog #810150**

⊗ Cy5-PE **Avanti Polar Lipids, Inc. Catalog #810345**

Solutions to prepare:

Liposome buffer:

	A	B
	HEPES (7.4)	20 mM
	KCl	100 mM
	TCEP	0.5 mM



Liposome preparation

1h 5m

- 1 Dissolve lipid mixtures with chloroform in glass vials in moles percent as follows:
 - **PS liposomes:** 99.5% brain PS:0.5% Rhod-PE.
 - **GC/PS nanotubes:** 39.5% Galactosylceramide:60% brain PS:0.5% Cy5-PE.
- 2 Evaporate chloroform under a stream of nitrogen gas to produce a lipid film on the glass surface.
- 3 Dry the lipid film in a vacuum oven for  01:00:00 . 
- 4 Rehydrate the dried lipid films in liposome buffer at a final concentration of  1 mg/mL (~  1.2 millimolar (mM)). 
- 5 For PS mixtures, form the liposomes by three freeze (liquid N₂)–thaw ( 37 °C water bath) cycles.
- 5.1 For GC/PS mixtures, form the lipid nanotubes by a brief vortexing instead of freeze–thaw cycles.
- 6 Remove large aggregates by a brief centrifugation ( 500 x g for  00:05:00) and store in the dark at  4 °C to avoid photooxidation. 
