

Nov 12, 2023

Liquid Chromatography with Tandem Mass Spectrometry (LC-MS/MS)

DOI

dx.doi.org/10.17504/protocols.io.14egn3nrml5d/v1

Leonardo A Parra-Rivas¹

¹University of California, San Diego

Metabolomics Protocols & Workflows
Tech. support email: bbmisraccb@gmail.com



Leonardo A Parra-Rivas

University of California, San Diego

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn3nrml5d/v1>

Protocol Citation: Leonardo A Parra-Rivas 2023. Liquid Chromatography with Tandem Mass Spectrometry (LC-MS/MS). [protocols.io https://dx.doi.org/10.17504/protocols.io.14egn3nrml5d/v1](https://dx.doi.org/10.17504/protocols.io.14egn3nrml5d/v1)

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: November 12, 2023

Last Modified: November 12, 2023

Protocol Integer ID: 90827

Keywords: liquid chromatography with tandem mass spectrometry, tandem mass spectrometry, liquid chromatography, mass spectrometry, m, lc

Abstract

Liquid Chromatography with Tandem Mass Spectrometry (LC-MS/MS)

Troubleshooting



- 1 The data acquisition by LC-MS/MS was adapted from a published procedure

Citation

Seyfried NT, Dammer EB, Swarup V, Nandakumar D, Duong DM, Yin L, Deng Q, Nguyen T, Hales CM, Wingo T, Glass J, Gearing M, Thambisetty M, Troncoso JC, Geschwind DH, Lah JJ, Levey AI (2017)
. A Multi-network Approach Identifies Protein-Specific Co-expression in Asymptomatic and Symptomatic Alzheimer's Disease..

<https://doi.org/10.1016/j.cels.2016.11.006>

LINK

- 2 Derived peptides were resuspended in a loading buffer (0.1% trifluoroacetic acid, TFA) and separated on a Waters' Charged Surface Hybrid (CSH) column (150 µm internal diameter x 15 cm; particle size: 1.7 µm).
- 3 The samples were run on an EVOSEP liquid chromatography system using the 15 samples per day preset gradient (88 min) and were monitored on a Q-Exactive Plus Hybrid Quadrupole-Orbitrap Mass Spectrometer (Thermo Scientific).
- 4 The mass spectrometer cycle was programmed to collect one full MS scan followed by 20 data-dependent MS/MS scans. The MS scans (400-1600 m/z range, 3×10^6 AGC target, 100 ms maximum ion time) were collected at a resolution of 70,000 at m/z 200 in profile mode.
- 5 The HCD MS/MS spectra (1.6 m/z isolation width, 28% collision energy, 1×10^5 AGC target, 100 ms maximum ion time) were acquired at a resolution of 17,500 at m/z 200. Dynamic exclusion was set to



exclude previously sequenced precursor ions for 30 seconds. Precursor ions with +1, and +7, +8, or higher charge states were excluded from sequencing.

Citations

Step 1

Seyfried NT, Dammer EB, Swarup V, Nandakumar D, Duong DM, Yin L, Deng Q, Nguyen T, Hales CM, Wingo T, Glass J, Gearing M, Thambisetty M, Troncoso JC, Geschwind DH, Lah JJ, Levey AI. A Multi-network Approach Identifies Protein-Specific Co-expression in Asymptomatic and Symptomatic Alzheimer's Disease.

<https://doi.org/10.1016/j.cels.2016.11.006>