

Feb 01, 2024

Live-cell imaging; Calcium imaging with Fura-2

DOI

<https://dx.doi.org/10.17504/protocols.io.261genzqwg47/v1>

Minee-Liane Choi¹

¹UCL Institute of Neurology



Minee L Choi

Korea Advanced Institute of Science & Technology, The Franci...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261genzqwg47/v1>

Protocol Citation: Minee-Liane Choi 2024. Live-cell imaging; Calcium imaging with Fura-2. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.261genzqwg47/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 31, 2022



Last Modified: May 31, 2024

Protocol Integer ID: 67954

Keywords: ASAPCRN, calcium imaging, cell imaging, calcium, imaging

Abstract

This protocol describes a method to perform calcium imaging with Fura-2 (live-cell imaging)



- 1 Cells are washed with HBSS (Cat No. 14025092, ThermofisherScientific)
- 2 Cells are loaded with 5uM of Fura 2-AM (F1221, ThermofisherScientific) in HBSS.

Note

Fura 2-AM: 340/380 ratiometric measurement. Upon binding calcium, Fura-2 exhibits 340 nm and does 380 nm in the rest of the condition.

- 3 CCD protocol is set up by selecting a501 nm wheel (both 340 nm and 380 nm wavelength are excited b 501 nm).