

May 29, 2025

Live Cell Imaging

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.36wgqqz4ygk5/v1>

Camille Goldman¹

¹Icahn School of Medicine at Mount Sinai



Camille Goldman

Icahn School of Medicine at Mount Sinai, University of Oxford...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.36wgqqz4ygk5/v1>

Protocol Citation: Camille Goldman 2025. Live Cell Imaging. protocols.io

<https://dx.doi.org/10.17504/protocols.io.36wgqqz4ygk5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 15, 2025

Last Modified: May 29, 2025

Protocol Integer ID: 218390

Keywords: ASAPCRN, preparing cell, cell, cell imaging, imaging

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-024297

NASA

Grant ID: 80ARC022CA004

NIH/NINDS

Grant ID: R01NS114239

NIH/NINDS

Grant ID: UH3NS115064

NIH/NIA

Grant ID: T32AG04968

NIH/NINDS

Grant ID: F31NS13090

New York State Department of Health

Grant ID: NYSTEM-C32561GG

Abstract

Protocol for preparing cells for live-cell imaging



Preparing Astrocytes

- 1 Refer to the following protocols for midbrain astrocyte differentiation:
[10.1016/j.xpro.2021.100463](https://www.protocols.io/view/astrocyte-extraction-from-brain-organoids-261ge364wl47/v2)
<https://www.protocols.io/view/astrocyte-extraction-from-brain-organoids-261ge364wl47/v2>
- 2 Refer to the following protocols for astrocyte lentiviral transduction:
<https://protocols.io/view/astrocyte-viral-transduction-gzhtbx36p>

Plating Astrocytes

- 3 No fewer than 3 days after lentiviral transduction, and after completion of antibiotic selection, passage astrocytes into plate for imaging
- 3.1 Coat desired wells of a 96-well, black plate with cover-slip thickness plastic (Greiner [655096](#) or similar) with 0.1% gelatin
- 3.2 Wash astrocytes with 1x PBS
- 3.3 Dissociate astrocytes in a trypsin enzyme solution (TrypLE Select, Gibco 12563011). For a 6-well plate of astrocytes, use 1mL per well
- 3.4 Transfer dissociated cells to a falcon tube with 5mL astrocyte media (AM; ScienCell #1801)
- 3.5 Centrifuge at 300rcf for 3 minutes at room temperature
- 3.6 Resuspend the cell pellet in 1mL AM and count cells
- 3.7 Plate 1×10^4 cells per well of the precoated 96-well plate

Imaging Astrocytes



- 4 Approximately 3 days later, incubate the cells with a live-cell compatible nuclear dye (such as Hoechst33342 or DRAQ5) for 10 minutes
- 5 Wash the cells 1x in PBS and replace with fresh AM media
- 6 Acquire images on a microscope equipped with an incubation unit and maintain the temperature at 37°C, CO₂ at 5%, and humidity at 60%