


May 28, 2025

🌐 Live-Imaging of GCase Activity in iPSC-Derived Dopaminergic Neurons Using the PFB-FDGlu Probe

DOI

<https://dx.doi.org/10.17504/protocols.io.4r3l2br2jl1y/v1>



Roxanne Larivière¹, Jace Jones-Tabah¹, Edward A. Fon¹

¹Montreal Neurological Institute - McGill University



Roxanne Larivière

Montreal Neurological Institute - McGill University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l2br2jl1y/v1>

External link: [https://www.thermofisher.com/order/catalog/product/P11947?](https://www.thermofisher.com/order/catalog/product/P11947?ef_id=Cj0KCQjw8zABhDKARIsAHXuD7aSDcnUtkKhkbBO8g5l9jYQZgQfBoWsJqbjj9kh4NqQuoL_c8Pwr_UaAiw9EALw_wcB:G:s&s_kwcid=AL!3652!3!447292198754!!!g!!!10506731179!109642167571&cid=bid_pca_iva_r01_co_cp1359_pjt0000_bid00000_Ose_gaw_dy_pur_con&gad_source=1&gbraid=0AAAAADxi_GGSu27ji9GnhWmeiP1C-ycML&gclid=Cj0KCQjw8zABhDKARIsAHXuD7aSDcnUtkKhkbBO8g5l9jYQZgQfBoWsJqbjj9kh4NqQuoL_c8Pwr_UaAiw9EALw_wcB)

[ef_id=Cj0KCQjw8zABhDKARIsAHXuD7aSDcnUtkKhkbBO8g5l9jYQZgQfBoWsJqbjj9kh4NqQuoL_c8Pwr_UaAiw9EALw_wcB:G:s&s_kwcid=AL!3652!3!447292198754!!!g!!!10506731179!109642167571&cid=bid_pca_iva_r01_co_cp1359_pjt0000_bid00000_Ose_gaw_dy_pur_con&gad_source=1&gbraid=0AAAAADxi_GGSu27ji9GnhWmeiP1C-ycML&gclid=Cj0KCQjw8zABhDKARIsAHXuD7aSDcnUtkKhkbBO8g5l9jYQZgQfBoWsJqbjj9kh4NqQuoL_c8Pwr_UaAiw9EALw_wcB](https://www.thermofisher.com/order/catalog/product/P11947?ef_id=Cj0KCQjw8zABhDKARIsAHXuD7aSDcnUtkKhkbBO8g5l9jYQZgQfBoWsJqbjj9kh4NqQuoL_c8Pwr_UaAiw9EALw_wcB:G:s&s_kwcid=AL!3652!3!447292198754!!!g!!!10506731179!109642167571&cid=bid_pca_iva_r01_co_cp1359_pjt0000_bid00000_Ose_gaw_dy_pur_con&gad_source=1&gbraid=0AAAAADxi_GGSu27ji9GnhWmeiP1C-ycML&gclid=Cj0KCQjw8zABhDKARIsAHXuD7aSDcnUtkKhkbBO8g5l9jYQZgQfBoWsJqbjj9kh4NqQuoL_c8Pwr_UaAiw9EALw_wcB)



Protocol Citation: Roxanne Larivière, Jace Jones-Tabah, Edward A. Fon 2025. Live-Imaging of GCase Activity in iPSC-Derived Dopaminergic Neurons Using the PFB-FDGlu Probe. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.4r3l2br2jl1y/v1>

Manuscript citation:

The Parkinson's disease risk gene cathepsin B promotes fibrillar alpha-synuclein clearance, lysosomal function and glucocerebrosidase activity in dopaminergic neurons, *Mol Neurodegener.* 2024 Nov 25;19(1):88. doi: 10.1186/s13024-024-00779-9.

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 01, 2025

Last Modified: May 28, 2025

Protocol Integer ID: 204067

Keywords: GCase, iPSC, dopaminergic neurons, PFB-FDGlu, live-imaging, cellular assay, imaging of gcase activity, cell analysis of gcase activity, assessing gcase activity, active gcase enzyme, fdglu gcase substrate, gcase activity, cell imaging, derived dopaminergic neuron, dopaminergic neuron, fluorescence detection, fdglu probe, cell imaging with the opera phenix, fluorescent dye, cell analysis, lysosomal compartments upon cellular uptake, fdglu, fdglu probe this protocol, lysosomal compartment, cellular uptake, cell viability, cell, neuron

Funders Acknowledgements:

Michael J. Fox Foundation for Parkinson's research



Grant ID: MJFF-020696



Abstract

This protocol describes the procedure for assessing GCase activity in iPSC-derived dopaminergic neurons using the PFB-FDGlu GCase substrate. The assessment is performed via live-cell imaging with the Opera Phenix high-content imaging system and analyzed using Harmony software.

PFB-FDGlu is a cell-permeable compound that localizes to endosomal and lysosomal compartments upon cellular uptake, without compromising cell viability. When cleaved by active GCase enzymes, it releases a green-fluorescent dye (PFB-F), enabling specific, single-cell analysis of GCase activity through fluorescence detection.



(A) Representative images of Cell Line 1 from the GCase activity assay in CBE-treated and untreated iPSC-derived dopaminergic neurons differentiated for 6 weeks. **(B)** Corresponding GCase activity curves.



Materials

Reagents	Company	Cat. Number
LysoTracker Deep Red	ThermoFisher Scientific	L12492
PFB-FDGlu (5-(Pentafluorobenzoylamino)Fluorescein Di- β -D-Glucopyranoside)	ThermoFisher Scientific	P11947
FluoroBrite DMEM	ThermoFisher Scientific	A1896701
CBE (Conduritol B epoxide)	Millipore Sigma	C5424
DMSO	FisherScientific	BP231-100

Safety warnings

- ⚠ Be extremely gentle when dispensing into the wells, as dopaminergic neurons are prone to detaching with forceful pipetting

Before start

iPSC-derived dopaminergic neurons are cultured on 96-well plates for up to 4 to 6 weeks of differentiation. We use black Costar 96-well assay plates (Corning cat. no. 3904), which are compatible with imaging using the Opera Phenix high-content screening system.

One hour prior to imaging, turn on the Opera Phenix high-content screening system. Set the incubation chamber temperature to 37 °C and adjust the CO₂ concentration to 5%.

CBE treatment

16h

- 1 For the negative control, incubate neurons  Overnight at  37 °C with  25 micromolar (μM) of the GCase inhibitor Conduritol B epoxide (CBE)

16h

Labeling lysosomes with LysoTracker Deep Red

30m

- 2 Prepare a LysoTracker Deep Red working solution at  50 nanomolar (nM) in media. Remove the existing media from the wells and add  60 μL of the LysoTracker solution. Incubate for  00:30:00 at  37 °C .

30m

Preparing GCase substrate stock solution

- 3 While the cells are incubating, prepare the GCase substrate (PFB-FDGlu) stock solution. Reconstitute one  5 mg vial of PFB-FDGlu in  154 μL of DMSO for a  37.5 millimolar (mM) stock, then aliquot into  10 μL single-use stock vials. Store aliquots at  -20 °C , protected from light, and avoid repeated freeze-thaw cycles.

Addition of the GCase substrate

- 4
 1. Prepare the PFB-FDGlu working solution by diluting the stock 1:200 in FluoroBrite DMEM.
 2. Carefully remove the LysoTracker Deep Red solution from the wells and gently add  50 μL of the PFB-FDGlu working solution to each well.
 3. Proceed with imaging

Imaging using the Opera Phenix high-content imaging system

2h 30m

- 5
 1. Insert the prepared plate into the Opera Phenix chamber.
 2. Utilize a 40 $\times$ water immersion objective. Capture images using the 488 nm and 647 nm laser channels. Acquire 15-18 fields of view per well every  00:15:00 for a

2h 30m



total duration of  02:30:00 . Maintain the plate under standard cell culture conditions ( 37 °C , 5% CO₂) throughout the imaging session.

Image Analysis Workflow in Harmony Software

- 6
 1. **Load Images**
Open the *Data Analysis* tab in the Harmony software and load the desired image dataset.
 2. Create an Analysis Sequence
 3. **Identify Cells**
Use the **Alexa 647** channel with the **P method** to detect cells. Parameters: Area > 100 µm², Splitting Sensitivity: 0.70, Common Threshold: 0.55). Save this population as "**Cells**".
 4. **Select Cell Population**
Apply Common Filters to the "Cells" population to refine it, creating a new population named "**Cells Selected**".
 5. **Calculate Morphological Properties**
For the "Cells Selected" population, use the Standard Method to calculate the **Area**.
 6. **Calculate Intensity Properties**
 - For the **Alexa 647** channel, use the Standard Method to calculate the **mean intensity** for the "Cells Selected" population.
 - For the **Alexa 488** channel, use the Standard Method to calculate the **mean intensity** for the "Cells Selected" population.
 7. **Define Results Output**
 - Output: **Cells Selected - Intensity Cell Alexa 647 Mean**
 - Output: **Cells Selected - Intensity Cell Alexa 488 Mean**
 8. **Run the Analysis**
 9. **Export Data**
Download the results as a **CSV file**.

Assess GCase activity

- 7
 - From the table generated by Harmony, calculate the ratio of the PFB-FDGlu mean fluorescence intensity to the LysoTracker Deep Red (647) baseline reading (at t=0) for each time point.
 - Enter these calculated values into Prism to generate the GCase activity curves.