

Dec 19, 2019

LIVE IMAGING OF i³NEURONS (Support Protocol 5)

 In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.5w7g7hn

iPSCs

Michael S. Fernandopulle¹, Ryan Prestil¹, Christopher Grunseich¹, Chao Wang², Li Gan², Michael E. Ward¹

¹National Institute of Neurological Disorders and Stroke, National Institutes of Health, Bethesda, Maryland;

²Gladstone Institute of Neurological Disease, Gladstone Institutes, San Francisco, California

Neurodegeneration Method Development Community

Tech. support email: ndcn-help@chanzuckerberg.com



Julia Rossmanith

protocols.io

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.5w7g7hn

External link: <https://doi.org/10.1002/cpcb.51>

Protocol Citation: Michael S. Fernandopulle, Ryan Prestil, Christopher Grunseich, Chao Wang, Li Gan, Michael E. Ward 2019. LIVE IMAGING OF i³NEURONS (Support Protocol 5) . **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5w7g7hn>

Manuscript citation:

Fernandopulle, M. S., Prestil, R., Grunseich, C., Wang, C., Gan, L., & Ward, M. E. (2018). Transcription-factor mediated differentiation of human iPSCs into neurons. Current Protocols in Cell Biology, e51. doi:<https://doi.org/10.1002/cpcb.51>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 27, 2019

Last Modified: December 19, 2019

Protocol Integer ID: 26303

Keywords: i3LMN, i3Neurons, iPSC, iPSC-derived neurons, transcription factor-mediated differentiation

Abstract

Live imaging permits visualization of molecular and organellar dynamics within the neuron. While a standard confocal microscope is sufficient for short imaging experiments, extended imaging applications (>1 hr) are best served by a 37 °C live imaging chamber outfitted onto the microscope. CM should also be changed to Hibernate A Low Fluorescence Medium (BrainBits LLC, [cat. no. SKU#HAPR](#)) for extended imaging.

Attachments





[fernandopulle2018.pd...](#)


1.7MB

Guidelines

Live imaging permits visualization of molecular and organellar dynamics within the neuron. While a standard confocal microscope is sufficient for short imaging experiments, extended imaging applications (>1 hr) are best served by a 37 °C live imaging chamber outfitted onto the microscope. CM should also be changed to Hibernate A Low Fluorescence Medium (BrainBits LLC, [cat. no. SKU#HAPR](#)) for extended imaging. Hibernate A permits long-term maintenance of neuronal cultures in ambient carbon dioxide levels (0.04 % vs. 5 % for standard cell culture incubators) and provides a better imaging environment by reducing autofluorescence from phenol red-containing medium. Finally, medium (either CM for short imaging or Hibernate A for extended imaging) should be supplemented with SOS neuronal supplement (Cell Guidance Systems, [M09-50](#)) instead of B27. SOS supplement does not contain phototoxic components present in B27 and other neuronal supplements. Imaging is best done on glass-bottom slides, such as Ibidi μ -slides (Ibidi, [cat. no. 80827](#)).

 Hibernate A Low Fluorescence **BrainBits Catalog #HALF**

 SOS® neuronal supplement **Catalog #M09-50**

 μ -Slide 8 Well Glass Bottom **Ibidi Catalog #80827**

Safety warnings

 Please see SDS (Safety Data Sheet) for hazards and safety warnings.

