

Jul 31, 2023

Live imaging to investigate mitophagy kinetics and NEMO recruitment in HeLa-M cells

DOI

dx.doi.org/10.17504/protocols.io.rm7vzy9w5lx1/v1

OLIVIA HARDING^{1,2}, Erika L.F. Holzbaur^{1,2}

¹Department of Physiology, Perelman School of Medicine, University of Pennsylvania, Philadelphia, PA 19104;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815

Liv



OLIVIA HARDING

Children's Hospital of Philadelphia

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzy9w5lx1/v1>

Protocol Citation: OLIVIA HARDING, Erika L.F. Holzbaur 2023. Live imaging to investigate mitophagy kinetics and NEMO recruitment in HeLa-M cells. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.rm7vzy9w5lx1/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 22, 2022

Last Modified: May 31, 2024

Protocol Integer ID: 65074

Keywords: Tissue culture, Mitophagy, Fluorescent ligands, Confocal microscopy, Live cell imaging, ASAPCRN, damaged mitochondria, live cell imaging, tagged mitochondria, mitochondrial depolarization, substitute for live cell imaging, mitochondria with less correlation, live imaging of nemo occupancy, microscopy, p62 during the cellular response, mitophagy kinetic, subcellular biology, live imaging, confocal microscopy system, imaging, live cell, cell, cellular response

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: Mechanisms of mitochondrial damage control by PINK1 and Parkin (ASAP-000350)

Abstract

There is no substitute for live cell imaging in the investigation of the kinetics of subcellular biology. Here, we enumerate a protocol to visualize fluorescently-tagged mitochondria, NEMO, OPTN, and p62 during the cellular response to mitochondrial depolarization. Because live cells are sensitive to photo-damage, we used low laser power and short exposure times in our confocal microscopy system. Even with minimally damaging parameters, we were able to collect high-content data that we subsequently analyzed. With this technique, we showed that NEMO and OPTN, despite containing highly similar domains, were recruited to damaged mitochondria with less correlation than NEMO and p62. Furthermore, live imaging of NEMO occupancy on damaged mitochondria was a necessary complement to our parallel fixation studies. Since fixation introduces artifacts, especially in samples with concentrated proteins like those in mitophagy, results from our live cell imaging corroborated our findings in fixation studies. Our reporting of these results would not have been possible without real-time, live cell imaging.

Attachments



[470-983.pdf](#)

251KB

Guidelines

- This protocol was developed with the HeLa subtype, HeLa-M. HeLa-M cells are flatter than standard HeLa cells, making them easier to image. They also uptake siRNA better than standard HeLa. Regardless, the protocol would be easily adaptable to standard HeLa cells or other cell culture lines.
- This protocol was created in order to investigate Parkin-dependent mitophagy. Parkin and several other fluorescently-tagged mitophagy components are intended for use in the protocol.



Materials

Materials/Reagents

✂ 1.5 mL capped tubes **Merck MilliporeSigma (Sigma-Aldrich) Catalog #EP022364120**

Reagents

- FBS (HyClone)
- ✂ Leibovitz's L-15 Medium **Thermo Fisher Catalog #11415064**
- ✂ Dimethyl sulfoxide **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D2650**
- Ethanol
- ✂ Antimycin A from *Streptomyces* sp. **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A8674**
- ✂ Oligomycin A **Merck MilliporeSigma (Sigma-Aldrich) Catalog #75351**

Equipment:

- Vacuum apparatus
- 🌡 37 °C water bath
- Confocal microscope with 60X objective, heated stage chamber, focus correction system, and associated software

Note

Our system is not equipped with a pressurized stage chamber to sustain 5% CO₂ conditions. Thus, we use L-15 in order to buffer the samples in atmospheric CO₂ conditions.

Troubleshooting



Before start

- The start point for this protocol is after cells grown on \pm 35 mm glass bottom dishes have been transfected with Mito-dsRed2, EGFP-NEMO, HaloOPTN, and Parkin for 18:00:00 - 24:00:00 and tagged with Halo ligand.
- Prepare 45 millimolar (mM) stock of Antimycin A by suspending 50 mg solid AntA in 2 mL ethanol.
- Prepare 10 millimolar (mM) stock Oligomycin A by suspending 5 mg solid OligA in 630 μ L DMSO.
- Prepare imaging media by making a 10% FBS solution in L-15 and warm in water bath.






Note

Will use ~ 3 mL imaging media per dish.

- Prepare working AntA/OligA solution by adding 1 μ L 10 millimolar (mM) OligA and 4 mL 45 millimolar (mM) AntA to 0.5 mL Imaging media in a 1.5 mL tube. Keep warm in bath or imaging chamber.
- Heat microscope imaging chamber to 37 $^{\circ}$ C .

Replace Standard Media with Imaging Media

10m

- 1 Wait until imaging chamber is heated to  37 °C .
- 2 Aspirate media from sample.
- 3 Repeat the following two steps for a total of 2 washes.
 - Add  200 μ L Imaging media gently to dish.
 - Aspirate.
- 4 Add  1.5 mL Imaging media and place dish in  37 °C imaging chamber.
- 5 Incubate sample in imaging conditions for at least  00:10:00 before imaging.



10m




Imaging set-up

15m

- 6 Place dish on microscope stage, moving stage adjusters if necessary to fit the dish.

Note

Remove the lid of the dish to maneuver it more easily without spilling its contents. Replace the lid after dish is firmly secured on the stage.

- 7 Raise objective so that there is an oil interface between the objective lens and glass bottom of the dish.
- 8 Allow dish to settle at least  00:05:00 in this position before imaging to minimize drift during imaging.
- 9 Using 60X objective and RFP epifluorescence, find the focal plane of fixed cells by looking for dsRed2-labeled mitochondria.

5m

**Note**

We find that 60X is sufficient magnification to collect several cells in each field of view with enough resolution to perform analysis and quantification. Other magnifications may be appropriate for various applications.

Note

Note the health of cells, transfection efficiency of the constructs, and brightness of fluorescence. If you will image multiple samples with varying conditions, you may want to observe all dishes before confirming the imaging parameters so that no images are overexposed or too dim.

- 10 Configure 488, 561, and 647 lasers and accompanying exposure times.


Note

Set parameters as low as possible while still detecting the signals in order to avoid phototoxicity and bleaching. For a benchmark, we use <10% power and <200 ms exposure for each channel. Save settings in order to maintain consistency among experiments.

STEP CASE**Imaging confocal (single z section)** 4 steps

- 11

Note

This option is used to collect information in single sections of the sample over the course of 1-  02:00:00 . Data can be used to reconstruct recruitment kinetics, however single section imaging is not ideal for assays that require a high signal to-noise ratio.

Choose a field of view with multiple cells that look healthy and express all three fluorescent constructs. Set focus on a confocal section approximately between 5 and 7, where 0 is the bottom, attached part of the cell and 10 is the top of the cell.

**Note**

If scope is equipped with multi-stage position feature, choose several fields that fulfill these requirements and add to sequence.

12 Begin imaging at acquisition rate of 1 frame per min.



13 After 5-10 frames (00:05:00 - 00:10:00), remove the lid of the imaging dish and add the working AntA/OligA solution for a total of 2 mL (5 micromolar (μ M) AntA/ 10 micromolar (μ M) OligA).

15m

Note

Record the frame number at which AntA/OligA was added.

14 Replace the lid to the dish and continue imaging for as long as desired.

**Note**

- OPTN recruitment is visible between 30 and 01:00:00 .
- After 01:00:00 imaging at 1 frame per minute, even low laser power settings will damage the cells. Take caution in analyzing data from samples after 01:00:00 imaging.