

Dec 03, 2019 Version 2

# 🌐 Loop L1 (odd level) Bsal type IIS cloning into pCk vectors V.2

🔗 Version 1 is forked from [Loop L1 \(odd level\) type IIS cloning - pCk-ye vectors](#)

DOI

[dx.doi.org/10.17504/protocols.io.92gh8bw](https://dx.doi.org/10.17504/protocols.io.92gh8bw)

Eftychis Frangedakis<sup>1</sup>, Susana Sauret-Gueto<sup>2</sup>, Anthony West<sup>3</sup>, Nicola Patron<sup>4</sup>, Marta Marta Tomaselli<sup>2</sup>, Marius Rebmann<sup>2</sup>, Jim Haseloff<sup>2</sup>

<sup>1</sup>University of Cambridge; <sup>2</sup>Plant Sciences, University of Cambridge, OpenPlant;


<sup>3</sup>previously at Earlham Institute, Norwich; <sup>4</sup>Earlham Institute, Norwich

OpenPlant Project



**Susana Sauret-Gueto**

Plant Sciences, University of Cambridge, OpenPlant



OPEN ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.92gh8bw](https://dx.doi.org/10.17504/protocols.io.92gh8bw)

**Protocol Citation:** Eftychis Frangedakis, Susana Sauret-Gueto, Anthony West, Nicola Patron, Marta Marta Tomaselli, Marius Rebmann, Jim Haseloff 2019. Loop L1 (odd level) Bsal type IIS cloning into pCk vectors. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.92gh8bw>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** December 03, 2019

**Last Modified:** December 03, 2019

**Protocol Integer ID:** 30504



## Abstract

Protocol based on:

**Pollak B, Cerda A, Delmans M, et al (2019) Loop assembly: a simple and open system for recursive fabrication of DNA circuits. New Phytol 222:628–640**

<https://doi.org/10.1111/nph.15625>

## Materials

### MATERIALS

⊗ Bsal - 5,000 units **New England Biolabs Catalog #R0535L**

⊗ T4 DNA Ligase - 20,000 units **New England Biolabs Catalog #M0202S**


⊗ Sterile water

⊗ BSA, molecular biology grade, 20 mg/ml **New England Biolabs Catalog # B9000S**

⊗ 10X NEB T4 DNA ligase buffer **New England Biolabs**

## Loop pCk vectors

1



### Loop vectors for nuclear transformation: pCks (A) and for chloroplast transformation pCkchlo (B).


Loop fusion sites in the pCk vectors to assemble different L0 parts into a L1 construct using a pCk vector and BsaI are: A (GGAG) and F (CGCT).

Loop fusion sites in the pCk vectors to assemble different L1 constructs into a L2 construct using a pCs vector and SmaI are: a (ATC), b (GCA), d (TAC), e (CAG) and o (GCT).

Left (LB) and right border (RB) repeats from nopaline C58 T-DNA for Agrobacterium-mediated nuclear transformation. KanR: kanamycin bacterial resistance cassette. LacZ: lacZa cassette for blue-white screening of colonies.

## Example of assembly of L0 parts into a transcription unit (L1)

2



Loop assembly of multiple L0 parts into a transcription unit (L1) using a pCk plasmid and BsaI.

## Protocol for assembly of L0 parts into a transcription unit (L1)

- Determine the concentrations of each DNA plasmid needed (L0 plasmids and pCk acceptor plasmid) by spectrophotometry (Nanodrop).  
In the example in step 2, determine concentration of plasmids L0\_PROM5, L0\_CDS, L0\_3TERM and pCk1.
- Prepare aliquots for each plasmid at a concentration of 15 nM for the L0 plasmids and of 7.5 nM for the acceptor pCk vector. With this final concentration, 1  $\mu$ L of each plasmid is added to the plasmids mix (see step 6).

To calculate the concentration in ng/ $\mu$ L:

- For a final concentration of 15 nM, the concentration in [ng/ $\mu$ L] equals N (the length in bp of the plasmid) divided by 110. This is an approximation of the formula:

$15 \cdot 10^{-9} \text{ mol/L} \times ((607.4 \times N) + 157.9) \text{ g/mol} \times 10^{-6} \text{ L}/\mu\text{L} \times 10^9 \text{ ng/g} = \text{concentration (ng}/\mu\text{L)}$



- For a final concentration of 7.5 nM, the concentration in [ng/ul] equals N divided by 220.

- 5 Prepare Loop assembly Level 1 reaction master mix (MM) according to Table , if four or less number of L0 parts are assembled into a pCk vector (otherwise see step 8)

Com pone nts	Volu me ( $\mu$ L)
Steril e water	3
10x T4 ligase buffer (NEB)	1
1 mg/mL bovin e seru m albu min (NEB)	0.5
T4 DNA ligase at 400 U/ $\mu$ L (NEB)	0.25
10 U/ $\mu$ L BsaI (NEB)	0.25
Final volum e	5



- 6 Prepare plasmids mix by adding in a 0.2 mL tube: 1  $\mu$ L of each L0 plasmid , 1  $\mu$ L of the pCk vector (see step 4), and sterile water up to 5  $\mu$ L. Mix well.  
  
When 4 L0 parts are assembled into a pCk plasmid, the volume of the plasmid mix is 5  $\mu$ L, and thus no volume of water is added.
- 7 Add 5  $\mu$ L of MM (step 5) to the 5  $\mu$ L of plasmids mix (step 6), to a final volume of 10  $\mu$ L. Mix well.
- 8 If more than 4 L0 parts are to be assembled into a pCk vector, reduce the water volume in the MM by 1  $\mu$ L (step 5) for each extra 1  $\mu$ L of DNA part added in the plasmids mix (step 6).
- 9 Place samples in a thermocycler and use the following program:  
Assembly: 26 cycles of 37 °C for 3 min and 16 °C for 4 min.  
Termination and enzyme denaturation: 50 °C for 5 min and 80 °C for 10 min.
- 10 Transform 20  $\mu$ L of chemically competent E. coli cells (transformation efficiency of  $1 \times 10^7$  transformants/ $\mu$ g plasmid DNA) using 2  $\mu$ L of the Loop assembly reaction and then plate on LB agar plates containing 50  $\mu$ g/mL kanamycin and 40  $\mu$ g/mL of X-gal for blue-white screening.
- 11 Incubate overnight at 37 °C.
- 12 Colonies with white color are likely to contain an L1 insert cloned into the pCk vector (In the example in step 2: PROM5:CDS-3TERM)  
Blue color colonies will contain undigested pCk vector with LacZ
- 13 Confirm the presence of the correct insert with Sanger sequencing using the primers pC\_F (GCAACGCTCTGTCATCGTTAC) and pC\_R (GTAAGTTAGGACTTGTGCGACATGTC) for pCk vectors, and pC\_F and pC\_R2 (CAATCTGCTCTGATGCCGCATAGTTAAG) for pCkchlo vectors.